Service Repair Sheet SRS68519

Contact Name

Company/ Hospital Name

Department

Position

Direct Phone

General Phone

Opera Account

Email

Order Number

Date Received

Booked in By

Main Company

Type Return

Andrew Tracey

Smit Mobile Equipment

Service Engineer

07917170368

0127664490

00020025

17/Jul/2023

Robert Connor

Viamed

Quote

VIAMEDclear Goods In Only Decontamination certificate provided by customer / Cleaned by Viamed, I if no declaration certificate from customer Signed: Date: Goods Out Only Cleaned by Viamed before returning to customer Signed: Date:

Notes 17/Jul/2023 Sophie Lines

17/Jul/2023 Sophie Lines

Has an MX300 which is displaying a fault - he is driving past today (17th July) and would like to drop it off for us to post once sorted. S/N: 332862.

No email address on system for him - going to phone and ask for it

17/Jul/2023 Sophie Lines

He couldn't remember his work email address - private one is a tracy@ntlworld.com (confirmed no e in tracy in email address)

17/Jul/2023 Robert Connor

Received 1 x MX300i s/n 332882, with 3 x AA battery.

Ready For quote

Canean 25-7-23 Canean 1-8-23

Repair Complete Signed

| SRN | Equipment | Stock Ref | Serial Number | Warranty |
|----------|----------------|-----------|---------------|----------|
| SRN35931 | Oxygen Monitor | 0111235 | 332882 | |

0140005 x 0.5 @ E63 per Mr. SIN, SRS, SRN 0131214 X 1 0 £2.40 SRS, SRN 710 R12.

Service Repair Sheet 68519

Contact Name

Company/ Hospital Name

Position
Direct Phone
General Phone
Opera Account
Order Number

Date Received

Booked in By

Andrew Tracey

Smit Mobile Equipment

Service Engineer 07917170368 0127664490 00020025

MP10055 17/Jul/2023 Robert Connor

Repair

SRN35931

Ref

S/N

0111235 332882

Equipment Type

Oxygen Monitor

Complete - Repaired Time :0.5 Hour(s)

Parts Replaced

Qty: 1 0131214 Qty: 0.5 0140005

The MX300 Oxygen monitor has been repaired. The sensor input connector has been replaced.



Smit Mobile Equipment (UK) Ltd.
Unit 17, Nelson Way
Tuscam Trading Estate
Camberley
Surrey GU15 3DH

Telephone: +44 (0) 1276 64490 Fax: +44 (0) 1276 63350

PURCHASE ORDER

| Name and Address of Supplier: | Purchase Order No.: MP 10055 |
|-------------------------------|------------------------------|
| VIAMED | Date: 27/07/2023 |
| | Project / ES Ref: No. |
| Attn: JANINE GILL | RELOC 2 |
| | |

Delivery Date: SUTNEST

Delivery Point (if different to above address): ANDREW TRACEY

34 KINTORE DRIVE, GREAT SANKEY, WARRING TON. WASSIF

| Qty | Description | Unit Price* | Total Price |
|-----|--------------------------------|-------------|-------------|
| | REPAIR OF MX300 OXYGEN MONITOR | | |
| | 040005 LABOUR | 31-50 | |
| | 0131214 DUPUT CONNECTOR | 2-40 | |
| | SHIPPING | 12-00 | |
| | REF: 58568519 | | |
| | Total Value | 1 | 45-90 |

| | Total Value 45-90 |
|--|---|
| Signed for Smit Mobile Equipment (UK) Ltd. | Name of Person Signing (in Capitals) MARK HALS |
| Please specify pric | e where possible |

Price exclusive of VAT, but includes all other charges. This order is subject to our Terms & Conditions. All packages, delivery notes and invoice must bear above order number.

Registered in England No. 4184642

VAT Reg. No. 775 461502