





Spectrometer:	Prema SPM 9001
Manufacturer of probe:	UDT
Part number:	P858RA
Serial number:	8A02877
Cable length:	8 ft
Connector:	9 pin D male - pins 4 & 8 missing
Connects to:	Extension cable
For use on:	Nellcor
Tested on:	Nellcor N-200
Red peak wavelength:	654
Infra red peak wavelength:	876.5
LED configuration:	2 lead



Spectrometer: Prema SPM 9001

Manufacturer of probe: UDT

Part number: P858RA

Serial number: 8D05614

Cable length: 8 ft

Connector: 9 pin D male - pins 4 & 8 missing

Connects to: Extension cable


For use on: Nellcor

Tested on: Nellcor N-200

Red peak wavelength: 658

Infra red peak wavelength: 873.3

LED configuration: 2 lead



Spectrometer: Prema SPM 9001


Manufacturer of probe: UDT
Part number: P858RA
8005615

Connects to: Extension cable

For use on: Nellcor
Tested on: Nellcor N-200

Red peak wavelength: 658
Infra red peak wavelength: 870

LED configuration: 2 lead



Spectrometer: Prema SPM 9001

Manufacturer of probe: UDT

Part number: P858RA

Serial number: 8D05616

Cable length: 8 ft

Connector: 9 pin D male - pins 4 & 8 missing

Connects to: Extension cable


For use on: Nellcor

Tested on: Nellcor N-200

Red peak wavelength: 658

Infra red peak wavelength: 873.3

LED configuration: 2 lead



Spectrometer: Prema SPM 9001

Manufacturer of probe: UDT

Part number: P858RA

Serial number: 8D05617

Cable length: 8 ft

Connector: 9 pin D male - pins 4 & 8 missing

Connects to: Extension cable

For use on: Nellcor

Tested on: Nellcor N-200

Red peak wavelength: 658

Infra red peak wavelength: 870

LED configuration: 2 lead