

Viamed repair

WATSON, Becky (ROYAL DEVON UNIVERSITY HEALTHCARE NHS FOUNDATION TRUST)
<becky.watson@nhs.net>
To: Janine Gill <office@viamed.co.uk>

19 January 2024 at
08:35

Good morning Janine.

The MX300 will not calibrate. When turned on, it gives a display of 00.0 and when going through calibration with O2, it counts down the calibration, but does not register a reading. The SN is 388702 .

The O2 cell has been replaced.

This handset has probably been left unused for a few years and came into the workshop with flat batteries.

Regards
Becky

Becky

Rebecca Watson
Senior Clinical Technologist - Specialist Services
Ex 2060

From: Main Account <viamedinbox@gmail.com> **On Behalf Of** Janine Gill
Sent: 17 January 2024 12:37
To: WATSON, Becky (ROYAL DEVON UNIVERSITY HEALTHCARE NHS FOUNDATION TRUST)
<becky.watson@nhs.net>
Subject: Viamed repair

You don't often get email from office@viamed.co.uk. [Learn why this is important](#)

This message originated from outside of NHSmail. Please do not click links or open attachments unless you recognise the sender and know the content is safe.

[Quoted text hidden]

This message may contain confidential information. If you are not the intended

recipient please:

- i) inform the sender that you have received the message in error before deleting it; and
 - ii) do not disclose, copy or distribute information in this e-mail or take any action in relation to its content (to do so is strictly prohibited and may be unlawful).
- Thank you for your co-operation.

NHSmial is the secure email, collaboration and directory service available for all NHS staff in England. NHSmial is approved for exchanging patient data and other sensitive information with NHSmial and other accredited email services.

For more information and to find out how you can switch visit [Joining NHSmial – NHSmial Support](#)