
pulse oximetry



MD300-C

Finger Pulse Oximeter

Ultra-compact, lightweight pulse oximeter with LED display

Low cost

Easy to use; 1-button operation

Automatic power-off when finger is removed

Economical pulse oximeter for general applications



Actual Size

The MD300-C can be used in a variety of applications, from use in hospitals, paramedic ambulances, clinics and dental surgeries, through to spot-checking before, during and after exercise.





MD300-C Finger Pulse Oximeter

Measuring Range

SpO₂ Pulse Rate 70 - 100% 30 - 250 bpm Accuracy SpO₂ (70 - 100%) SpO₂ Resolution ± 2 digits 1% Pulse Rate; beats per minute (30 - 99 bpm) ± 2 bpm Pulse Rate; beats per minute (100 - 250 bpm) ± 2% Pulse Rate Resolution 1 bpm

Display

Features LED display

Data displayed Oxygen saturation, pulse rate, pulse indicator, bar graph

Indicators Signal strength, battery status

Environmental Conditions

5°C to 40°C Operating temperature -20°C to +55°C Storage temperature

Regulatory Approvals

Standards CE Class IIa

Assembly **Dimensions** Weight Power Supply Battery Life

56 mm (L) x 32 mm (W) x 30 mm (H) Approx. 50 g (including batteries) 2 x AAA/LR03/MN2400 1.5V alkaline batteries Approximately 16 hours of continuous operation. Automatically powers off when finger is removed for longer than 8 seconds.

Part Numbers

MD300-C EAN-13 (GTIN-13) Barcode: 5051826003046 2810010

Standard accessories included:

Lanyard

2 x AAA Batteries

Operator Manual

Optional Accessories	
Soft Carry Case, with belt loop, colour - black	0022190
Soft Carry Case, with belt loop, colour - orange	0022191
Soft Carry Case, with belt loop, colour - green	0022192
Soft Carry Case, with belt loop, colour - blue	0022193

Clinical trials carried out in accordance with BS EN ISO 9919:2005 standards. Materials tested in accordance with ISO 10993-, -5, -10.







Optional carrying cases with belt loop and D-ring

C€0123

Specifications Subject to Change

Distributed by:



VIAMED

Viamed Limited · 15 Station Road · Cross Hills Keighley · West Yorkshire · BD20 7DT · United Kingdom Tel: +44 (0)1535 634 542 Fax: +44 (0)1535 635 582 Email: info@viamed.co.uk Website: www.viamed.co.uk





ISO 13485:2003



Part Number: 2890099