

3-FEB-00 DO 7:35

00410319514673
ING. NUFER AG, GUEMLIGEN

00410319514673

PL ✓ S. 01



Ing.Nufer AG
Schürmattstrasse 6
Postfach 125
CH-3073 Gümligen

Tel. ++41 31.958 66 66 Fax ++41 31 951 46 73
E-Mail: info@nufer-medical.ch
Internet: <http://www.nufer-medical.ch>
UBS Bern
MWSt-Nr. 121 519



FAX MESSAGE

Gümligen, 2 February 2000	Pages: 1 incl. cover
To: Viamed Limited, Keighley	
Fax No.: 0044 1535 635 582	
Attention: Mr Peter Lamb	

Dear Peter,

Further to your fax of today we would like to confirm that the mobile stand has been tested to remain stable in any position on a 10 degree slope.

The worst position is of course when the warmer itself is swung out in a 45° angle with the horizontal arm of the stand swung out by 45° in the same direction but the stand does remain stable on a 10° slope even in this position. (However, in this position all three ends of the stand would bear a high risk of stumbling over which we consider a far greater risk to anyone involved.)

We do of course recommend that the warmer is used in the straight or near straight position because of the design of the stand and to swing out the horizontal arm only in order to position the warmer parallel or across the area to be heated (to achieve this position one usually has to swing the horizontal arm in one direction and the warmer itself in the other).

In any case the stand has been approved to remain stable on a 10 degree slope (but we would like to stress that it is naturally still possible to make it tumble on purpose, e.g. by pushing it).

We hope to have answered your question satisfactorily but remain at your disposal for any further information which you may require.

Best regards,

ING. NUFER AG

Dieter Indermühle
Product Manager

PS: Thank you for displaying the Ceratherm at the Egypt Med 2000. Will you be visiting other fairs, too. Please let us know. Thank you