

Viamed Limited - 15 Station Road - Cross Hills Keighley · West Yorkshire · BD20 7DT · United Kingdom Tel: +44 (0)1535 634 542 · Fax: +44 (0)1535 635 582 Email: info@viamed.co.uk Website: www.viamed.co.uk

Service Repair Sheet 68742

Contact Name

Company/ Hospital Name Department

Position Direct Phone General Phone Opera Account

Email

Order Number Date Received Booked in By

David Watson

**Darlington Memorial Hospital** Clinical Engineering Department Clinical Engineering Technician

01325743048 01325380100 00001350

david.watson28@nhs.net

CDD0005996494 02/Apr/2024 Robert Connor

Repair

S/N SRN36350 0310030 VMD105

Equipment Type Tom Thumb

Complete - Repaired Time: 0 Hour(s)

Parts Replaced

0380000 Qty: 1

Tom Thumb TT480 VMD105 has been stripped and cleaned then reassembled with new O rings. The two valves have been checked and reset before QA.