



Derek Lamb <liquidgands@gmail.com>

Fwd: Lifespan of Products

1 message

Ryan Swaine <ryan.swaine@viamed.co.uk>
To: Derek Lamb <derek.lamb@viamed.co.uk>

27 March 2015 at 13:51

Ryan Swaine
International Sales Manager
VIAMED Ltd.
<http://www.viamed-online.com>
Email ryan.swaine@viamed.co.uk
Tel: +44 (0)1535 634542
Fax: +44 (0)1535 635582
Mob: +44 (0)7803 907117

----- Forwarded message -----

From: **Froehlich, Heidi** <h.froehlich@bluepoint-medical.com>
Date: 27 March 2015 at 11:20
Subject: AW: Lifespan of Products
To: Ryan Swaine <ryan.swaine@viamed.co.uk>
Cc: Steve Nixon <steve.nixon@viamed.co.uk>, "Lindner, Bernd" <b.lindner@bluepoint-medical.com>

Hi Ryan,

Do you mean with "lifespan" the "service life"? See definition according to 60601-1 3rd ed and also obligatory to be stated in the User Manual for Home Care use below:

60601-1 © IEC:2005

- 53 -

3.28**EXPECTED SERVICE LIFE**

maximum period of useful life as defined by the MANUFACTURER

Just confirmed by Bernd the following Service life will apply:

- 5 years for all monitors (also VM-2500-S, VM-2500-M, VM-2160, VM-2101 etc.)

- 5 years also for the reusable SpO2 sensors
- Spirotrue A flow sensor are disposable. Shelf life according to the label

Let me know if you have questions

Best regards,

Heidi

Von: viamed.ryan.swaine@gmail.com [mailto:viamed.ryan.swaine@gmail.com] **Im Auftrag von** Ryan Swaine

Gesendet: Freitag, 27. März 2015 11:41

An: Froehlich, Heidi

Betreff: Lifespan of Products

Hi Heidi

I hope that you can help me or point me in the right direction.

For registration purposes I need the specified lifespan of the following products:

VM-2500-S

VM-2500-M

VM-2160

All the reusable SPO2 sensors for use with the above monitors

VM-2101

Spirotrue A flow sensors

I look forward to your reply.

Best regards

Ryan

Ryan Swaine
International Sales Manager
VIAMED Ltd.
<http://www.viamed-online.com>
Email ryan.swaine@viamed.co.uk
Tel: +44 (0)1535 634542
Fax: +44 (0)1535 635582
Mob: +44 (0)7803 907117