

MG Full, CS Full



Post market Surveillance and Risk Assesment Re-evaluation

Tom Thumb

[Last Reviewed by Derek Lamb on 03/05/24](#)
[\(templink print file description\)](#)

Associated Stock References

Stock Reference	Description
0310027	Tom Thumb Infant Resuscitator (Air).
0310029	Tom Thumb Infant Resuscitator (Air).
0310030	Tom Thumb Infant Resuscitator.
0310031	Tom Thumb Infant Resuscitator (Oxygen).
0310032	Tom Thumb Infant Resuscitator (Oxygen).
0310033	Tom Thumb Infant Resuscitator (Oxygen).
0310034	Tom Thumb Infant Resuscitator (Oxygen).
0310072	Tom Thumb Infant Resuscitator.
0310074	Tom Thumb Infant Resuscitator.
0310075	Tom Thumb Infant Resuscitator.
0310076	Tom Thumb Infant Resuscitator.
0310080	Tom Thumb Infant Resuscitator (RO).
0310082	Tom Thumb Infant Resuscitator (RO).
0310092	Tom Thumb Infant Resuscitator (TG).
0310093	Tom Thumb Infant Resuscitator.
0310094	Tom Thumb Infant Resuscitator.
0310096	Tom Thumb Oxygen Vent with Pole Mount
0310097	Tom Thumb Infant Resuscitator (TG).
0310098	Tom Thumb Oxygen Flow.
0310099	Tom Thumb Oxygen Vent.
0310250	Tom Thumb Conversion Kit.
0320102	Tom Thumb Rectangular Mounting Block.
0320103	Tom Thumb Rectangular Mounting Block.
0330185	Wing Nut Screw Black 25mm Thread
0330190	Universal clamp - 44mm horizontal
0330191	Clamp Spacer
0330193	Pressure gauge: 63mm dial -5 to +60 cmWG
0330200	TT460 MAIN Body Block
0330201	TT480/490 MAIN Body Block
0330202	TT495 MAIN Body Block
0330203	Pressure gauge for `Tom Thumb` Infant
0330204	Flowmeter (air) 0-15 Lpm
0330205	Flowmeter Meter Bolt
0330206	Blocking Bolt And Nut
0330207	Blanking Bolt
0330208	Inlet 10 mm
0330209	Outlet 15mm
0330210	Precision Valve
0330211	Adjustable Valve
0330212	Gauge Seal
0330213	O Ring Thick
0330214	O Ring Thin
0330215	Clamp Spacer
0330216	Screws For Clamp M5
0330217	Oxygen hose 1m. C/W indirect probe &
0330218	Oxygen Hose 3m. C/W indirect probe &
0330219	Right Angle For Hose
0330220	Fomblin OT20N Grease 100g tube
0330221	Lock Tight
0330222	Precision Valve Spring
0330223	Adjustable Valve Spring
0330224	Precision Valve Body
0330225	Precision/Adjustable Valve Seating
0330227	Precision/Adjustable Valve Seat
0330228	Precision Valve Adjustment Screw
0330229	Adjustable Valve Body

0330230	Adjustable Valve Screw
0330231	Adjustable Valve Internal Stop
0330232	Adjustable Valve Collar
0330233	Adjustable Valve Centre Screw
0330234	Adjustable Valve Side Screw
0330235	Inlet Restrictor
0330236	Replacement Gauge Face.
0330241	TT480/490 MAIN Body Block
0330242	Medirail Mounting Bracket
0330243	Adjustable Valve Side Screw. Bag 100
0330246	Blocking Bolt And Nut
0330247	Blanking Bolt
0330248	Inlet 10 mm
0330249	Outlet 15mm
0330250	Oxygen hose 1.5m. C/W DIN fitting:
0330251	Oxygen hose 3.0m. C/W DIN fitting:
0330263	Oxygen Hose 2m. C/W indirect probe &
0330268	Air Hose 3m. C/W indirect probe &
0330291	Inlet 10mm Thumbelina Inlet
0330292	Outlet 15mm Thumbelina
0330295	Outlet 10mm do not use
0330296	Blanking Bolt
0330297	Clamp Spacer - dove tailed.
0330300	Key Ring chain 1 inch ball
0330800	Label - pressure control for Tom Thumb
0340000	Labour - Tom Thumb infant resuscitator.
0340001	Training - Tom Thumb infant resuscitator
0340002	Training - Tom Thumb infant resuscitator
0350120	Packing box for Tom Thumb TT480.
0350121	Packing box for Tom Thumb TT490 Series.
0380000	Annual service - Tom Thumb
0380030	Tom Thumb Infant Resuscitator.
0390003	Tom Thumb Instruction Manual
0390004	Leaflet - Tom Thumb Infant Resuscitator.
0390005	Tom Thumb Training Manual
0390006	Tom Thumb Instruction Manual
0390007	Tom Thumb Instruction Manual
0390008	Tom Thumb Serial Number Label
0390009	`Tom Thumb` Label
0390011	Tom Thumb CE Label
0390015	Tom Thumb label
0390016	Address label for Tom Thumb gauge.
0390017	Gauge label - Tom Thumb
0390020	Tom Thumb Service Manual
0390021	Tom Thumb Service Manual TT460
0390022	Tom Thumb Service Manual TT480
0390023	Tom Thumb Service Manual TT490-15
0390024	Tom Thumb Service Manual TT496
0390060	Tom Thumb Instruction Manual - TT495

Comments on Stock References

part numbers upto date

Further Action Required ☐

Supplier Information Main and Component Parts

Stock_ Reference	Description	Supplier AC	Supplier PN	Supplier Name	Rating	CRM	Reviewed
0330210	Precision Valve	00009117	9901-4000009117	JNK Tooling	C		No longer Supplier (2015)
0330211	Adjustable Valve	00009117	9901-1800009117	JNK Tooling	C		No longer Supplier (2015)
0330214	O Ring Thin	00009123	00009123	RS Components Ltd	B		No longer Supplier (2020)
0390004	Leaflet - Tom Thumb Infant Resuscitator.	00009360	00009360	Print Centre	C		No longer Supplier (2017)
0330201	TT480/490 MAIN Body Block	00010160	033020100010160	KSJ Engineering Limited	B		No longer Supplier (2016)
0330216	Screws For Clamp M5	00009123	483-961600009123	RS Components Ltd	B		No longer Supplier (2020)

0330030	Dovetail (male)	00009150	450500009150	Therapy Equipment Ltd	B		No longer Supplier (2019)
0330203	Pressure gauge for `Tom Thumb` Infant	00010395	00010395	Wika Instruments Ltd	B		No longer Supplier (2020)
0330205	Flowmeter Meter Bolt	00010160	033020500010160	KSJ Engineering Limited	B		No longer Supplier (2016)
0330209	Outlet 15mm	00010160	00010160	KSJ Engineering Limited	B		No longer Supplier (2016)
0330215	Clamp Spacer	00010160	9901-2200010160	KSJ Engineering Limited	B		No longer Supplier (2016)
0330212	Gauge Seal	00009150	00009150	Therapy Equipment Ltd	B		No longer Supplier (2019)
0330218	Oxygen Hose 3m. C/W indirect probe &	00009150	00009150	Therapy Equipment Ltd	B		No longer Supplier (2019)
0320060	Flowmeter (Oxygen) 0-15 Lpm.	00009150	650100009150	Therapy Equipment Ltd	B		No longer Supplier (2019)
0330220	Fomblin OT20N Grease 100g tube	00009123	00009123	RS Components Ltd	B		No longer Supplier (2020)
0330221	Lock Tight	00009123	514-51500009123	RS Components Ltd	B		No longer Supplier (2020)
0330219	Right Angle For Hose	00010770	202000010770	Chelmer Pneumatics Ltd	B		No longer Supplier (2010)
0330201	TT480/490 MAIN Body Block	00010725	033020100010725	MAJ Hi-Spec Ltd	B		No longer Supplier (2016)
0330205	Flowmeter Meter Bolt	00010725	033020500010725	MAJ Hi-Spec Ltd	B		No longer Supplier (2016)
0330209	Outlet 15mm	00010725	00010725	MAJ Hi-Spec Ltd	B		No longer Supplier (2016)
0330234	Adjustable Valve Side Screw	00009117	033023400009117	JNK Tooling	C		No longer Supplier (2015)
0330235	Inlet Restrictor	00009117	033023500009117	JNK Tooling	C		No longer Supplier (2015)
0330200	TT460 MAIN Body Block	00010160	033020000010160	KSJ Engineering Limited	B		No longer Supplier (2016)
0330202	TT495 MAIN Body Block	00010160	033020200010160	KSJ Engineering Limited	B		No longer Supplier (2016)
0330204	Flowmeter (air) 0-15 Lpm	00009150	650200009150	Therapy Equipment Ltd	B		No longer Supplier (2019)
0330200	TT460 MAIN Body Block	00010725	033020000010725	MAJ Hi-Spec Ltd	B		No longer Supplier (2016)
0330202	TT495 MAIN Body Block	00010725	033020200010725	MAJ Hi-Spec Ltd	B		No longer Supplier (2016)
0330055	Universal clamp - 44mm horizontal	00009150	200200009150	Therapy Equipment Ltd	B		No longer Supplier (2019)
0330208	Inlet 10 mm	00010160	00010160	KSJ Engineering Limited	B		No longer Supplier (2016)
0330236	Replacement Gauge Face.	00010395	139168200010395	Wika Instruments Ltd	B		No longer Supplier (2020)
0330207	Blanking Bolt	00009123	033020700009123	RS Components Ltd	B		No longer Supplier (2020)
0330208	Inlet 10 mm	00010725	00010725	MAJ Hi-Spec Ltd	B		No longer Supplier (2016)
0330800	Label - pressure control for Tom Thumb	00009930	00009930	Renlim Labels	B		No longer Supplier (2019)
0330222	Precision Valve Spring	00012028	C640900012028	Harlow Springs Limited	C		No longer Supplier (2018)
0330223	Adjustable Valve Spring	00012028	996800012028	Harlow Springs Limited	C		No longer Supplier (2018)
0330217	Oxygen hose 1m. C/W indirect probe &	00009150	00009150	Therapy Equipment Ltd	B		No longer Supplier (2019)
0330213	O Ring Thick	00009123	00009123	RS Components Ltd	B		No longer Supplier (2020)
0330207	Blanking Bolt	00010725	033020700010725	MAJ Hi-Spec Ltd	B		No longer Supplier (2016)
0330206	Blocking Bolt And Nut	00010160	00010160	KSJ Engineering Limited	B		No longer Supplier (2016)
0390015	Tom Thumb label	00010665	TT IR00010665	Pendle Signs & Plastics Ltd.	B		No longer Supplier (2018)
0390016	Address label for Tom Thumb gauge.	00010665	TT ADR00010665	Pendle Signs & Plastics Ltd.	B		No longer Supplier (2018)
0390017	Gauge label - Tom Thumb	00010665	TT VIA00010665	Pendle Signs & Plastics Ltd.	B		No longer Supplier (2018)
0330216	Screws For Clamp M5	00009061	483-961600009061	Farnell	B		No longer Supplier (2020)

0330224	Precision Valve Body	00009117	00009117	JNK Tooling	C		No longer Supplier (2015)
0330228	Precision Valve Adjustment Screw	00009117	00009117	JNK Tooling	C		No longer Supplier (2015)
0330227	Precision/Adjustable Valve Seat	00009117	00009117	JNK Tooling	C		No longer Supplier (2015)
0330295	Outlet 10mm do not use	00010160	033029500010160	KSJ Engineering Limited	B		No longer Supplier (2016)
0330300	Key Ring chain 1 inch ball	00012395	BC600012395	Seton	B		No longer Supplier (2010)
3210071	Single use silicone round facemask	00009721	00009721	Marshall Products Limited	B		No longer Supplier (2012)
3210011	NeoPEEP neonatal single patient use	00009721	00009721	Marshall Products Limited	B		No longer Supplier (2012)
0330250	Oxygen hose 1.5m. C/W DIN fitting:	00009150	00009150	Therapy Equipment Ltd	B		No longer Supplier (2019)
0330251	Oxygen hose 3.0m. C/W DIN fitting:	00009150	00009150	Therapy Equipment Ltd	B		No longer Supplier (2019)
0320103	Tom Thumb Rectangular Mounting Block.	00009022	PP869400009022	Tanda Engineering Ltd	B		No longer Supplier (2020)
0330206	Blocking Bolt And Nut	00010725	00010725	MAJ Hi-Spec Ltd	B		No longer Supplier (2016)
3210071	Single use silicone round facemask	00009602	00009602	Bunzl Healthcare	B		No longer Supplier (2020)
3210011	NeoPEEP neonatal single patient use	00009602	00009602	Bunzl Healthcare	B		No longer Supplier (2020)
0330241	TT480/490 MAIN Body Block	00010160	9901-0100010160	KSJ Engineering Limited	B		No longer Supplier (2016)
0330246	Blocking Bolt And Nut	00010160	9901-0600010160	KSJ Engineering Limited	B		No longer Supplier (2016)
0330248	Inlet 10 mm	00010160	4701-6100010160	KSJ Engineering Limited	B		No longer Supplier (2016)
0330249	Outlet 15mm	00010160	9901-3400010160	KSJ Engineering Limited	B		No longer Supplier (2016)
0350120	Packing box for Tom Thumb TT480.	00009353	RPT10B00009353	Rajapack Ltd	B		No longer Supplier (2019)
0350121	Packing box for Tom Thumb TT490 Series.	00009353	RPT17B00009353	Rajapack Ltd	B		No longer Supplier (2019)
0390006	Tom Thumb Instruction Manual	00009345	00009345	Duffield Printers	A		No longer Supplier (2014)
0390007	Tom Thumb Instruction Manual	00009345	00009345	Duffield Printers	A		No longer Supplier (2014)
0330234	Adjustable Valve Side Screw	00009073	033023400009073	Modern Screws	B		No longer Supplier (2014)
0330268	Air Hose 3m. C/W indirect probe &	00009150	ZZ HOSE 3M00009150	Therapy Equipment Ltd	B		No longer Supplier (2019)
0330243	Adjustable Valve Side Screw. Bag 100	00009073	02086507020602000009073	Modern Screws	B		No longer Supplier (2014)
0330296	Blanking Bolt	00009123	367-609800009123	RS Components Ltd	B		No longer Supplier (2020)
3210011	NeoPEEP neonatal single patient use	FL39B01	FL39B01		Needs Grading		No longer Supplier (1970)
0320221	Pole Clamp without slide bracket.	00007317	IHCPC1-SPECIAL00007317	Inspiration Healthcare Limited	B		No longer Supplier (2020)
0330190	Universal clamp - 44mm horizontal	00009032	00009032	R & L Enterprises Ltd	B		No longer Supplier (2019)
0330212	Gauge Seal				B		No longer Supplier (2020)
0390004	Leaflet - Tom Thumb Infant Resuscitator.	00009363	00009363	ACI Print Group Impress Print	B		No longer Supplier (2016)
0330209	Outlet 15mm	00009032	00009032	R & L Enterprises Ltd	B		No longer Supplier (2019)
0330208	Inlet 10 mm	00009032	00009032	R & L Enterprises Ltd	B		No longer Supplier (2019)
0330206	Blocking Bolt And Nut	00009032	00009032	R & L Enterprises Ltd	B		No longer Supplier (2019)
0330224	Precision Valve Body	00009032	00009032	R & L Enterprises Ltd	B		No longer Supplier (2019)
0330228	Precision Valve Adjustment Screw	00009032	00009032	R & L Enterprises Ltd	B		No longer Supplier (2019)
0330225	Precision/Adjustable Valve Seating	00009032	00009032	R & L Enterprises Ltd	B		No longer Supplier (2019)
0330227	Precision/Adjustable Valve Seat	00009032	00009032	R & L Enterprises Ltd	B		No longer Supplier (2019)

0330229	Adjustable Valve Body	00009032	00009032	R & L Enterprises Ltd	B		No longer Supplier (2019)
0330234	Adjustable Valve Side Screw	00009032	033023400009032	R & L Enterprises Ltd	B		No longer Supplier (2019)
0330233	Adjustable Valve Centre Screw	00009032	00009032	R & L Enterprises Ltd	B		No longer Supplier (2019)
0330205	Flowmeter Meter Bolt	00009032	033020500009032	R & L Enterprises Ltd	B		No longer Supplier (2019)
0330230	Adjustable Valve Screw	00009032	00009032	R & L Enterprises Ltd	B		No longer Supplier (2019)
0330232	Adjustable Valve Collar	00009032	00009032	R & L Enterprises Ltd	B		No longer Supplier (2019)
0330231	Adjustable Valve Internal Stop	00009032	00009032	R & L Enterprises Ltd	B		No longer Supplier (2019)
0330201	TT480/490 MAIN Body Block	00009032	033020100009032	R & L Enterprises Ltd	B		No longer Supplier (2019)
0330263	Oxygen Hose 2m. C/W indirect probe &	00009150	ZZ HOSE 2M00009150	Therapy Equipment Ltd	B		No longer Supplier (2019)
0330193	Pressure gauge: 63mm dial -5 to +60 cmWG	00009097	7465600009097	GNW - Gauges North West Ins	D		No longer Supplier (2017)
0330191	Clamp Spacer	00009032	A491700009032	R & L Enterprises Ltd	B		No longer Supplier (2019)
0390006	Tom Thumb Instruction Manual	00009366	00009366	Warwick Printing	B		No longer Supplier (2020)
0330242	Medirail Mounting Bracket	00009032	A509000009032	R & L Enterprises Ltd	B		No longer Supplier (2019)
0390003	Tom Thumb Instruction Manual	00009366	00009366	Warwick Printing	B		No longer Supplier (2020)
0330230	Adjustable Valve Screw	00009117	00009117	JNK Tooling	C		No longer Supplier (2015)
0330185	Wing Nut Screw Black 25mm Thread	00009032	333500009032	R & L Enterprises Ltd	B		No longer Supplier (2019)
0330217	Oxygen hose 1m. C/W indirect probe &				B		No longer Supplier (2020)
0330203	Pressure gauge for `Tom Thumb` Infant				B		No longer Supplier (2020)
3210071	Single use silicone round facemask	FL39B01	FL39B01		Needs Grading		No longer Supplier (1970)
0330218	Oxygen Hose 3m. C/W indirect probe &				B		No longer Supplier (2020)
0330220	Fomblin OT20N Grease 100g tube				B		No longer Supplier (2020)
0330218	Oxygen Hose 3m. C/W indirect probe &	00009211	00009211	MEC Medical Limited	B		No longer Supplier (2020)

Comments on Suppliers

no issues with suppliers

Further Action Required ☐

Sales Information

Stock Ref		2001	2002	2003	2004	2005	2006	2007	2008	2009	2010	2011	2012	2013	2014	2015	2016	2017	2018	2019
0310027	Tom Thumb Infant Resuscitator (Air).																			
0310029	Tom Thumb Infant Resuscitator (Air).																			
0310030	Tom Thumb Infant Resuscitator.	12	6	3	31	44	29	33	102	58	80	121	259	94	129	67	48	52	58	10
0310031	Tom Thumb Infant Resuscitator (Oxygen).	3																		
0310032	Tom Thumb Infant Resuscitator (Oxygen).	11	2	26	6	9	5	2			6			3						
0310033	Tom Thumb Infant Resuscitator (Oxygen).	3																		
0310034	Tom Thumb Infant Resuscitator (Oxygen).	5	16	4	1	6	1	1		4		1		1	2	1				
0310072	Tom Thumb Infant Resuscitator.								11											

0310074	Tom Thumb Infant Resuscitator.																		
0310075	Tom Thumb Infant Resuscitator.																		
0310076	Tom Thumb Infant Resuscitator.																		
0310080	Tom Thumb Infant Resuscitator (RO).					12							2						
0310082	Tom Thumb Infant Resuscitator (RO).																		
0310092	Tom Thumb Infant Resuscitator (TG).								2										
0310093	Tom Thumb Infant Resuscitator.											2							
0310094	Tom Thumb Infant Resuscitator.										4		9		6	1	1		
0310096	Tom Thumb Oxygen Vent with Pole Mount												1						
0310097	Tom Thumb Infant Resuscitator (TG).												5	7	18	16	1	6	
0310098	Tom Thumb Oxygen Flow.																		
0310099	Tom Thumb Oxygen Vent.					1			5				1	1					
0310250	Tom Thumb Conversion Kit.											6	6	28	3	33		1	1
0320102	Tom Thumb Rectangular Mounting Block.																		
0320103	Tom Thumb Rectangular Mounting Block.								1										
0330185	Wing Nut Screw Black 25mm Thread																		
0330190	Universal clamp - 44mm horizontal												2	100	50		43		10
0330191	Clamp Spacer												1				1	1	
0330193	Pressure gauge: 63mm dial -5 to +60 cmWG														18	22	17	13	2
0330200	TT460 MAIN Body Block																		
0330201	TT480/490 MAIN Body Block				11	2	5	1				1					1		
0330202	TT495 MAIN Body Block												1						
0330203	Pressure gauge for 'Tom Thumb' Infant	6	10	16	16	18	17	19	17	30	20	13	10	21	9		1		24
0330204	Flowmeter (air) 0-15 Lpm																		
0330205	Flowmeter Meter Bolt																		
0330206	Blocking Bolt And Nut											16	12	10				2	
0330207	Blanking Bolt	1																	
0330208	Inlet 10 mm			2	1							46	12	10		5		1	
0330209	Outlet 15mm	0						5								15	7	5	15
0330210	Precision Valve	2	2	-1		2	13	1		2	2	4		3	2	3	2	8	4
0330211	Adjustable Valve	1	1	1	3	-3	5	4	1		2	1	2	3	3	5	7	7	13
0330212	Gauge Seal	6		3	31	16	13	19	34	29	55	31	40	28	28	69	55	81	69
0330213	O Ring Thick	10	20	24	120	26	26	52	12	40	13	77	8	34	193	90	95	92	77
0330214	O Ring Thin	25	20	24	160	98	125	142	151	235	135	358	48	227	137	534	541	674	97
0330215	Clamp Spacer			2	2	10			4	1									
0330216	Screws For Clamp M5					20			28		4							2	
0330217	Oxygen hose 1m. C/W indirect probe &			1	1	1	1			8	2			6	6	2	4	16	5
0330218	Oxygen Hose 3m. C/W indirect probe &			1	10	3		12	6	28	10	15	2	7	17	12		3	7
0330219	Right Angle For Hose					3					1	2							10

0330220	Fomblin OT20N Grease 100g tube														1				1	
0330221	Lock Tight																			
0330222	Precision Valve Spring											1								
0330223	Adjustable Valve Spring				5		5					2		1						
0330224	Precision Valve Body																			
0330225	Precision/Adjustable Valve Seating						5													
0330227	Precision/Adjustable Valve Seat				10		4			1		2	1		2	3		9		
0330228	Precision Valve Adjustment Screw						4													
0330229	Adjustable Valve Body						1						1			1				
0330230	Adjustable Valve Screw														4					
0330231	Adjustable Valve Internal Stop																			
0330232	Adjustable Valve Collar																			
0330233	Adjustable Valve Centre Screw																			
0330234	Adjustable Valve Side Screw													3						
0330235	Inlet Restrictor		3				1	4			5					3				
0330236	Replacement Gauge Face.		5	18	2	14	2	20	5	9	17	6	5	13	13	24	8	19	3	6
0330241	TT480/490 MAIN Body Block																			
0330242	Medirail Mounting Bracket																			
0330243	Adjustable Valve Side Screw. Bag 100																			
0330246	Blocking Bolt And Nut																			
0330247	Blanking Bolt																			
0330248	Inlet 10 mm																			
0330249	Outlet 15mm																			
0330250	Oxygen hose 1.5m. C/W DIN fitting:																			
0330251	Oxygen hose 3.0m. C/W DIN fitting:																			
0330263	Oxygen Hose 2m. C/W indirect probe &												3	12	2		2	5	2	
0330268	Air Hose 3m. C/W indirect probe &																			
0330291	Inlet 10mm Thumbelina Inlet																			
0330292	Outlet 15mm Thumbelina																			
0330295	Outlet 10mm do not use									2						8				
0330296	Blanking Bolt																4			
0330297	Clamp Spacer - dove tailed.										5									
0330300	Key Ring chain 1 inch ball																			
0330800	Label - pressure control for Tom Thumb																			
0340000	Labour - Tom Thumb infant resuscitator.																			
0340001	Training - Tom Thumb infant resuscitator																			
0340002	Training - Tom Thumb infant resuscitator																			

0350120	Packing box for Tom Thumb TT480.																			
0350121	Packing box for Tom Thumb TT490 Series.																			
0380000	Annual service - Tom Thumb																			
0380030	Tom Thumb Infant Resuscitator.																			
0390003	Tom Thumb Instruction Manual																			
0390004	Leaflet - Tom Thumb Infant Resuscitator.									25		6			50		10			
0390005	Tom Thumb Training Manual																			
0390006	Tom Thumb Instruction Manual						49				6	3			20					
0390007	Tom Thumb Instruction Manual										6	3								
0390008	Tom Thumb Serial Number Label																			
0390009	Tom Thumb Label																			
0390011	Tom Thumb CE Label																			
0390015	Tom Thumb label																5		2	
0390016	Address label for Tom Thumb gauge.																			
0390017	Gauge label - Tom Thumb																			
0390020	Tom Thumb Service Manual																			
0390021	Tom Thumb Service Manual TT460																			
0390022	Tom Thumb Service Manual TT480																			
0390023	Tom Thumb Service Manual TT490-15																			
0390024	Tom Thumb Service Manual TT496																			
0390060	Tom Thumb Instruction Manual - TT495																			
		85	85	124	410	282	305	316	379	476	370	710	411	541	711	952	849	1029	1046	14

Comments on Sales Information

just spare parts going out, no issues

Further Action Required ☐

New Countrys

Country	2001	2002	2003	2005	2004	2006	2008	2011	2007	2009	2010	2013	2012	2014	2015	2016	2018	2017	2019	2020	2023	2
Belgium																						
Canada																						
China																						
Czech Republic																						
Finland																						
France																						
Germany																						
IN India																						
Italy																						
Nigeria																						
Portugal																						
Russia																						
Saudi Arabia																						
Scotland																						
Sri Lanka																						
Switzerland																						

Turkey																						
UK																						
UK United Kingdom																						
United Kingdom																						

Comments on Sales Countrys

no issues with countries

Further Action Required ☐

Sales Countrys - New Risks?

Further Action Required ☐

Returns Information

Info	Stock Ref	Fault	2006	2007	2008	2009	2010	2011	2012	2013	2014	2015	2016	2017	2018	2019	2020	2021	2022	2023	2024
	0310030	Needs Evaluating												1							
	0310030	Adjusted Factory Set Valves							3			1	2	1	2						
	0310030	Calibrated Gauge									1										
	0310030	Converted To TT480		1				1							1	2					
	0310030	Converted to TT480 and service																	4	5	
	0310030	Customer Declined Quotation										2									
	0310030	Damaged Pressure Gauge		1	2	4	2			1	1		2	1	1						
	0310030	N/a	1																		
	0310030	New							1		2								2		
	0310030	No Fault Found			1							2		2	2						
	0310030	No Fault Return To Stock			1																
	0310030	Renovated															3	13	3		
	0310030	Repaired Clamp And Serviced										1									
	0310030	Replace Gauge And Service								1	1	4	3	2	2	11	2	4	6	2	
	0310030	Replaced Adjustable Valve													1	1					
	0310030	Replaced Gauge													2	1		1	1		
	0310030	Replaced Gauge and Adjustable Valve													2						
	0310030	Replaced gauge face													1	1					

	0310030	Replaced Gauge Face And Service											1			1	2	1	3	1	1
	0310030	Replaced precision valve and serviced																	1		
	0310030	Serviced	1		2	3							1	10	14	3	4		15	23	3
	0310030	Serviced / Adjusted Factory Set Valves				2	1		3	11	11	10	12	3	7	15	5	10	2	2	22
	0310030	Serviced and Upgraded															1				
	0310030	Stuck Gauge												1	1						
	0310030	Tested				1						1	1	1	1					3	
	0310031	New																		1	
	0310031	Serviced											2								
	0310032	Adjusted Factory Set Valves					1														
	0310032	Calibrated Gauge				1															
	0310032	Converted To TT480						1													
	0310032	Damaged Pressure Gauge					1														
	0310032	New			39																
	0310032	Replaced Gauge Face And Service												2		1					
	0310032	Serviced	2		1									1	3	1					
	0310032	Serviced / Adjusted Factory Set Valves				1	2			1				1							
	0310032	Serviced and Upgraded												4							
	0310032	Tested				1				1											
	0310033	Replace Gauge And Service										1									
	0310033	Serviced		1																	
	0310034	Adjusted Factory Set Valves				2															
	0310034	Converted To TT480										1									
	0310034	Converted to TT480 and service																			1
	0310034	Damaged Pressure Gauge			1																
	0310034	New																			1
	0310034	Replace Gauge And Service														1					

	0310034	Replaced hose and service.																			1
	0310034	Serviced	1	1					1				1			2	1				
	0310034	Serviced / Adjusted Factory Set Valves					2	1	1	2		5		2		2		2			
	0310034	Unchecked - Returned To Stock					12														
	0310076	New				1															
	0310080	Replace Gauge And Service									1										
	0310094	Converted To TT480								9											
	0310094	Not Delivering More Than 10cm												1							
	0310094	Replace Gauge And Service												1							
	0310094	Replaced Gauge Face And Service												1					1		
	0310097	Serviced												1	1		1		1	1	
	0310097	Serviced / Adjusted Factory Set Valves									1	1	1	1		1					1
	0310098	Serviced / Adjusted Factory Set Valves									3	1									
	0310099	Serviced / Adjusted Factory Set Valves									2										
	0330193	No Fault Found													1						
	0330193	Unrepaired													2						
	0330201	Needs Evaluating										10									
	0330217	Unchecked - Returned To Stock					6				8										
	0330218	Tested for leaks															4				
	0380030	Needs Evaluating											3			1					
	0380030	New															3	3			
	0380030	Replaced gauge face																	1		
	0380030	Service exchange unit													3						
	0380030	Tested														3				3	
	0390023	Serviced		1																	

Comments on Returns

nothing new
just normal servicing

Further Action Required ☐

Returns / New Re-Work Risks?

Further Action Required ☐

Design Changes / Product Changes

Showing Documents Filed in Y 14 Design Changes Red lines are New Since last review : 03/05/24

Added On	Added By	ID	Document Description	
09/04/18	John Lamb		25816	Tom Thumb Design Changes Mounting screws
12/03/18	Derek Lamb		25578	Tom Thumb Design Changes QC28 B Tom Thumb Outlet
13/11/17	John Lamb		23872	Tom Thumb Design Changes Tom Thumb Darlington Memorial
13/11/17	John Lamb		23867	Tom Thumb Design Changes BRI 2005 Tom Thumb
08/11/17	Derek Lamb		23803	Tom Thumb Product Changes Discontinuation of Products in the range
02/11/17	John Lamb		23703	BRI Documentation 2005
09/03/11	Derek Lamb		7914	Tom Thumb Drawings 0330296
07/03/11	John Lamb		7911	Tom Thumb Design Changes TT480 modification
07/03/11	John Lamb		7909	Tom Thumb Design Changes 15mm outlets to 10mm taper
07/03/11	John Lamb		7905	Tom Thumb Design Changes Dual Outlet 15mm to 10mm Taper
14/02/11	Derek Lamb		7658	0310030 TT480 labelling modification 140211
19/01/11	john lamb		7452	Tom Thumb Product Changes G&R Hydroparts Block rounded corners
19/01/11	john lamb		7449	Tom Thumb Product Changes 2005
31/01/08	Derek Lamb		2995	Tom Thumb Design Changes 2005 dual outlet taper_G & R
31/01/08	Derek Lamb		2994	Tom Thumb Design Changes 2005 Design Change Rationale
31/01/08	Derek Lamb		2974	Tom Thumb Product Changes 2004
19/06/07	Derek Lamb		2451	Tom Thumb Design Changes QC 28
19/06/07	Derek Lamb		2221	Tom Thumb Product Changes Product changes 2000 Restrictor
19/06/07	Derek Lamb		2220	Tom Thumb Product Changes Product Changes 1998 Flow meter
19/06/07	Derek Lamb		2219	Tom Thumb Product Changes

Comments on Design Changes / Product Changes

no design changes

Further Action Required ☐

Risks Associated with Design Changes / Product Changes or
Do any Notified bodies MHRA require notification of any changes

Further Action Required / Yes requires notification ☐

Current User / Instruction

Showing Documents Filed in F 5 User Instructions for use Red lines are New Since last review : 03/05/24

Added On	Added By	ID	Document Description
01/09/15	Derek Lamb	15572	Tom Thumb Instructions for Use / User Manual TT490-15
01/09/15	Derek Lamb	15571	Tom Thumb Instructions for Use User Manual TT480

01/09/15	Derek Lamb	15570	Tom Thumb Instructions for Use User Manual 0310027
24/08/15	Derek Lamb	15511	Tom Thumb Instructions for Use User Manual TT480 Instructions Dutch High Quality 5mm bleed
24/08/15	Derek Lamb	15510	Tom Thumb Instructions for Use User Manual TT480 Instructions Dutch High Quality
24/08/15	Derek Lamb	15509	Tom Thumb Instructions for Use User Manual TT480 Instructions Dutch Email
14/10/13	Derek Lamb	13047	Tom Thumb Instructions for Use User Manual French
30/11/12	Derek Lamb	12211	Tom Thumb Instructions for Use User Manual TT495-15
23/03/12	John Lamb	11732	Tom Thumb Instructions for Use User Manual
06/07/11	John Lamb	8512	Tom Thumb Instructions for Use User Manual NeoPEEP
06/07/11	John Lamb	8506	Instructions For Use Tom Thumb TT495
23/02/09	Derek Lamb	5107	Instructions For Use Tom Thumb TT480

Comments on Instruction Manuals- Still fit for Purpose

no instruction updes

Further Action Required ☐

Risks Associated with Instruction Manuals

Further Action Required ☐

Current Labels

Showing Documents Filed in F 7 / F 8 Labels Red lines are New Since last review : 03/05/24

Added On	Added By	ID	Document Description
07/12/18	Derek Lamb	28357	Tom Thumb Labels Original labels as at Dec 2018
17/07/17	Derek Lamb	20970	Tom Thumb Labels Top Box Sticker example
19/10/16	Catrin Hird	17722	0330193 Pressure gauge: 63mm dial -5 to +60 cmWG (1416 / 1437)
27/03/12	Derek Lamb	11748	Tom Thumb Labels 0390016 Address Sticker Round outside Gauge
27/03/12	Derek Lamb	11747	Tom Thumb Labels Caution Label 0390018 Top Block
27/03/12	Derek Lamb	11746	Tom Thumb Labels 0390017 Front of Gauge Sticker
27/03/12	Derek Lamb	11745	Tom Thumb Labels 0390015 Front Block Logo Sticker
27/03/12	Derek Lamb	11744	Tom Thumb Labels 0390010 Adjustable Valve Sticker
27/03/12	Derek Lamb	11743	Tom Thumb Labels Part number to Label
27/03/12	Derek Lamb	11742	Tom Thumb Labels TT491 Serial
27/03/12	Derek Lamb	11741	Tom Thumb Labels TT490 Serial
27/03/12	Derek Lamb	11740	Tom Thumb Labels Serial TT480
14/02/11	Derek Lamb	7658	0310030 TT480 labelling modification 140211
14/05/09	Derek Lamb	5922	Tom Thumb Labels mediprema labels
19/06/07	Derek Lamb	2254	Tom Thumb Parts List Adj label 04
19/06/07	Derek Lamb	2204	Tom Thumb Labels Tomb Thumb Infant Resuscitator

Comments on Labels- Still fit for Purpose

no label changes

Further Action Required ☐

Risks Associated with Labels

Further Action Required ☐

All other Documentation

Showing New / Updated Documents Red lines are New Since last review : 03/05/24

Added On	Added By	ID	Document Description
03/08/23	Derek Lamb	125906	VM3COP50.19 TT480 service manual
13/04/23	Derek Lamb	116149	Tom Thumb Post Market Surveillance Periodic Safety Update Report
25/10/21	Steve Nixon	73112	WIKA ISO 9001:2015 Certificate

15/10/21	Steve Hardaker	72371	Tom Thumb Material Safety Data Sheet 2004 MH052_MSDS_20050405
30/09/19	Derek Lamb	31154	VM3COP50.13 Quality Control Tom Thumb
30/09/19	Derek Lamb	31150	VM3COP50.01 Tom Thumb TT490 (0310031 / 0310033)
30/09/19	Derek Lamb	31146	VM3COP50.02 Tom Thumb production procedures TT480 (0310030)
30/09/19	Derek Lamb	31142	VM3COP50.03.2 Tom Thumb TT490_15_Restricted
30/09/19	Derek Lamb	31138	VM3COP50.04 Tom Thumb TT495
30/09/19	Derek Lamb	31134	VM3COP50.05 Tom Thumb Adjustable Valve
30/09/19	Derek Lamb	31130	VM3COP50.06 Tom Thumb Precision Valve
30/09/19	Derek Lamb	31126	VM3COP50.09 Tom Thumb TT490 - TT480 Conversion Procedure
30/09/19	Derek Lamb	31122	VM3COP50.16 TT480 Top Gauge 0310097
30/09/19	Derek Lamb	31118	VM3COP50.18 TT Valve Setup and Service
30/09/19	Derek Lamb	31001	Tom Thumb Risk analysis reports Risk assessment
30/09/19	Derek Lamb	30993	Tom Thumb Risk analysis reports Risks of Assembling or Servicing a Tom Thumb
29/05/19	Derek Lamb	29753	Viamed CE Certificate CE01389 Internal Use Only Cancellation Letter
04/04/19	Steve Hardaker	29246	Tom Thumb Data Sheets - MSDS - 7114898 Klinger PTFE tape_04_04_19
13/12/18	Derek Lamb	28381	Tom Thumb Parts List Relief valve adjustable inc part numbers
13/12/18	Derek Lamb	28379	Tom Thumb Parts List Relief valve inc part numbers
07/12/18	Derek Lamb	28357	Tom Thumb Labels Original labels as at Dec 2018
20/09/18	Derek Lamb	27474	VM3COP02.02 Viamed Company Responsibility's organisation chart structure
04/09/18	Derek Lamb	27256	Tom Thumb VM3COP Procedures Rework of Tom Thumbs
29/08/18	Derek Lamb	27092	Tom Thumb Literature reviews Resuscitation Guidelines 2000
20/08/18	Derek Lamb	27040	Tom Thumb Service Manual Parts required for Annual Service of TT480 and TT490 Tom Thumb v1.2
27/07/18	Derek Lamb	26802	Tom Thumb Training Information Pressure Pipes World Wide
19/07/18	John Lamb	26720	Tom Thumb Literature reviews Tom Thumb Papers
15/06/18	John Lamb	26487	Tom Thumb Manufacturing route
04/06/18	John Lamb	26301	Tom Thumb Legal Actions and Litigation reports Shrewsbury
30/05/18	Derek Lamb	26245	Tom Thumb QA Procedures QA Records
23/04/18	John Lamb	25973	Tom Thumb Leakage and Liquid damage ingress
09/04/18	John Lamb	25814	Tom Thumb Description of Device Variants discontinued
06/04/18	Derek Lamb	25770	Tom Thumb Promotional Leaflets - Master Art work Gauge Wika
06/04/18	John Lamb	25768	Tom Thumb Classification rationale
05/04/18	John Lamb	25758	Tom Thumb Leakage and Liquid damage Ingress Tom Thumb
03/04/18	John Lamb	25744	Tom Thumb Legal Actions and Litigation reports Sunderland Royal Inf.
03/04/18	John Lamb	25742	Tom Thumb Legal Actions and Litigation reports
13/03/18	Derek Lamb	25596	Tom Thumb Test and Design Reports GNW Gauge Fault Reports
06/03/18	Derek Lamb	25530	Tom Thumb Manufacturing route Current Suppliers 2017-2018
20/02/18	Derek Lamb	25356	VM3COP50.14 Calibration / Test Sheet Blank Tom Thumbs
13/02/18	John Lamb	25282	Tom Thumb Description of Accessories (8)
04/01/18	Derek Lamb	24820	Tom Thumb Tooling List Tom Thumb valve collar insertion tool - 0360011 Tom Thumb Drawings Tom Thumb valve collar insertion tool - 0360011
18/12/17	John Lamb	24672	Tom Thumb Flammability and Fire Risks
11/12/17	John Lamb	24623	Tom Thumb Product usability Misuse
27/11/17	Derek Lamb	24197	Tom Thumb Clinical Trials Reports Reviews and Post Market Surveillance Tom Thumb Probability
23/11/17	John Lamb	24149	Tom Thumb Packaging Trials and validation Tom Thumb Background Notes for Tom Thumb
22/11/17	Derek Lamb	24117	Tom Thumb Leakage and Liquid damage Ingress statistics
22/11/17	Derek Lamb	24101	Tom Thumb Leakage and Liquid damage TOM Thumb Ingress Issue 108582
14/11/17	John Lamb	23922	Tom Thumb Clinical Trials Reports Reviews and Post Market Surveillance Clinical Evaluation format Tom Thumb
14/11/17	John Lamb	23912	Tom Thumb Clinical Trials Research JSL Tom Thumb Clinical Papers Tom Thumb Clinical papers
14/11/17	John Lamb	23908	Tom Thumb Residual Risks Estimation of Frequency of Use of Tom Thumb Resuscitators 20171108
14/11/17	John Lamb	23904	Tom Thumb Clinical Trials Reports Reviews and Post Market Surveillance comparison of T piece occluder
14/11/17	John Lamb	23900	Tom Thumb Pre-Clinical Trials Data Tom Thumb summation annex X 1.1b
14/11/17	John Lamb	23896	Tom Thumb Residual Risks Foreseeable Misuse of Product Tom Thumb
13/11/17	John Lamb	23878	Oxygen Tents/Hoods Tanda Design Input Device Design Theory Oxygen delivery systems for infants and children

13/11/17	John Lamb	23874	Tom Thumb Purchase specifications Wika 6111.10 gauge Tom Thumb Specifications of Materials 1997 Wika
10/11/17	Derek Lamb	23855	Tom Thumb Specifications 2017
10/11/17	Derek Lamb	23849	Tom Thumb Essential requirements
09/11/17	Derek Lamb	23831	Tom Thumb Risk analysis reports Risk analysis 2017
09/11/17	John Lamb	23811	Tom Thumb Clinical Trials Reports Reviews and Post Market Surveillance Tom Thumb Probability of Incident
09/11/17	John Lamb	23809	Tom Thumb Clinical Trials Reports Reviews and Post Market Surveillance Estimated Frequency of use Probability
08/11/17	Derek Lamb	23803	Tom Thumb Product Changes Discontinuation of Products in the range
02/11/17	John Lamb	23705	BRI Documentation 2005 10mm adapter Report
02/11/17	John Lamb	23703	BRI Documentation 2005
27/10/17	Derek Lamb	23511	Tom Thumb Specifications Wika quote gauge cmH2o
26/10/17	Derek Lamb	23423	VM3COP50.20 TT Component Inspection
19/10/17	Derek Lamb	22982	Tom Thumb Description of Accessories NeoPeep File Location
18/10/17	John Lamb	22910	Tom Thumb Product Life Background Notes for Tom Thumb Life expectancy
18/10/17	John Lamb	22906	Tom Thumb Design Input Device History Tom Thumb Design
18/10/17	John Lamb	22902	Tom Thumb Design Input Device History Background Notes for Tom Thumb Design SN
18/10/17	John Lamb	22898	Tom Thumb Description of Accessories Background Notes for Tom Thumb Accessories
18/10/17	John Lamb	22860	Tom Thumb Location of responsibilities JK. Location of Responsibilities_Tom Thumb
12/10/17	John Lamb	22531	Tom Thumb Sterilisation Tom Thumb Sterility
11/10/17	John Lamb	22519	Tom Thumb Expected Life of Products Tom Thumb life Rationale (Internal document)
10/10/17	John Lamb	22494	Tom Thumb Purchase specifications Flowmeter
10/10/17	John Lamb	22492	Tom Thumb Purchase specifications Fomblin RT15
10/10/17	John Lamb	22488	Tom Thumb Specifications of Materials Nitrile O rings
10/10/17	Derek Lamb	22480	Tom Thumb Literature reviews NewbornLifeSupportClinicalGuidelinesFor
26/09/17	Derek Lamb	22240	Tom Thumb Work Instructions and Tests GNW 21_09_2017 (4)
26/09/17	Derek Lamb	22238	Tom Thumb Work Instructions and Tests GNW 21_09_2017 (3)
26/09/17	Derek Lamb	22236	Tom Thumb Work Instructions and Tests GNW 21_09_2017 (2)
26/09/17	Derek Lamb	22234	Tom Thumb Work Instructions and Tests GNW 21_09_2017 (1)
11/09/17	Derek Lamb	21995	Tom Thumb Product Life User feedback old tom thumbs
04/09/17	Derek Lamb	21738	Tom Thumb PPE Equipment
04/09/17	Derek Lamb	21722	Tom Thumb Notified Body Reports 8348484 Final Assessment E-Report
04/09/17	Derek Lamb	21720	Tom Thumb Machinery Equipment
01/09/17	Derek Lamb	21716	Tom Thumb Design Input Device History BS EN 60601-1-6 2015 Non Compliance
01/09/17	Derek Lamb	21712	Tom Thumb Design Input Device History BS EN 62366-1 2015 non compliance
31/08/17	Derek Lamb	21694	Product Control Tom Thumb
27/08/17	Derek Lamb	21625	Tom Thumb QA Procedures Example QA Report Barcode 1009081
27/08/17	Derek Lamb	21623	Tom Thumb QA Procedures Example QA Records Tom Thumb 0310030
17/07/17	Derek Lamb	20970	Tom Thumb Labels Top Box Sticker example
17/07/17	Derek Lamb	20968	Tom Thumb Expected Life of Products Proof2
17/07/17	Derek Lamb	20966	Tom Thumb Expected Life of Products Proof1
11/07/17	Derek Lamb	20916	QC 33b Tom Thumb TT480 Service Sheet
07/04/17	Derek Lamb	19538	Tom Thumb Purchase specifications GNW Instrumentation - Tom Thumb pressure gauge Tom Thumb Product Changes gauge
06/04/17	Derek Lamb	19534	Tom Thumb EC Declaration of Conformance
29/11/16	Derek Lamb	18139	Tom Thumb Work Instructions and Tests NWG test certificates 2
29/11/16	Derek Lamb	18138	Tom Thumb Work Instructions and Tests NWG test certificates
21/11/16	Derek Lamb	18115	Tom Thumb Purchase specifications Tom Thumb parts cost comparisons 17-11-2016
19/10/16	Catrin Hird	17722	0330193 Pressure gauge: 63mm dial -5 to +60 cmWG (1416 / 1437)
06/09/16	Derek Lamb	17411	Tom Thumb Clinical Papers Annex X 1.1d Tom Thumb Post Market Clinical Evaluation Follow-up
03/09/16	Derek Lamb	17382	Tom Thumb Clinical Papers Clinical Papers Annex X 1.1d supporting docs
02/09/16	Derek Lamb	17379	Tom Thumb Data Approach Database 2008 vma9001 service
02/09/16	Derek Lamb	17378	Tom Thumb Data Approach Database VMA9001 2007 service
02/09/16	Derek Lamb	17377	Tom Thumb Data Approach Database Queen elizabeth 1997 Invoice
02/09/16	Derek Lamb	17376	Tom Thumb Data Approach Database sunderland pre ce mark 1998

02/09/16	Derek Lamb	17375	Tom Thumb Data Approach Database proof of sale March 1999
02/09/16	Derek Lamb	17374	Tom Thumb Data Approach Database vma9001 2006.
02/09/16	Derek Lamb	17373	Tom Thumb Data Approach Database TT480 Tom Thumb Test Sheets VMA9001 - Feb 2000
02/09/16	Derek Lamb	17372	Tom Thumb Data Approach Database Service record vma9001 2005
02/09/16	Derek Lamb	17371	Tom Thumb Data Approach Database 2016 vma9001 service
02/09/16	Derek Lamb	17370	Tom Thumb Data Approach Database 2010_2016 vma9001 services
02/09/16	Derek Lamb	17369	Tom Thumb Data Approach Database 2003 vma9001 service
02/09/16	Derek Lamb	17368	Tom Thumb Data Approach Database 1999 vma9001 service
02/09/16	Derek Lamb	17367	Tom Thumb Data Approach Database Sales Tom Thumbs 1998_1999
02/09/16	Derek Lamb	17366	Tom Thumb Data Approach Database 1994_TT490_Sales
26/08/16	John Lamb	17336	Tom Thumb Clinical Papers
26/08/16	Derek Lamb	17334	Resuscitation Guide Neonatal Life support 2005 Tom Thumb Literature reviews Resuscitation Council (UK) Newborn Life Support
26/08/16	Derek Lamb	17333	Tom Thumb Literature reviews Royal Berkshire GL443 Neonatal Guidelines
26/08/16	Derek Lamb	17332	Tom Thumb Literature reviews Nottingham Neonatal Service - Clinical Guidelines
26/08/16	Derek Lamb	17331	Tom Thumb Literature reviews Royal Cornwall Newborn guide lines
26/08/16	Derek Lamb	17330	Tom Thumb Literature reviews Scottish Neonatal Resuscitation Course Manual
23/08/16	Derek Lamb	17307	Tom Thumb Literature reviews T-piece or self inflating bag for positive pressure ventilation during delivery room resuscitation_ An RCT - Resuscitation
23/08/16	Derek Lamb	17306	Tom Thumb Literature reviews Comparison of the T-piece resuscitator with other neonatal manual ventilation devices: A qualitative review
26/11/15	Derek Lamb	16043	Tom Thumb Clinical Trials Reports Reviews and Post Market Surveillance Risk benefits
30/09/15	Derek Lamb	15756	Tom Thumb Residual Risks to Benefits Residual risks 2015
30/09/15	Derek Lamb	15754	Tom Thumb Clinical Trials Reseach
24/09/15	Derek Lamb	15734	Tom Thumb Clinical Trials Reports Reviews and Post Market Surveillance University of Kentucky / UK HealthCare Policy and Procedure RC04-04P
24/09/15	Derek Lamb	15733	Tom Thumb Clinical Trials Reports Reviews and Post Market Surveillance Neonatal resuscitation – a practical approach. The experience of one UK tertiary neonatal unit
24/09/15	Derek Lamb	15732	Tom Thumb Clinical Trials Reports Reviews and Post Market Surveillance T-piece_resuscitator_2
24/09/15	Derek Lamb	15731	Tom Thumb Clinical Trials Reports Reviews and Post Market Surveillance T-piece_resuscitator_1
21/09/15	Derek Lamb	15681	Service Manual Tom Thumb Valves v1.1
21/09/15	Derek Lamb	15679	Tom Thumb Email Service Manual Tom Thumb TT490-15 v2.0
16/09/15	John Lamb	15647	Tom Thumb Training Information Assessment Form Tom Thumb Training Information Training assessment V1.0 2011
16/09/15	John Lamb	15644	Tom Thumb Training Information 2009
16/09/15	John Lamb	15641	Tom Thumb Training Information Training V1.1
16/09/15	John Lamb	15639	Tom Thumb Labels Original label sizes
11/09/15	Derek Lamb	15634	Tom Thumb Description of Accessories (4)
10/09/15	John Lamb	15624	Tom Thumb Components Breakdown 9901-02 Flowmeter bolt
10/09/15	John Lamb	15622	Tom Thumb Components Breakdown 9901-01 Main Block original
10/09/15	John Lamb	15621	Tom Thumb Description of Accessories (2)
01/09/15	Derek Lamb	15572	Tom Thumb Instructions for Use / User Manual TT490-15
01/09/15	Derek Lamb	15571	Tom Thumb Instructions for Use User Manual TT480
01/09/15	Derek Lamb	15570	Tom Thumb Instructions for Use User Manual 0310027
24/08/15	Derek Lamb	15522	Tom Thumb Competitors Information and Patients Neopuff technical manual
24/08/15	Derek Lamb	15521	Tom Thumb Competitors Information and Patients Neopuff specifications
24/08/15	Derek Lamb	15520	Tom Thumb Competitors Information and Patients Neopuff solutions
24/08/15	Derek Lamb	15519	Tom Thumb Competitors Information and Patients Neopuff operating instructions
24/08/15	Derek Lamb	15518	Tom Thumb Competitors Information and Patients Neopuff masks operating instructions
24/08/15	Derek Lamb	15517	Tom Thumb Competitors Information and Patients Neopuff key references
24/08/15	Derek Lamb	15516	Tom Thumb Competitors Information and Patients Neopuff catalogue
24/08/15	Derek Lamb	15515	Tom Thumb Competitors Information and Patients Neopuff benefits
24/08/15	Derek Lamb	15512	Tom Thumb General Certificates Tom Thumb TT480 Instruction Booklet Dutch Project Completion certificate Omnilingua
24/08/15	Derek Lamb	15511	Tom Thumb Instructions for Use User Manual TT480 Instructions Dutch High Quality 5mm bleed

24/08/15	Derek Lamb	15510	Tom Thumb Instructions for Use User Manual TT480 Instructions Dutch High Quality
24/08/15	Derek Lamb	15509	Tom Thumb Instructions for Use User Manual TT480 Instructions Dutch Email
18/08/15	John Lamb	15497	Tom Thumb Literature reviews Clinical papers
14/08/15	Derek Lamb	15461	Tom Thumb Packaging Trials and validation
12/08/15	John Lamb	15454	Tom Thumb Literature reviews Airmed leaflet PIP & PEEP against bag
11/08/15	Derek Lamb	15439	Tom Thumb Risk / benefits assessment and Clinical Re evaluation
06/08/15	Derek Lamb	15390	VM3COP50.17 Tom Thumb 30CM limited QA Sheet
06/08/15	Derek Lamb	15385	VM3COP50.15 Packaging Procedures Tom Thumbs
05/08/15	Derek Lamb	15367	VM3COP50.12 Quality Control / Service Checks Tom Thumb
21/07/15	Derek Lamb	15183	No risks disposal
26/06/15	Derek Lamb	15118	Tom Thumb Trade Mark Certificate Tom thumb registration trade mark Image
26/06/15	Derek Lamb	15117	Tom Thumb Trade Mark Certificate Tom thumb registration trade mark TEXT
05/11/14	Derek Lamb	14506	Tom Thumb Trade Mark Certificate Trade Mark Company
05/11/14	Derek Lamb	14504	Tom Thumb Trade Mark Certificate
03/10/14	Derek Lamb	14437	Tom Thumb Promotional Leaflets 0330237 Viamed MH052 Leaflet
03/10/14	Derek Lamb	14433	Resuscitation Unit and TC400 Promotional Leaflets Tom Thumb Infant Resus Leaflet
15/07/14	Derek Lamb	14099	Tom Thumb Drawings 0330193 dial artwork
03/07/14	Derek Lamb	14091	Tom Thumb Possible replacement Parts Specification Sheets 0330193 Quotation
10/06/14	Derek Lamb	14046	Tom Thumb Description of Device
04/06/14	Derek Lamb	14011	Tom Thumb Clinical Trials Reports
04/06/14	John Lamb	14007	Tom Thumb Design Input Device History IMI positive pressure ventilator
04/06/14	John Lamb	14006	Tom Thumb Quotations IMI Pre Tom Thumb
06/05/14	Derek Lamb	13680	Tom Thumb Specifications Temperature Range
22/04/14	Derek Lamb	13654	Tom Thumb Competitor Instructions for use for similar products Fisher Pake NeoPuff IFU_22_04_14
31/10/13	Derek Lamb	13136	Medicinal Substances None
31/10/13	Derek Lamb	13135	No Animal Tissues No Animal Substances
14/10/13	Derek Lamb	13047	Tom Thumb Instructions for Use User Manual French
30/04/13	Derek Lamb	12348	Tom Thumb Preliminary drawings 0330191 Drawing possible replacement spacer.
30/04/13	Derek Lamb	12347	Tom Thumb Preliminary drawings 0330190 Drawing possible replacement clamp
14/03/13	Derek Lamb	12301	VM3COP50.10 Tom Thumb TT480 (0310094)
14/01/13	Derek Lamb	12231	Tom Thumb Drawings 0121184 (Note incorrectly refers to 0121198 in drawing)
14/01/13	Derek Lamb	12230	Tom Thumb Drawings 0121182 (note incorrectly refers to 0121201 in drawing)
14/01/13	Derek Lamb	12229	Tom Thumb Drawings 0121181 Drawing (note incorrectly refers to 0121201)
02/01/13	Derek Lamb	12226	Tom Thumb Preliminary drawings 0330055 14-001
02/01/13	Derek Lamb	12225	Tom Thumb Preliminary drawings 0330055 14-004
30/11/12	Derek Lamb	12211	Tom Thumb Instructions for Use User Manual TT495-15
23/11/12	Derek Lamb	12201	VM3COP50.091 Tom Thumb TT490 - TT480 Conversion Procedure Customer Reply Notification
24/10/12	john Lamb	12128	Tom Thumb Location of Documents Drawing register
24/10/12	john Lamb	12127	Tom Thumb Drawings Precision valve 0330210
24/10/12	john Lamb	12126	Tom Thumb Drawings Adjustable valve body 0330229
13/09/12	Derek Lamb	11903	VM3COP50.08 0310250 Convert TT490 to TT480
26/04/12	Derek Lamb	11775	Tom Thumb Expected Life of Products
02/04/12	steve nixon	11763	Tom Thumb Purchase specifications 0330203 Pressure Gauge
27/03/12	Derek Lamb	11748	Tom Thumb Labels 0390016 Address Sticker Round outside Gauge
27/03/12	Derek Lamb	11747	Tom Thumb Labels Caution Label 0390018 Top Block
27/03/12	Derek Lamb	11746	Tom Thumb Labels 0390017 Front of Gauge Sticker
27/03/12	Derek Lamb	11745	Tom Thumb Labels 0390015 Front Block Logo Sticker
27/03/12	Derek Lamb	11744	Tom Thumb Labels 0390010 Adjustable Valve Sticker
27/03/12	Derek Lamb	11743	Tom Thumb Labels Part number to Label
27/03/12	Derek Lamb	11742	Tom Thumb Labels TT491 Serial
27/03/12	Derek Lamb	11741	Tom Thumb Labels TT490 Serial
27/03/12	Derek Lamb	11740	Tom Thumb Labels Serial TT480
23/03/12	John Lamb	11732	Tom Thumb Instructions for Use User Manual
22/03/12	John Lamb	11669	Tom Thumb Product usability
20/03/12	John Lamb	11656	Tom Thumb Literature reviews Use of self-inflating bags

07/03/12	Derek Lamb	9781	0330240 Thumb Screw for Clamp
18/10/11	Derek Lamb	9533	Tom Thumb drawings Adjusting ring 0330232
18/10/11	Derek Lamb	9476	Tom Thumb Test Reports Dual Tom Thumb Flowrates
18/10/11	Derek Lamb	9475	Tom Thumb Tooling List
18/10/11	Derek Lamb	9474	Tom Thumb Leakage and Liquid damage
18/10/11	Derek Lamb	9456	Generic CE File Attached to All Ionising Radiation None
18/10/11	Derek Lamb	9455	Carcinogenic Substances
18/10/11	Derek Lamb	9441	Tom Thumb Essential requirements 200747ec
18/10/11	Derek Lamb	9358	Tom Thumb Pre-Clinical Trials Data
18/10/11	Derek Lamb	9140	Notified Body Certificate Not Required
18/10/11	Derek Lamb	9098	No 510k Required
12/10/11	Derek Lamb	8664	Service Manual Tom Thumb TT490-5 v1.4
06/07/11	John Lamb	8512	Tom Thumb Instructions for Use User Manual NeoPEEP
06/07/11	John Lamb	8506	Instructions For Use Tom Thumb TT495
04/07/11	John Lamb	8489	Tom Thumb EN ISO 14971:2001 Annex D
04/07/11	John Lamb	8486	Tom Thumb Work Instructions and Tests Work Instructions 2004
07/06/11	John Lamb	8315	Human Blood & derivatives Statement of non use
31/03/11	John Lamb	8211	Tom Thumb Service Manual Archive Service Manual 400 series Originals
16/03/11	John Lamb	8043	Machinery Directive Not Applicable
09/03/11	Derek Lamb	7916	Tom Thumb Service Procedures Repositioning the patient outlet on the Tom Thumb
09/03/11	Derek Lamb	7914	Tom Thumb Drawings 0330296
07/03/11	John Lamb	7912	Tom Thumb Preliminary drawings Drager TT495 Blender
07/03/11	John Lamb	7910	Tom Thumb Preliminary drawings TT85.
07/03/11	John Lamb	7908	Tom Thumb Work Instructions and Tests
07/03/11	John Lamb	7907	Tom Thumb Drawings 10 mm Inlet 0330208 0330248
07/03/11	John Lamb	7906	Tom Thumb Drawings 15mm Outlet 0330209 0330249
07/03/11	John Lamb	7902	Tom Thumb Preliminary drawings 9901-06 Blocking bolt
07/03/11	John Lamb	7901	Tom Thumb Preliminary drawings Adjustable valve
07/03/11	John Lamb	7900	Tom Thumb Preliminary drawings Preliminary Dwg Issue 4
07/03/11	John Lamb	7899	Tom Thumb Preliminary drawings Preliminary Dwg Issue 3
07/03/11	John Lamb	7898	Tom Thumb Preliminary drawings Preliminary Dwg Issue 2
07/03/11	John Lamb	7897	Tom Thumb Preliminary drawings Preliminary Dwg Issue 1
07/03/11	John Lamb	7896	Tom Thumb Design Reviews QC24
07/03/11	John Lamb	7894	Tom Thumb Drawings Drager 490
07/03/11	John Lamb	7893	Tom Thumb Drawings Drager TT460
07/03/11	John Lamb	7892	Tom Thumb Drawings Drager Gas control box
07/03/11	John Lamb	7891	Tom Thumb Preliminary drawings Modura Rail
07/03/11	John Lamb	7890	Tom Thumb Preliminary drawings Preliminary sketches
07/03/11	John Lamb	7889	Tom Thumb Parts List
07/03/11	John Lamb	7888	Tom Thumb Final Costings
07/03/11	John Lamb	7887	Tom Thumb Specifications of Materials Therapy Equipment Flow meter
07/03/11	John Lamb	7884	Tom Thumb Specifications of Materials TT component specifications all
07/03/11	John Lamb	7882	OLD DO NOT USE Tom Thumb Service Procedures
07/03/11	John Lamb	7880	Tom Thumb Training Information
24/02/11	Derek Lamb	7726	Tom Thumb Packaging Trials and validation TT495
24/02/11	Derek Lamb	7724	Tom Thumb Packaging Trials and validation Packagine TT480
15/02/11	John Lamb	7672	Tom Thumb Pre-Design Risk analysis report
14/02/11	Derek Lamb	7658	0310030 TT480 labelling modification 140211
14/02/11	Derek Lamb	7656	0310030 TT480 packaging modification 140211
21/01/11	Derek Lamb	7492	Tom Thumb Drawings Precision Valve Adjustable Screw 0330228
21/01/11	Derek Lamb	7491	Tom Thumb Drawings 1994 Archived complete drawings set
20/01/11	john lamb	7473	Tom Thumb EN Standards 2004
20/01/11	john lamb	7472	Tom Thumb Specification Brief Drager specs
20/01/11	john lamb	7471	Tom Thumb Full inspection and test 1997 QA TT495
20/01/11	john lamb	7469	Tom Thumb Preliminary drawings Flowmeter bolt 9901-02
20/01/11	john lamb	7468	Tom Thumb Drawings Drager Oxygen Air
20/01/11	john lamb	7467	Tom Thumb Drawings adjusting screw 0330230
20/01/11	john lamb	7465	Tom Thumb Preliminary drawings spacer Block 4 hole mount
20/01/11	john lamb	7464	Tom Thumb Preliminary drawings Dual outlet
20/01/11	john lamb	7463	Tom Thumb Quotations Buckingham fastners
20/01/11	john lamb	7462	Tom Thumb Quotations GR Hydro 1996
19/01/11	john lamb	7460	Tom Thumb Project Progress QC25 1997 Drager project
19/01/11	john lamb	7459	Tom Thumb Work Logs Drager project
19/01/11	john lamb	7458	Tom Thumb Design Input Device History valve theory

19/01/11	john lamb	7457	Tom Thumb Design Input Device History Air Oxygen calculations
19/01/11	john lamb	7456	Tom Thumb Purchase specifications 1997
19/01/11	john lamb	7455	Tom Thumb Tooling List pressure tester
19/01/11	john lamb	7454	Tom Thumb Tooling List manometer
19/01/11	john lamb	7453	Tom Thumb Components Breakdown parts Total photograph
19/01/11	john lamb	7452	Tom Thumb Product Changes G&R Hydroparts Block rounded corners
19/01/11	john lamb	7451	Tom Thumb Instructions for Use User Manual Application of Tom Thumb
19/01/11	john lamb	7450	Tom Thumb Analysis of complaints location of files
18/01/11	john lamb	7440	Tom Thumb Software
18/01/11	john lamb	7437	Tom Thumb Type Examinations
07/10/09	Derek Lamb	6584	Tom Thumb Possible replacement Parts Specification Sheets Guage2 Wiki
07/10/09	Derek Lamb	6583	Tom Thumb Possible replacement Parts Specification Sheets Guage1 Wiki
17/06/09	Steve Hardaker	6022	Tom Thumb User Feedback Blenders & Oxygen monitors
04/06/09	Derek Lamb	6009	Tom Thumb Final Drawings 10mm Tom Thumb outlet 20050605 (special)
14/05/09	Derek Lamb	5922	Tom Thumb Labels mediprema labels
12/03/09	Derek Lamb	5141	Tom Thumb Preliminary drawings 9901-34 Final Drawings 15mm Outlet Tom Thumb outlet 0330209
26/02/09	Derek Lamb	5120	Tom Thumb Clinical Information / Data
23/02/09	Derek Lamb	5107	Instructions For Use Tom Thumb TT480
15/12/08	Derek Lamb	4936	Tom Thumb Final Drawings TT0100_15_Retaining Screw 0330234 9901-14
15/07/08	Derek Lamb	4298	Tom Thumb Design Reviews QC24 Thumbelina Mediprema FeedBack 2008 June
09/07/08	John Lamb	4155	Tom Thumb Final Drawings Modura rail 2002 CE FILE Resuscitation Unit and TC400
30/06/08	Derek Lamb	3942	Tom Thumb Final Drawings TT0100_27_10mm Inlet 0330208
30/06/08	Derek Lamb	3940	Tom Thumb Final Drawings 9901-06 Blocking bolt and Nut 0330206
16/04/08	John Lamb	3593	Tom Thumb Final Drawings TT0100-Sp_02 special body block
28/03/08	John Lamb	3544	Tom Thumb Viamed Leaflets
19/03/08	John Lamb	3477	Tom Thumb Specifications of Materials TT0100-25 Blocking Bolt & Nut 0330206
31/01/08	Derek Lamb	3010	Tom Thumb Specifications of Materials 1998 GSM 4 001
31/01/08	Derek Lamb	3009	Tom Thumb Specifications of Materials 1997 Wika 4 001
31/01/08	Derek Lamb	3007	Tom Thumb Specifications of Materials 1997 Therapy Equipment 3 001
31/01/08	Derek Lamb	3005	Tom Thumb Specifications of Materials 1997 Harlow 001
31/01/08	Derek Lamb	3004	Tom Thumb Specifications of Materials 1997 G R Hydroparts 001
31/01/08	Derek Lamb	3003	Tom Thumb Specifications of Materials 1997 Buckingham fastners 001
31/01/08	Derek Lamb	3002	Tom Thumb Full inspection and test Drager Test Sheet
31/01/08	Derek Lamb	3001	Tom Thumb Full inspection and test Drager TT460 Inspection
31/01/08	Derek Lamb	3000	Tom Thumb Test Reports Flow checks Test Results
31/01/08	Derek Lamb	2999	Tom Thumb Test Reports TT Peep Results
31/01/08	Derek Lamb	2998	Tom Thumb Design Reviews QC24 Drager Design review
31/01/08	Derek Lamb	2997	Tom Thumb Validation design compliance
31/01/08	Derek Lamb	2996	Tom Thumb Validation 2005 validation
31/01/08	Derek Lamb	2993	Tom Thumb Working Drawings
31/01/08	Derek Lamb	2992	Tom Thumb Preliminary drawings Drager ALL Hand Drawings - Pre Cad
31/01/08	Derek Lamb	2990	Tom Thumb Quality Assurance Quality Agreements
31/01/08	Derek Lamb	2989	Tom Thumb Quality Assurance Newsham Engineering
31/01/08	Derek Lamb	2988	Tom Thumb Quality Assurance MAJ
31/01/08	Derek Lamb	2987	Tom Thumb Quality Assurance G R HydroParts
31/01/08	Derek Lamb	2986	Tom Thumb Photos Parts Photographs
31/01/08	Derek Lamb	2985	Tom Thumb Photos Midwife version
31/01/08	Derek Lamb	2984	Tom Thumb Photos TT490-15
31/01/08	Derek Lamb	2983	Tom Thumb Photos 1995 Originals
31/01/08	Derek Lamb	2982	Tom Thumb Literature reviews Neopuff
31/01/08	Derek Lamb	2981	Tom Thumb User Feedback 1997 Sunderland Clinical evaluation
31/01/08	Derek Lamb	2980	Tom Thumb Analysis of complaints 2001 Drager PEEP
31/01/08	Derek Lamb	2979	Tom Thumb Analysis of complaints 1997 QE Service Reports
31/01/08	Derek Lamb	2978	Tom Thumb Analysis of complaints 2006 James Cook excess pressures
31/01/08	Derek Lamb	2977	Tom Thumb Analysis of complaints 2004 Igi leaking oxygen
31/01/08	Derek Lamb	2973	Tom Thumb Training Information Training Manual Training Information

31/01/08	Derek Lamb	2970	Tom Thumb EN1998 Risk Analysis
31/01/08	Derek Lamb	2969	Tom Thumb MHRA Correspondence 2005 MHRA Statement 10mm connectors EN stds
31/01/08	Derek Lamb	2967	Tom Thumb Classification rationale 2004
11/07/07	Derek Lamb	2534	MH052 TDS UK Permabond
19/06/07	Derek Lamb	2475	Tom Thumb Compatibility Trials Sunderland 1997
19/06/07	Derek Lamb	2474	Tom Thumb Bio-compatibility
19/06/07	Derek Lamb	2471	Tom Thumb Test Reports
19/06/07	Derek Lamb	2470	Tom Thumb Final Costings GR Hydroparts
19/06/07	Derek Lamb	2469	Tom Thumb Final Costings
19/06/07	Derek Lamb	2466	Tom Thumb Clinical Trials Sunderland 1993
19/06/07	Derek Lamb	2465	Tom Thumb Specification Brief Specification_QC22 1993
19/06/07	Derek Lamb	2464	Tom Thumb Specification Brief Specification
19/06/07	Derek Lamb	2463	Tom Thumb Specification Brief Drager spec
19/06/07	Derek Lamb	2462	Tom Thumb Test Reports Wika calibration
19/06/07	Derek Lamb	2461	Tom Thumb Test Reports Wika cal p02
19/06/07	Derek Lamb	2460	Tom Thumb Test Reports Flush rays 2000
19/06/07	Derek Lamb	2459	Tom Thumb Validation Project Validation
19/06/07	Derek Lamb	2458	Tom Thumb Validation Compliance
19/06/07	Derek Lamb	2456	Tom Thumb Design Reviews QC24 TT Restrictor Year 2000
19/06/07	Derek Lamb	2455	Tom Thumb Design Reviews QC24 Product Review_03 March 2004
19/06/07	Derek Lamb	2454	Tom Thumb Design Reviews QC24 Product Review Restrictor 02 March 2000
19/06/07	Derek Lamb	2453	Tom Thumb Design Reviews QC24 Product Review 02 March 2000
19/06/07	Derek Lamb	2452	Tom Thumb Design Reviews QC24 Product Review 01 March 2000
19/06/07	Derek Lamb	2449	Tom Thumb Final Drawings TT0100_05_Precision Valve Assembly 0330210
19/06/07	Derek Lamb	2448	Tom Thumb Final Drawings TT0100_04_Valve Body 0330224
19/06/07	Derek Lamb	2447	Tom Thumb Final Drawings TT0100_03_Internal Valve 0330225
19/06/07	Derek Lamb	2446	Tom Thumb Final Drawings TT0100_01_Adjustment Screw 0330228
19/06/07	Derek Lamb	2442	Tom Thumb Final Drawings Y_11. TT0100_20 Body
19/06/07	Derek Lamb	2440	Tom Thumb Final Drawings TT0100_26_Flowmeter Bolt 0330205
19/06/07	Derek Lamb	2439	Tom Thumb Final Drawings TT0100_20_Body Block TT 480_490 0330201
19/06/07	Derek Lamb	2438	Tom Thumb Final Drawings TT0100_18_Spring 0330223
19/06/07	Derek Lamb	2435	Tom Thumb Final Drawings TT0100_16_Adjustable Valve Assembly 0330211
19/06/07	Derek Lamb	2433	Tom Thumb Final Drawings TT0100_14_Side Screw 0330234
19/06/07	Derek Lamb	2432	Tom Thumb Final Drawings TT0100_13_Internal Stop Nut 0330231
19/06/07	Derek Lamb	2431	Tom Thumb Final Drawings TT0100_10_Valve Body 0330229
19/06/07	Derek Lamb	2429	Tom Thumb Final Drawings TT0100_02_Valve Seat 0330227
19/06/07	Derek Lamb	2428	Tom Thumb Final Drawings Precision Valve Service Layout1 0330210
19/06/07	Derek Lamb	2427	Tom Thumb Final Drawings Final Drawings
19/06/07	Derek Lamb	2426	Tom Thumb Final Drawings Final diagrams 2004
19/06/07	Derek Lamb	2424	Tom Thumb Final Drawings Drager Resus
19/06/07	Derek Lamb	2422	Tom Thumb Final Drawings 9901-26
19/06/07	Derek Lamb	2421	Tom Thumb Final Drawings 9901-08
19/06/07	Derek Lamb	2420	Tom Thumb Final Drawings 9901-07
19/06/07	Derek Lamb	2419	Tom Thumb Final Drawings 9901-04
19/06/07	Derek Lamb	2417	Tom Thumb Final Drawings 0330215 Clamp Spacer 0330215
19/06/07	Derek Lamb	2416	Tom Thumb Final Drawings 9901-34
19/06/07	Derek Lamb	2415	Tom Thumb Final Drawings 9901-13
19/06/07	Derek Lamb	2414	Tom Thumb Final Drawings 9901-12
19/06/07	Derek Lamb	2413	Tom Thumb Final Drawings 9901-11
19/06/07	Derek Lamb	2412	Tom Thumb Final Drawings Indicator nut 9901-10
19/06/07	Derek Lamb	2411	Tom Thumb Final Drawings 9901-09
19/06/07	Derek Lamb	2410	Tom Thumb Final Drawings 9901-06
19/06/07	Derek Lamb	2409	Tom Thumb Final Drawings 9901-05
19/06/07	Derek Lamb	2408	Tom Thumb Final Drawings 9901-03
19/06/07	Derek Lamb	2407	Tom Thumb Preliminary drawings TT490-15 1999
19/06/07	Derek Lamb	2406	Tom Thumb Preliminary drawings TT490 v1.1 1993
19/06/07	Derek Lamb	2405	Tom Thumb Preliminary drawings TT480 modification
19/06/07	Derek Lamb	2404	Tom Thumb Preliminary drawings TT179
19/06/07	Derek Lamb	2403	Tom Thumb Preliminary drawings Sketches 1993
19/06/07	Derek Lamb	2402	Tom Thumb Preliminary drawings screw
19/06/07	Derek Lamb	2401	Tom Thumb Preliminary drawings Saturday, December 01, 2001 (2)

19/06/07	Derek Lamb	2400	Tom Thumb Layout Designer
19/06/07	Derek Lamb	2399	Tom Thumb Preliminary drawings Relief valve v2
19/06/07	Derek Lamb	2398	Tom Thumb Preliminary drawings Preliminary Drawings
19/06/07	Derek Lamb	2397	Tom Thumb Preliminary drawings Prelim Drawings 1993
19/06/07	Derek Lamb	2396	Tom Thumb Preliminary drawings Modura clamp
19/06/07	Derek Lamb	2395	Tom Thumb Preliminary drawings Main bolt
19/06/07	Derek Lamb	2394	Tom Thumb Preliminary drawings Gas outlet
19/06/07	Derek Lamb	2393	Tom Thumb Preliminary drawings Flow meter & clamp
19/06/07	Derek Lamb	2392	Tom Thumb Preliminary drawings Drager TT490
19/06/07	Derek Lamb	2391	Tom Thumb Preliminary drawings Drager TT460 a
19/06/07	Derek Lamb	2390	Tom Thumb Preliminary drawings Drager PEEP problems 2001
19/06/07	Derek Lamb	2389	Tom Thumb Preliminary drawings Drager mixer Prelim
19/06/07	Derek Lamb	2388	Tom Thumb Preliminary drawings Drager DRF001
19/06/07	Derek Lamb	2387	Tom Thumb Preliminary drawings Adjustable valve 2001
19/06/07	Derek Lamb	2386	Tom Thumb Preliminary drawings 9901-39
19/06/07	Derek Lamb	2384	Tom Thumb Preliminary drawings 9901-26
19/06/07	Derek Lamb	2383	Tom Thumb Preliminary drawings 9901-13
19/06/07	Derek Lamb	2382	Tom Thumb Preliminary drawings 9901-12
19/06/07	Derek Lamb	2381	Tom Thumb Preliminary drawings valve body 9901-11
19/06/07	Derek Lamb	2380	Tom Thumb Preliminary drawings 9901-10
19/06/07	Derek Lamb	2379	Tom Thumb Preliminary drawings 9901-09
19/06/07	Derek Lamb	2378	Tom Thumb Preliminary drawings 9901-08b
19/06/07	Derek Lamb	2377	Tom Thumb Preliminary drawings 9901-06
19/06/07	Derek Lamb	2376	Tom Thumb Preliminary drawings 9901-05
19/06/07	Derek Lamb	2375	Tom Thumb Preliminary drawings 9901-04b
19/06/07	Derek Lamb	2374	Tom Thumb Preliminary drawings 9901-03
19/06/07	Derek Lamb	2373	Tom Thumb Preliminary drawings 14 September 2001b
19/06/07	Derek Lamb	2372	Tom Thumb Quotations Quotations 1993
19/06/07	Derek Lamb	2371	Tom Thumb Expenditure Purchases 1996
19/06/07	Derek Lamb	2370	Tom Thumb Expenditure Expenditure details 1993
19/06/07	Derek Lamb	2369	Tom Thumb Project Diary Project Progress
19/06/07	Derek Lamb	2368	Tom Thumb Design Input Time Scale Time scales 2004
19/06/07	Derek Lamb	2367	Tom Thumb Design Input Time Scale Time scale 1997
19/06/07	Derek Lamb	2366	Tom Thumb Work Logs Job Progress
19/06/07	Derek Lamb	2365	Tom Thumb Design Input Device History Restrictor calculations
19/06/07	Derek Lamb	2364	Tom Thumb Design Input Device History Original Specification Relief pressures
19/06/07	Derek Lamb	2363	Tom Thumb Design Input Device History Flow calculations
19/06/07	Derek Lamb	2362	Tom Thumb Design Input Device History Device Theory 2005
19/06/07	Derek Lamb	2361	Tom Thumb Design Input Device History Device History & Theory 2004
19/06/07	Derek Lamb	2360	Tom Thumb Design Input Device History Calculations of Flow
19/06/07	Derek Lamb	2359	Tom Thumb Design Input OEM Files Flow Rates
19/06/07	Derek Lamb	2358	Tom Thumb Design Input OEM Files Documentation Checklist 1995
19/06/07	Derek Lamb	2357	Tom Thumb Design Input OEM Files Checklist 1993
19/06/07	Derek Lamb	2356	Tom Thumb Purchase specifications R S
19/06/07	Derek Lamb	2355	Tom Thumb Purchase specifications Wika
19/06/07	Derek Lamb	2354	Tom Thumb Purchase specifications Purchase Specifications 2004
19/06/07	Derek Lamb	2353	Tom Thumb Purchase specifications Newsham Engineering
19/06/07	Derek Lamb	2351	Tom Thumb Purchase specifications Pressure Tester
19/06/07	Derek Lamb	2350	Tom Thumb Purchase specifications M A J Hi-spec
19/06/07	Derek Lamb	2349	Tom Thumb Purchase specifications Intersurgical
19/06/07	Derek Lamb	2348	Tom Thumb Purchase specifications G & R Hydroparts
19/06/07	Derek Lamb	2347	Tom Thumb Purchase specifications Fromlim
19/06/07	Derek Lamb	2346	Tom Thumb Purchase specifications CPC
19/06/07	Derek Lamb	2345	Tom Thumb Purchase specifications Chelmer Pneumatics
19/06/07	Derek Lamb	2344	Tom Thumb Purchase specifications Buckingham fastners
19/06/07	Derek Lamb	2343	Tom Thumb Specifications of Materials TE guage 1997 9901-15
19/06/07	Derek Lamb	2342	Tom Thumb Specifications of Materials TE Flowmeters
19/06/07	Derek Lamb	2341	Tom Thumb Specifications of Materials T E Copper washers
19/06/07	Derek Lamb	2340	Tom Thumb Specifications of Materials O rings 1997
19/06/07	Derek Lamb	2338	Tom Thumb Specifications of Materials Fomblin grease
19/06/07	Derek Lamb	2337	Tom Thumb Specifications of Materials Dial a flow 1994
19/06/07	Derek Lamb	2336	Tom Thumb Components Breakdown THE 2002 Modurail Clamp
19/06/07	Derek Lamb	2335	Tom Thumb Components Breakdown Spacer THE2002
19/06/07	Derek Lamb	2334	Tom Thumb Components Breakdown O rings

19/06/07	Derek Lamb	2333	Tom Thumb Components Breakdown Modura rail screw
19/06/07	Derek Lamb	2332	Tom Thumb Components Breakdown 99091-06 nut for seating
19/06/07	Derek Lamb	2331	Tom Thumb Components Breakdown 99091-06 nut Nut for seating
19/06/07	Derek Lamb	2330	Tom Thumb Components Breakdown 9901-40 precision vavle breakdown
19/06/07	Derek Lamb	2328	Tom Thumb Components Breakdown 9901-38b Adjustment screw
19/06/07	Derek Lamb	2326	Tom Thumb Components Breakdown 9901-34 adapter 15mm male
19/06/07	Derek Lamb	2325	Tom Thumb Components Breakdown 9901-31 Blanking Bolt
19/06/07	Derek Lamb	2324	Tom Thumb Components Breakdown 9901-28 Copper sealing washer (diamond)
19/06/07	Derek Lamb	2323	Tom Thumb Components Breakdown 9901-27 Main valve spring
19/06/07	Derek Lamb	2322	Tom Thumb Components Breakdown 9901-26 Main Valve spring
19/06/07	Derek Lamb	2319	Tom Thumb Components Breakdown 9901-20 Rt Angle connector for O2
19/06/07	Derek Lamb	2318	Tom Thumb Components Breakdown 9901-19 Flowmeter bolt with hole
19/06/07	Derek Lamb	2317	Tom Thumb Components Breakdown 9901-18
19/06/07	Derek Lamb	2316	Tom Thumb Components Breakdown 9901-17
19/06/07	Derek Lamb	2315	Tom Thumb Components Breakdown 9901-16 O2 Hose 1mtr
19/06/07	Derek Lamb	2314	Tom Thumb Components Breakdown 9901-15 Wika guage
19/06/07	Derek Lamb	2313	Tom Thumb Components Breakdown 9901-14 Screw
19/06/07	Derek Lamb	2312	Tom Thumb Components Breakdown 9901-13 Main valve side screw
19/06/07	Derek Lamb	2310	Tom Thumb Components Breakdown 9901-12 Main valve seat
19/06/07	Derek Lamb	2307	Tom Thumb Components Breakdown 9901-11 (2) Main valve internal
19/06/07	Derek Lamb	2305	Tom Thumb Components Breakdown 9901-10b Main valve internal stop
19/06/07	Derek Lamb	2304	Tom Thumb Components Breakdown 9901-10a Main valve internal stop
19/06/07	Derek Lamb	2303	Tom Thumb Components Breakdown 9901-10 Main valve internal stop
19/06/07	Derek Lamb	2302	Tom Thumb Components Breakdown 9901-09 (2) Main valve adjusting screw
19/06/07	Derek Lamb	2300	Tom Thumb Components Breakdown 9901-08a Main valve housing
19/06/07	Derek Lamb	2298	Tom Thumb Components Breakdown 9901-06 (2) Safety Valve seat
19/06/07	Derek Lamb	2296	Safety valve seat Tom Thumb Components Breakdown 9901-05-b Safety valve seat
19/06/07	Derek Lamb	2294	Tom Thumb Components Breakdown 9901-04 (2) Main valve internal
19/06/07	Derek Lamb	2293	Tom Thumb Components Breakdown 9901-04 Main valve adjusting ring
19/06/07	Derek Lamb	2292	Tom Thumb Components Breakdown 9901-03-c Safety valve housing
19/06/07	Derek Lamb	2291	Tom Thumb Components Breakdown 9901-03 Variable valve housing
19/06/07	Derek Lamb	2289	Tom Thumb Components Breakdown 9901-01 Main block
19/06/07	Derek Lamb	2288	Tom Thumb Components Breakdown 9901-00 Nut for seating
19/06/07	Derek Lamb	2287	Tom Thumb Components Breakdown 6575 TE Flowmeters b
19/06/07	Derek Lamb	2285	Tom Thumb Components Breakdown 4701-61 O2 Inlet Christmas tree
19/06/07	Derek Lamb	2284	Tom Thumb Components Breakdown 4450-12 Blanking nut
19/06/07	Derek Lamb	2281	Tom Thumb Accessories Photographs DSCF1001
19/06/07	Derek Lamb	2280	Tom Thumb Accessories Photographs DSCF1000
19/06/07	Derek Lamb	2279	Tom Thumb Accessories Photographs DSCF0999
19/06/07	Derek Lamb	2278	Tom Thumb Accessories Photographs DSCF0998
19/06/07	Derek Lamb	2277	Tom Thumb Accessories Photographs DSCF0997
19/06/07	Derek Lamb	2276	Tom Thumb Accessories Photographs DSCF0996
19/06/07	Derek Lamb	2275	Tom Thumb Complete Unit Photographs TT490
19/06/07	Derek Lamb	2274	Tom Thumb Complete Unit Photographs TT480
19/06/07	Derek Lamb	2273	Tom Thumb Complete Unit Photographs Originals
19/06/07	Derek Lamb	2272	Tom Thumb Complete Unit Photographs Midwife
19/06/07	Derek Lamb	2271	Tom Thumb Complete Unit Photographs Complete Units 2004
19/06/07	Derek Lamb	2270	Tom Thumb Complete Unit Photographs Complete Units 2002
19/06/07	Derek Lamb	2269	Tom Thumb Components Breakdown
19/06/07	Derek Lamb	2265	Tom Thumb Parts List Pictorial Info Tom Thumb Parts List TT490
19/06/07	Derek Lamb	2264	Tom Thumb Parts List Parts List 2005
19/06/07	Derek Lamb	2262	Tom Thumb Parts List Parts LIST 2001
19/06/07	Derek Lamb	2261	Tom Thumb Parts List Parts List 1994
19/06/07	Derek Lamb	2260	Tom Thumb Parts List Parts list 1993 original
19/06/07	Derek Lamb	2258	Tom Thumb Parts List O_Rings02
19/06/07	Derek Lamb	2256	Tom Thumb Parts List Copper Washer_03
19/06/07	Derek Lamb	2255	Tom Thumb Parts List Block
19/06/07	Derek Lamb	2254	Tom Thumb Parts List Adj label 04
19/06/07	Derek Lamb	2253	Tom Thumb Sterilisation

19/06/07	Derek Lamb	2251	Tom Thumb Quality Plan Quality Plan 2004
19/06/07	Derek Lamb	2250	Tom Thumb Quality Plan N. MAJ_HI Spec
19/06/07	Derek Lamb	2249	Tom Thumb Quality Plan N. G & R Hydroparts
19/06/07	Derek Lamb	2248	Tom Thumb Quality Plan N. D. Newsham
19/06/07	Derek Lamb	2247	Tom Thumb Specifications Technical Specification 1997
19/06/07	Derek Lamb	2246	Tom Thumb Photos Tom thumb in use
19/06/07	Derek Lamb	2245	Tom Thumb Packaging Trials and validation Packaging Photographs 2004
19/06/07	Derek Lamb	2244	Tom Thumb Packaging Trials and validation Packaging Trial 2004
19/06/07	Derek Lamb	2243	Tom Thumb Tooling List
19/06/07	Derek Lamb	2242	Tom Thumb Test Method QA Spring 9901-26
19/06/07	Derek Lamb	2240	Tom Thumb Work Instructions and Tests Assembly TT480
19/06/07	Derek Lamb	2239	Tom Thumb Work Instructions and Tests Assembly valve 9901-40
19/06/07	Derek Lamb	2238	Tom Thumb Work Instructions and Tests Assembly valve 9901-18
19/06/07	Derek Lamb	2237	Tom Thumb Work Instructions and Tests Assembly TT495
19/06/07	Derek Lamb	2236	Tom Thumb Work Instructions and Tests Assembly TT490
19/06/07	Derek Lamb	2235	Tom Thumb Work Instructions and Tests Assembly TT460
19/06/07	Derek Lamb	2233	Tom Thumb Location of responsibilities
19/06/07	Derek Lamb	2230	Tom Thumb Clinical Trials RVI 1997
19/06/07	Derek Lamb	2229	Tom Thumb Clinical Trials (R) Clinical Trials 1997
19/06/07	Derek Lamb	2228	Tom Thumb User Feedback Training Manual
19/06/07	Derek Lamb	2227	Tom Thumb Analysis of complaints Complaints RVI Outlet tube size
19/06/07	Derek Lamb	2226	Tom Thumb Analysis of complaints Repair Reports 1997-a
19/06/07	Derek Lamb	2225	Tom Thumb Analysis of complaints Complaints Sunderland pressure Gauge
19/06/07	Derek Lamb	2223	Tom Thumb Analysis of complaints Complaints Internal 1994
19/06/07	Derek Lamb	2222	Tom Thumb Analysis of complaints Complaint Stoke City General Pressure relief valve
19/06/07	Derek Lamb	2221	Tom Thumb Product Changes Product changes 2000 Restrictor
19/06/07	Derek Lamb	2220	Tom Thumb Product Changes Product Changes 1998 Flow meter
19/06/07	Derek Lamb	2215	Tom Thumb Product Life
19/06/07	Derek Lamb	2212	Tom Thumb Maintenance Service & Care of Tom Thumb 1994
19/06/07	Derek Lamb	2211	Tom Thumb Maintenance Oxygen Hose Maintenance
19/06/07	Derek Lamb	2210	Tom Thumb Maintenance medical-devices.gov
19/06/07	Derek Lamb	2209	Tom Thumb Maintenance
19/06/07	Derek Lamb	2204	Tom Thumb Labels Tomb Thumb Infant Resuscitator
19/06/07	Derek Lamb	2188	Tom Thumb Description of Accessories Datex-Ohmeda Giraffe Incubator
19/06/07	Derek Lamb	2187	Tom Thumb Viamed Data Sheets
19/06/07	Derek Lamb	2186	Tom Thumb Viamed Leaflets 2
19/06/07	Derek Lamb	2185	Tom Thumb Viamed Leaflets 1
19/06/07	Derek Lamb	2183	Tom Thumb EN ISO 14971:2001 Annex A
19/06/07	Derek Lamb	2180	Tom Thumb Risk analysis certificate
19/06/07	Derek Lamb	2178	Tom Thumb MHRA Correspondence MDA_2003_016
19/06/07	Derek Lamb	2177	Tom Thumb MHRA Correspondence MDA 11_03_03 Medical gas hoses
19/06/07	Derek Lamb	2176	Tom Thumb EMC rationale
19/06/07	Derek Lamb	2175	Tom Thumb Classification rationale
19/06/07	Derek Lamb	2174	Tom Thumb EN Standards
19/06/07	Derek Lamb	2172	Tom Thumb Essential requirements
19/06/07	Derek Lamb	2168	Tom Thumb Location of Documents Index
19/06/07	Derek Lamb	2167	Tom Thumb Location of Documents
19/06/07	Derek Lamb	2164	Tom Thumb CD Amendment Control
21/05/07	Derek Lamb	1663	Tom Thumb Description of Accessories silicone adaptor

Comments on Document Changes

no issues with documentations

Further Action Required ☐**Risks Associated with Design Changes / Product Changes**

Further Action Required ☐**Issues containing References Created / Updated since Review Date****Comments on Issues**

mainly just servicing issues, no problems

Further Action Required ☐**Risks Associated with Issues**Further Action Required ☐

Clinical Re-Assessment:

Search Google for [Clinical Reports](#)**Suggest Advanced search and narrow results to last 12 months**

Do any of the Results indicate a Risk / Problem with Tom Thumb

☐ Yes ☒ No

Do any of the Results indicate Tom Thumb is outdated Technology

☐ Yes ☒ No

Comments on Googles responses (include links to main pages)

FDA Reports:

Search Google for [FDA Reports](#)**Suggest Advanced search and narrow results to last 12 months**

Do any of the Results indicate a Risk / Problem with Tom Thumb

☐ Yes ☒ No

Do any of the Results indicate Tom Thumb is outdated Technology

☐ Yes ☒ No

Comments on Googles responses (include links to main pages)

Sign off / Store review Derek Lamb