

Description of Device.

Generic name of device: Timer

Model name: Apgar Timer

Model number: 0310100

Brief description of device and other devices designed to be used in combination with it:

It is intended for use with Resuscitation cabinets in Delivery rooms, Operating theatres, and Special Baby Care Units etc.

It is a Start / Freeze / Reset clock used in the compilation of “Apgar Scoring” on newborn infants based on their vital signs, and response to stimulus.

It is designed to be mounted within the cabinet on standard medi-rails via a mounting plate, alternatively it can be stand alone, by the use of the 2 in built camera mount studs.

The scoring of “Apgar” is a recognized technique and is therefore not explained in this document.