
VIAMED Management Meeting - CE Files Review

D4 EMC

Issue ID #62543

Date Created 08/Sep/2015

Issued To: Derek Lamb

Completed Status: Still Outstanding

Notes:

Added by John Lamb sent to Derek Lamb

This may need to be re-visited and expanded.

1) It can only be tested when operated.

2) A operator can not be in the room whilst EMC tests are being carried out.

3) The leads are attached to electrodes to the ULnar nerve which will act as an aerial. The pulses are up to 100v very fast but very short.

4) any electrical connections to internally circuitry to switch on the unit will effect EMC readings

4) Only way to test would be to design & manufacture a remote controlled finger. This would have to be electromechanical and could also effect the EMC results.
