

Viamed Limited · 15 Station Road · Cross Hills Keighley · West Yorkshire · BD20 7DT · United Kingdom Tel: +44 (0)1535 634 542 · Fax: +44 (0)1535 635 582 Email: info@viamed.co.uk Website: www.viamed.co.uk

Service Repair Sheet 68879

Contact Name

Company/ Hospital Name

Direct Phone General Phone Opera Account

Email

Date Received Booked in By

Clay Harper

Sage Services Group

8432845844 8772817243 00006261

harper@sageservicesgroup.com

28/Aug/2024 Robert Connor

Repair SRN37066

Ref

S/N

1410000 PR01883A18

Equipment Type

Foetal Simulator

Complete - Repaired Time: 0 Hour(s)

Parts Replaced

Qty: 1 1480000 Qty: 1 1430309

V1000 has had transducer interface cushion replaced and calibration has been checked.