

Viamed Limited + 15 Station Road - Cross Hills Keighley · West Yorkshire · BD20 7DT · United Kingdom Tel: +44 (0)1535 634 542 - Fax: +44 (0)1535 635 582 Email: info@viamed.co.uk Website: www.viamed.co.uk

Service Repair Sheet 68890

Contact Name

Company/ Hospital Name

Position Direct Phone Opera Account

Email

Order Number Date Received Booked in By

Derya Delmont

Queens Hospital/Siemens Medical Solutions

Technical Administrator

01708435482 00000264

derya.delmont.ext@siemens-healthineers.com

4515414434 09/Sep/2024 Robert Connor

Repair SRN37073 1410000 PR065A25

Ref

S/N

Equipment Type Foetal Simulator

Complete - Repaired Time: 0 Hour(s)

Parts Replaced

1480000 Qty: 1 1430309 Qty: 1

V1000 has had transducer interface cushion replaced and calibration has been checked.