

Viamed Limited · 15 Station Road · Cross Hills Keighley · West Yorkshire · BD20 7DT · United Kingdom Tel: +44 (0)1535 634 542 · Fax: +44 (0)1535 635 582 Email: Info@viamed.co.uk Website: www.viamed.co.uk

Service Repair Sheet 68923 Contact Name Company/ Hospital Name Position Direct Phone General Phone

Opera Account Email Date Received Booked in By Abbey McAnaney Vet Direct Services Ltd Technical Support RVN 01912642705 01912642705 00012824 tech@vet-direct.com 07/Oct/2024 Robert Connor

Repair Ref S/N SRN37126 0111275 A106337 Equipment Type Oxygen Analyser Under Warranty Complete - Repaired Time :0 Hour(s)

Device was reading 86.7% in air and couldn't be adjusted lower.

I have performed a reset and recentred the device on 20.9%

Repair Ref S/N SRN37127 0110660 A111126 Equipment Type Oxygen Sensor

No Fault Found Time :0 Hour(s)

The sensor is within spec. No fault found.