



Neuromuscular Monitor

Quantitatively monitors the level of neuromuscular transmission during surgery or in the intensive care unit, as a means of assessing the depth or degree of muscle relaxation



- Economical, uses standard low cost ECG electrodes
- Twitch response measured by acceleromyography (AMG)
- Displayed values
- Trend data





The new TOF3D is easy to use and, in addition to the standard TOF and PTC modes, there are alternative objective and subjective modes of stimulation available.

The monitor's functions and modes of stimulation can be quickly customised to meet your specific requirements.

+44 (0)1535 634542 www.viamed.co.uk