



Kate Griffiths <viamed.kate.griffiths@gmail.com>

RE: MAX 550E

2 messages

sales@medisyst.co.ke <sales@medisyst.co.ke>

6 November 2024 at 10:53

To: Kate Griffiths <kate.griffiths@viamed.co.uk>

Cc: williamomolo@medisyst.co.ke, watanga@icloud.com, kennedy@medisyst.co.ke

Hello Kate,

What I need is the sensor cell MAX 550E. Yes we have ordered it previously for use with the blenders that we had bought.

Regards,

Kennedy

From: Kate Griffiths <kate.griffiths@viamed.co.uk>**Sent:** Wednesday, November 6, 2024 1:39 PM**To:** sales@medisyst.co.ke**Cc:** williamomolo@medisyst.co.ke; watanga@icloud.com; kennedy@medisyst.co.ke**Subject:** Re: MAX 550E

Good morning Kennedy,

Thank you for your email.

Do you require part number 0110454 - Maxtec Oxygen Sensor MAX-550E with specific flow diverter for use with Maxtec MaxBlend range of blenders? You have ordered this previously.

There are two versions of the MAX-550E.

0110454 (R140P02-001) is a Max-550E for use with MaxBlend devices (MaxBlend 2, MaxBlend Lite and Blender Buddy 2).

It is supplied complete with a special flow diverter which is narrower in diameter and has a small bleed hole (0120086).