



VIA MED

Pulse Oximetry Sensors 4000 Serie

Reusable and disposable pulse oximetry sensors for continuous and spot-check monitoring of functional arterial oxygen saturation and pulse rate.



Class IIb

Assesment Route Annex II
NBOG MD1302

Carried out by Derek Lamb
20 / 04 / 16



VIAMED

-  Stock References Review
-  Supplier Review
-  Sales Review
-  Countries Review
-  Returns Review
-  Design Changes Review
-  User Instructions Review
-  Labels Review
-  Document Updates Review
-  Internal Issues Review
-  Clinical / FDA Incidents online search
-  Risk ISO 14971 : 2012 Review

Stock References Review

Stock Reference	Description
0014590	SpO2 Silicone Finger Sensor - Adult.
0014591	SpO2 Silicone Finger Sensor - Paediatric
0014592	SpO2 Wrap Sensor- Paediatric
0014595	SpO2 Silicone Finger Sensor - Adult.
0014596	SpO2 Silicone Finger Sensor - Paediatric
0014601	SpO2 Finger Sensor - Adult.
0014602	SpO2 Finger Sensor - Adult.
0014603	SpO2 Finger Sensor - Adult.
0014604	SpO2 Finger Sensor - Adult.
0014605	SpO2 Finger Sensor - Adult.
0014606	SpO2 Finger Sensor - Adult.
0014607	SpO2 Finger Sensor - Adult.
0014608	SpO2 Finger Sensor - Adult.
0014609	SpO2 Finger Sensor - Adult.
0014610	SpO2 Finger Sensor - Adult.
0014611	SpO2 Finger Sensor - Adult.
0014612	SpO2 Finger Sensor - Adult.
0014613	SpO2 Finger Sensor - Adult.
0014614	SpO2 Finger Sensor - Adult.
0014615	SpO2 Finger Sensor - Adult.
0014617	SpO2 Finger Sensor - Adult.
0014618	SpO2 Finger Sensor - Adult.
0014619	SpO2 Finger Sensor - Adult.
0014620	SpO2 Finger Sensor - Adult.
0014621	SpO2 Finger Sensor - Adult.
0014623	SpO2 Finger Sensor - Adult.
0014625	SpO2 Finger Sensor - Adult.
0014649	SpO2 Finger Sensor - Adult.
0014650	SpO2 Finger Sensor - Adult.
0014701	SpO2 Silicone Finger Sensor - Adult.
0014702	SpO2 Silicone Finger Sensor - Adult.
0014703	SpO2 Silicone Finger Sensor - Adult.
0014704	SpO2 Silicone Finger Sensor - Adult.
0014705	SpO2 Silicone Finger Sensor - Adult.
0014706	SpO2 Silicone Finger Sensor - Adult.
0014707	SpO2 Silicone Finger Sensor - Adult.
0014708	SpO2 Silicone Finger Sensor - Adult.
0014709	SpO2 Silicone Finger Sensor - Adult.
0014710	SpO2 Silicone Finger Sensor - Adult.
0014711	SpO2 Silicone Finger Sensor - Adult.
0014712	SpO2 Silicone Finger Sensor - Adult.
0014713	SpO2 Silicone Finger Sensor - Adult.

Stock Reference	Description
0014714	SpO2 Silicone Finger Sensor - Adult.
0014715	SpO2 Silicone Finger Sensor - Adult.
0014717	SpO2 Silicone Finger Sensor - Adult.
0014718	SpO2 Silicone Finger Sensor - Adult.
0014719	SpO2 Silicone Finger Sensor - Adult.
0014720	SpO2 Silicone Finger Sensor - Adult.
0014721	SpO2 Silicone Finger Sensor - Adult.
0014722	SpO2 Silicone Finger Sensor - Adult.
0014723	SpO2 Silicone Finger Sensor - Adult.
0014724	SpO2 Silicone Finger Sensor - Adult.
0014725	SpO2 Silicone Finger Sensor - Adult.
0014740	SpO2 Silicone Finger Sensor - Adult.
0014750	SpO2 Silicone Finger Sensor - Adult.
0014751	SpO2 Silicone Finger Sensor - Paediatric
0014755	SpO2 Silicone Finger Sensor - Adult.
0014756	SpO2 Silicone Finger Sensor - Paediatric
0014761	SpO2 Silicone Finger Sensor - Paediatric
0014762	SpO2 Silicone Finger Sensor - Paediatric
0014763	SpO2 Silicone Finger Sensor - Paediatric
0014764	SpO2 Silicone Finger Sensor - Paediatric
0014765	SpO2 Silicone Finger Sensor - Paediatric
0014766	SpO2 Silicone Finger Sensor - Paediatric
0014767	SpO2 Silicone Finger Sensor - Paediatric
0014768	SpO2 Silicone Finger Sensor - Paediatric
0014769	SpO2 Silicone Finger Sensor - Paediatric
0014770	SpO2 Silicone Finger Sensor - Paediatric
0014771	SpO2 Silicone Finger Sensor - Paediatric
0014772	SpO2 Silicone Finger Sensor - Paediatric
0014773	SpO2 Silicone Finger Sensor - Paediatric
0014774	SpO2 Silicone Finger Sensor - Paediatric
0014775	SpO2 Silicone Finger Sensor - Paediatric
0014777	SpO2 Silicone Finger Sensor - Paediatric
0014778	SpO2 Silicone Finger Sensor - Paediatric
0014779	SpO2 Silicone Finger Sensor - Paediatric
0014780	SpO2 Silicone Finger Sensor - Paediatric
0014781	SpO2 Silicone Finger Sensor - Paediatric
0014782	SpO2 Silicone Finger Sensor - Paediatric
0014785	SpO2 Silicone Finger Sensor - Paediatric
0014790	SpO2 Silicone Finger Sensor - Paediatric
0014801	Pulse Oximetry Ear Sensor
0014802	Pulse Oximetry Ear Sensor
0014803	Pulse Oximetry Ear Sensor
0014804	Pulse Oximetry Ear Sensor
0014805	Pulse Oximetry Ear Sensor
0014806	Pulse Oximetry Ear Sensor

Stock Reference	Description
0014807	Pulse Oximetry Ear Sensor
0014808	Pulse Oximetry Ear Sensor
0014816	Pulse Oximetry Ear Sensor
0014820	Pulse Oximetry Wrap Sensor
0014821	Pulse Oximetry Wrap Sensor
0014822	Pulse Oximetry Wrap Sensor
0014823	Pulse Oximetry Wrap Sensor
0014824	Pulse Oximetry Wrap Sensor
0014825	Pulse Oximetry Ear Sensor
0014826	Pulse Oximetry Wrap Sensor
0014827	Pulse Oximetry Wrap Sensor
0014828	Pulse Oximetry Wrap Sensor
0014829	Pulse Oximetry Wrap Sensor
0014830	Pulse Oximetry Wrap Sensor
0014831	Pulse Oximetry Wrap Sensor
0014832	Pulse Oximetry Wrap Sensor
0014833	Pulse Oximetry Wrap Sensor
0014834	Pulse Oximetry Wrap Sensor
0014835	Pulse Oximetry Wrap Sensor
0014836	Pulse Oximetry Wrap Sensor
0014837	Pulse Oximerty Wrap Sensor
0014838	Pulse Oximetry Wrap Sensor
0014840	Pulse Oximetry Wrap Sensor
0014841	Pulse Oximetry Wrap Sensor
0014842	Pulse Oximetry Wrap Sensor
0014843	Pulse Oximetry Wrap Sensor
0014844	Pulse Oximetry Wrap Sensor
0014845	Pulse Oximetry Multi-Site Y Sensor
0014848	Pulse Oximetry Wrap Sensor
0014895	Pulse Oximetry Extension Cable 1.2 m
0014896	Pulse Oximetry Extension Cable 2.4 m
0014901	Disposable SpO2 Sensors - Adult.
0014902	Disposable SpO2 Sensors - Paediatric.
0014903	Disposable SpO2 Sensors - Infant.
0014904	Disposable SpO2 Sensors - Neonatal.
0014905	Disposable SpO2 Sensors - Adult.
0014906	Disposable SpO2 Sensors - Paediatric.
0014909	Disposable SpO2 Sensors - Adult.
0014910	Disposable SpO2 Sensors - Paediatric.
0014911	Disposable SpO2 Sensors - Infant.
0014912	Disposable SpO2 Sensors - Neonatal.
0014913	Disposable SpO2 Sensors - Adult.
0014914	Disposable SpO2 Sensors - Paediatric.
0014915	Disposable SpO2 Sensors - Infant.
0014916	Disposable SpO2 Sensors - Neonatal.

Stock Reference	Description
0014917	Disposable SpO2 Sensors - Adult.
0014918	Disposable SpO2 Sensors - Paediatric.
0014921	Disposable SpO2 Sensors - Adult.
0014922	Disposable SpO2 Sensors - Paediatric.
0014923	Disposable SpO2 Sensors - Infant.
0014924	Disposable SpO2 Sensors - Neonatal.
0014925	Disposable SpO2 Sensors - Adult.
0014926	Disposable SpO2 Sensors - Paediatric.
0014927	Disposable SpO2 Sensors - Infant.
0014928	Disposable SpO2 Sensors - Neonatal.
0014929	Disposable SpO2 Sensors - Adult.
0014930	Disposable SpO2 Sensors - Paediatric.
0014933	Disposable SpO2 Sensors - Adult.
0014934	Disposable SpO2 Sensors - Paediatric.
0014935	Disposable SpO2 Sensors - Infant.
0014936	Disposable SpO2 Sensors - Neonatal.
0014937	Disposable SpO2 Sensors - Adult.
0014938	Disposable SpO2 Sensors - Paediatric.
0014939	Disposable SpO2 Sensors - Infant.
0014940	Disposable SpO2 Sensors - Neonatal.
0014941	Disposable SpO2 Sensors - Adult.
0014942	Disposable SpO2 Sensors - Paediatric.
0014945	Disposable SpO2 Sensors - Adult.
0014946	Disposable SpO2 Sensors - Paediatric.
0014947	Disposable SpO2 Sensors - Infant.
0014948	Disposable SpO2 Sensors - Neonatal.
0014949	Disposable SpO2 Sensors - Adult.
0014950	Disposable SpO2 Sensors - Paediatric.
0014951	Disposable SpO2 Sensors - Infant.
0014952	Disposable SpO2 Sensors - Neonatal.
0014953	Disposable SpO2 Sensors - Adult.
0014954	Disposable SpO2 Sensors - Paediatric.
0014957	Disposable SpO2 Sensors - Adult.
0014958	Disposable SpO2 Sensors - Paediatric.
0014959	Disposable SpO2 Sensors - Infant.
0014960	Disposable SpO2 Sensors - Neonatal.
0014961	Disposable SpO2 Sensors - Adult.
0014962	Disposable SpO2 Sensors - Paediatric.
0014963	Disposable SpO2 Sensors - Infant.
0014964	Disposable SpO2 Sensors - Neonatal.
0014965	Disposable SpO2 Sensors - Adult.
0014966	Disposable SpO2 Sensors - Paediatric.
0014969	Disposable SpO2 Sensors - Adult.
0014970	Disposable SpO2 Sensors - Paediatric.
0014971	Disposable SpO2 Sensors - Infant.

Stock Reference	Description
0014972	Disposable SpO2 Sensors - Neonatal.
0014973	Disposable SpO2 Sensors - Adult.
0014974	Disposable SpO2 Sensors - Paediatric.
0014975	Disposable SpO2 Sensors - Infant.
0014976	Disposable SpO2 Sensors - Neonatal.
0014977	Disposable SpO2 Sensors - Adult.
0014978	Disposable SpO2 Sensors - Paediatric.
0014981	Disposable SpO2 Sensors - Adult.
0014982	Disposable SpO2 Sensors - Paediatric.
0014983	Disposable SpO2 Sensors - Infant.
0014984	Disposable SpO2 Sensors - Neonatal.
0014985	Disposable SpO2 Sensors - Adult.
0014986	Disposable SpO2 Sensors - Paediatric.
0014987	Disposable SpO2 Sensors - Infant.
0014988	Disposable SpO2 Sensors - Neonatal.
0014989	Disposable SpO2 Sensors - Adult.
0014990	Disposable SpO2 Sensors - Paediatric.
0014993	Disposable SpO2 Sensors - Adult.
0014994	Disposable SpO2 Sensors - Paediatric.
0014995	Disposable SpO2 Sensors - Infant.
0014996	Disposable SpO2 Sensors - Neonatal.
0019001	Nellcor DEC-4 compatible SpO2 extension
0019002	Nellcor DEC-8 compatible SpO2 extension

Comments on Stock references review:

part numbers upto date

Supplier Review

Stock Ref.	Description	Supplier A/C	Supplier P/N	Supplier Name	Rating
0014740	SpO2 Silicone Finger Sen	00009299		Bluepoint Medical	A
0014790	SpO2 Silicone Finger Sen	00009299		Bluepoint Medical	A
0014901	Disposable SpO2 Sensors	00009299	1020131122	Bluepoint Medical	A
0014902	Disposable SpO2 Sensors	00009299	1020131125	Bluepoint Medical	A
0014903	Disposable SpO2 Sensors	00009299	1020131128	Bluepoint Medical	A
0014904	Disposable SpO2 Sensors	00009299	1020131130	Bluepoint Medical	A
0014905	Disposable SpO2 Sensors	00009299	1020131123	Bluepoint Medical	A
0014906	Disposable SpO2 Sensors	00009299	1020131126	Bluepoint Medical	A
0014909	Disposable SpO2 Sensors	00009299	1020131124	Bluepoint Medical	A
0014910	Disposable SpO2 Sensors	00009299	1020131127	Bluepoint Medical	A
0014911	Disposable SpO2 Sensors	00009299	1020131129	Bluepoint Medical	A
0014912	Disposable SpO2 Sensors	00009299	1020131131	Bluepoint Medical	A
0014925	Disposable SpO2 Sensors	00009299	1020131142	Bluepoint Medical	A
0014926	Disposable SpO2 Sensors	00009299	1020131145	Bluepoint Medical	A
0014927	Disposable SpO2 Sensors	00009299	1020131148	Bluepoint Medical	A
0014928	Disposable SpO2 Sensors	00009299	1020131150	Bluepoint Medical	A
0014957	Disposable SpO2 Sensors	00009299	1020131164	Bluepoint Medical	A
0014958	Disposable SpO2 Sensors	00009299	1020131167	Bluepoint Medical	A
0014959	Disposable SpO2 Sensors	00009299	1020131069	Bluepoint Medical	A
0014960	Disposable SpO2 Sensors	00009299	1020131071	Bluepoint Medical	A
0014590	SpO2 Silicone Finger Sen	00009300	302132142	Bluepoint Medical	A
0014591	SpO2 Silicone Finger Sen	00009300	3020132143	Bluepoint Medical	A
0014592	SpO2 Wrap Sensor- Paedia	00009300	5020132008	Bluepoint Medical	A
0014601	SpO2 Finger Sensor - Adu	00009300	1020132179	Bluepoint Medical	A
0014602	SpO2 Finger Sensor - Adu	00009300	1020132180	Bluepoint Medical	A
0014603	SpO2 Finger Sensor - Adu	00009300	1020132181	Bluepoint Medical	A
0014604	SpO2 Finger Sensor - Adu	00009300	1020132182	Bluepoint Medical	A
0014605	SpO2 Finger Sensor - Adu	00009300	1020132183	Bluepoint Medical	A
0014606	SpO2 Finger Sensor - Adu	00009300	1020132184	Bluepoint Medical	A
0014607	SpO2 Finger Sensor - Adu	00009300	1020132185	Bluepoint Medical	A
0014608	SpO2 Finger Sensor - Adu	00009300	1020132186	Bluepoint Medical	A
0014609	SpO2 Finger Sensor - Adu	00009300	1020132187	Bluepoint Medical	A
0014610	SpO2 Finger Sensor - Adu	00009300	1020132188	Bluepoint Medical	A
0014611	SpO2 Finger Sensor - Adu	00009300	1020132189	Bluepoint Medical	A
0014612	SpO2 Finger Sensor - Adu	00009300	1020132190	Bluepoint Medical	A
0014617	SpO2 Finger Sensor - Adu	00009300	1020132195	Bluepoint Medical	A
0014618	SpO2 Finger Sensor - Adu	00009300	1020132196	Bluepoint Medical	A
0014625	SpO2 Finger Sensor - Adu	00009300	4020132005	Bluepoint Medical	A
0014650	SpO2 Finger Sensor - Adu	00009300	1020132200	Bluepoint Medical	A
0014701	SpO2 Silicone Finger Sen	00009300	1020132223	Bluepoint Medical	A
0014702	SpO2 Silicone Finger Sen	00009300	1020132224	Bluepoint Medical	A
0014703	SpO2 Silicone Finger Sen	00009300	1020132225	Bluepoint Medical	A

Stock Ref.	Description	Supplier A/C	Supplier P/N	Supplier Name	Rating
0014704	SpO2 Silicone Finger Sen	00009300	1020132226	Bluepoint Medical	A
0014705	SpO2 Silicone Finger Sen	00009300	1020132227	Bluepoint Medical	A
0014706	SpO2 Silicone Finger Sen	00009300	1020132228	Bluepoint Medical	A
0014707	SpO2 Silicone Finger Sen	00009300	1020132229	Bluepoint Medical	A
0014708	SpO2 Silicone Finger Sen	00009300	1020132230	Bluepoint Medical	A
0014709	SpO2 Silicone Finger Sen	00009300	1020132231	Bluepoint Medical	A
0014710	SpO2 Silicone Finger Sen	00009300	1020132232	Bluepoint Medical	A
0014712	SpO2 Silicone Finger Sen	00009300	1020132234	Bluepoint Medical	A
0014713	SpO2 Silicone Finger Sen	00009300	1020132235	Bluepoint Medical	A
0014714	SpO2 Silicone Finger Sen	00009300	1020132236	Bluepoint Medical	A
0014715	SpO2 Silicone Finger Sen	00009300	1020132237	Bluepoint Medical	A
0014721	SpO2 Silicone Finger Sen	00009300	1020132243	Bluepoint Medical	A
0014725	SpO2 Silicone Finger Sen	00009300	4020132001	Bluepoint Medical	A
0014750	SpO2 Silicone Finger Sen	00009300	1020132245	Bluepoint Medical	A
0014751	SpO2 Silicone Finger Sen	00009300	1020132302	Bluepoint Medical	A
0014755	SpO2 Silicone Finger Sen	00009300	102132301	Bluepoint Medical	A
0014756	SpO2 Silicone Finger Sen	00009300	2020132102	Bluepoint Medical	A
0014761	SpO2 Silicone Finger Sen	00009300	2020132060	Bluepoint Medical	A
0014762	SpO2 Silicone Finger Sen	00009300	2020132061	Bluepoint Medical	A
0014763	SpO2 Silicone Finger Sen	00009300	2020132082	Bluepoint Medical	A
0014764	SpO2 Silicone Finger Sen	00009300	2020132083	Bluepoint Medical	A
0014766	SpO2 Silicone Finger Sen	00009300	2020132085	Bluepoint Medical	A
0014769	SpO2 Silicone Finger Sen	00009300	2020132088	Bluepoint Medical	A
0014770	SpO2 Silicone Finger Sen	00009300	2020132089	Bluepoint Medical	A
0014771	SpO2 Silicone Finger Sen	00009300	2020132090	Bluepoint Medical	A
0014774	SpO2 Silicone Finger Sen	00009300	2020132093	Bluepoint Medical	A
0014775	SpO2 Silicone Finger Sen	00009300		Bluepoint Medical	A
0014778	SpO2 Silicone Finger Sen	00009300	2020132097	Bluepoint Medical	A
0014785	SpO2 Silicone Finger Sen	00009300		Bluepoint Medical	A
0014790	SpO2 Silicone Finger Sen	00009300		Bluepoint Medical	A
0014801	Pulse Oximetry Ear Senso	00009300	1020132255	Bluepoint Medical	A
0014802	Pulse Oximetry Ear Senso	00009300	1020132256	Bluepoint Medical	A
0014804	Pulse Oximetry Ear Senso	00009300	102132258	Bluepoint Medical	A
0014805	Pulse Oximetry Ear Senso	00009300	1020132259	Bluepoint Medical	A
0014806	Pulse Oximetry Ear Senso	00009300	1020132251	Bluepoint Medical	A
0014807	Pulse Oximetry Ear Senso	00009300	1020132261	Bluepoint Medical	A
0014816	Pulse Oximetry Ear Senso	00009300	6020132001	Bluepoint Medical	A
0014820	Pulse Oximetry Wrap Sens	00009300		Bluepoint Medical	A
0014822	Pulse Oximetry Wrap Sens	00009300	2020132114	Bluepoint Medical	A
0014823	Pulse Oximetry Wrap Sens	00009300	202132115	Bluepoint Medical	A
0014824	Pulse Oximetry Wrap Sens	00009300	2020132123	Bluepoint Medical	A
0014825	Pulse Oximetry Ear Senso	00009300	1020132263	Bluepoint Medical	A
0014827	Pulse Oximetry Wrap Sens	00009300	2020132107	Bluepoint Medical	A
0014828	Pulse Oximetry Wrap Sens	00009300	2020132108	Bluepoint Medical	A
0014829	Pulse Oximetry Wrap Sens	00009300	2020132121	Bluepoint Medical	A

Stock Ref.	Description	Supplier A/C	Supplier P/N	Supplier Name	Rating
0014830	Pulse Oximetry Wrap Sens	00009300	2020132122	Bluepoint Medical	A
0014831	Pulse Oximetry Wrap Sens	00009300	0014031	Bluepoint Medical	A
0014835	Pulse Oximetry Wrap Sens	00009300	2020132008	Bluepoint Medical	A
0014848	Pulse Oximetry Wrap Sens	00009300		Bluepoint Medical	A
0014895	Pulse Oximetry Extension	00009300	1020132288	Bluepoint Medical	A
0014896	Pulse Oximetry Extension	00009300	1020132296	Bluepoint Medical	A
0014901	Disposable SpO2 Sensors	00009300	1020131122	Bluepoint Medical	A
0014902	Disposable SpO2 Sensors	00009300	1020131125	Bluepoint Medical	A
0014903	Disposable SpO2 Sensors	00009300	1020131128	Bluepoint Medical	A
0014904	Disposable SpO2 Sensors	00009300	1020131130	Bluepoint Medical	A
0014905	Disposable SpO2 Sensors	00009300	1020131123	Bluepoint Medical	A
0014906	Disposable SpO2 Sensors	00009300	1020131126	Bluepoint Medical	A
0014909	Disposable SpO2 Sensors	00009300	1020131124	Bluepoint Medical	A
0014910	Disposable SpO2 Sensors	00009300	1020131127	Bluepoint Medical	A
0014911	Disposable SpO2 Sensors	00009300	1020131129	Bluepoint Medical	A
0014912	Disposable SpO2 Sensors	00009300	1020131131	Bluepoint Medical	A
0014921	Disposable SpO2 Sensors	00009300	1020131134	Bluepoint Medical	A
0014922	Disposable SpO2 Sensors	00009300	1020131137	Bluepoint Medical	A
0014923	Disposable SpO2 Sensors	00009300	1020131139	Bluepoint Medical	A
0014924	Disposable SpO2 Sensors	00009300	1020131141	Bluepoint Medical	A
0014940	Disposable SpO2 Sensors	00009300	1020131160	Bluepoint Medical	A
0014957	Disposable SpO2 Sensors	00009300	1020131164	Bluepoint Medical	A
0014958	Disposable SpO2 Sensors	00009300	1020131167	Bluepoint Medical	A
0014959	Disposable SpO2 Sensors	00009300	1020131069	Bluepoint Medical	A
0014960	Disposable SpO2 Sensors	00009300	1020131071	Bluepoint Medical	A
0014964	Disposable SpO2 Sensors	00009300		Bluepoint Medical	A
0014981	Disposable SpO2 Sensors	00009300		Bluepoint Medical	A
0014982	Disposable SpO2 Sensors	00009300		Bluepoint Medical	A
0014983	Disposable SpO2 Sensors	00009300		Bluepoint Medical	A
0014984	Disposable SpO2 Sensors	00009300		Bluepoint Medical	A
0014990	Disposable SpO2 Sensors	00009300	1020131206	Bluepoint Medical	A
0014993	Disposable SpO2 Sensors	00009300	1020131204	Bluepoint Medical	A
0014994	Disposable SpO2 Sensors	00009300	1020131207	Bluepoint Medical	A
0014995	Disposable SpO2 Sensors	00009300	1020131209	Bluepoint Medical	A
0014996	Disposable SpO2 Sensors	00009300	1020131211	Bluepoint Medical	A
0019001	Nellcor DEC-4 compatible	00009300		Bluepoint Medical	A
0019002	Nellcor DEC-8 compatible	00009300		Bluepoint Medical	A
0031238	Connector kit Philips/H	00012399	KT-087	Unimed	B
0031256	Connector kit Nihon Koh	00012399	KT-084	Unimed	B

Comments on Suppliers:

supplier review upto date

Sales Information

Stock Reference	Description	2010	2011	2012	2013	2014	2015	2016
0014590	SpO2 Silicone Finger Sen	1						1
0014591	SpO2 Silicone Finger Sen			1				3
0014592	SpO2 Wrap Sensor- Paedia				1	11	6	1
0014595	SpO2 Silicone Finger Sensor - Adult.							
0014596	SpO2 Silicone Finger Sensor - Paediatric							
0014601	SpO2 Finger Sensor - Adu	14	82	6	20	25	4	7
0014602	SpO2 Finger Sensor - Adu	151	7	0	2	226	245	1
0014603	SpO2 Finger Sensor - Adu		2					
0014604	SpO2 Finger Sensor - Adu	4			11	10		3
0014605	SpO2 Finger Sensor - Adu	2	1	5	4			
0014606	SpO2 Finger Sensor - Adu	56	118	54	50	42	42	4
0014607	SpO2 Finger Sensor - Adu	100	82	21	7	6		0
0014608	SpO2 Finger Sensor - Adu	21	28	22	16	24	12	3
0014609	SpO2 Finger Sensor - Adu	6	27	5	5	12	6	
0014610	SpO2 Finger Sensor - Adu		11	8	6	6		
0014611	SpO2 Finger Sensor - Adu	2	9	3	7	1	28	
0014612	SpO2 Finger Sensor - Adu	3	49	14	11	18	5	2
0014613	SpO2 Finger Sensor - Adult.			3				
0014614	SpO2 Finger Sensor - Adult.							
0014615	SpO2 Finger Sensor - Adult.							
0014617	SpO2 Finger Sensor - Adu				3			
0014618	SpO2 Finger Sensor - Adu				5	15	5	4
0014619	SpO2 Finger Sensor - Adult.							
0014620	SpO2 Finger Sensor - Adult.							
0014621	SpO2 Finger Sensor - Adult.		10					
0014623	SpO2 Finger Sensor - Adult.				8			
0014625	SpO2 Finger Sensor - Adu		60	28	32	66	15	
0014649	SpO2 Finger Sensor - Adult.				45	5	1	
0014650	SpO2 Finger Sensor - Adu	8	195	228	202	165	189	82
0014701	SpO2 Silicone Finger Sen	25	48	17	2	14	1	
0014702	SpO2 Silicone Finger Sen	10	14	1	1	20	5	
0014703	SpO2 Silicone Finger Sen	21	7	18		1		
0014704	SpO2 Silicone Finger Sen	159	117	77	76	94	40	9
0014705	SpO2 Silicone Finger Sen			2		1		
0014706	SpO2 Silicone Finger Sen	6	50	13	4			0
0014707	SpO2 Silicone Finger Sen				15			0
0014708	SpO2 Silicone Finger Sen			1				
0014709	SpO2 Silicone Finger Sen		2				1	
0014710	SpO2 Silicone Finger Sen			2				
0014711	SpO2 Silicone Finger Sensor - Adult.					2	8	
0014712	SpO2 Silicone Finger Sen				10	1	4	
0014713	SpO2 Silicone Finger Sen					2		

Stock Reference	Description	2010	2011	2012	2013	2014	2015	2016
0014714	SpO2 Silicone Finger Sen							
0014715	SpO2 Silicone Finger Sen	91	66	12	37	23	24	
0014717	SpO2 Silicone Finger Sensor - Adult.							
0014718	SpO2 Silicone Finger Sensor - Adult.							
0014719	SpO2 Silicone Finger Sensor - Adult.							
0014720	SpO2 Silicone Finger Sensor - Adult.							
0014721	SpO2 Silicone Finger Sen							
0014722	SpO2 Silicone Finger Sensor - Adult.							
0014723	SpO2 Silicone Finger Sensor - Adult.	44	28	30	23	45	25	8
0014724	SpO2 Silicone Finger Sensor - Adult.	4		11	8	15	28	7
0014725	SpO2 Silicone Finger Sen	19	106	3	8	6	5	2
0014740	SpO2 Silicone Finger Sen							
0014750	SpO2 Silicone Finger Sen	342	422	409	654	503	305	143
0014751	SpO2 Silicone Finger Sen	93	74	47	51	43	66	18
0014755	SpO2 Silicone Finger Sen					2	3	
0014756	SpO2 Silicone Finger Sen							
0014761	SpO2 Silicone Finger Sen	11	13		22	11		
0014762	SpO2 Silicone Finger Sen	27	70	9	25	9	18	
0014763	SpO2 Silicone Finger Sen			1				
0014764	SpO2 Silicone Finger Sen		5		3	3		
0014765	SpO2 Silicone Finger Sensor - Paediatric							
0014766	SpO2 Silicone Finger Sen		1	5	2			
0014767	SpO2 Silicone Finger Sensor - Paediatric							
0014768	SpO2 Silicone Finger Sensor - Paediatric	1						
0014769	SpO2 Silicone Finger Sen		1					
0014770	SpO2 Silicone Finger Sen			1				
0014771	SpO2 Silicone Finger Sen		4					
0014772	SpO2 Silicone Finger Sensor - Paediatric							
0014773	SpO2 Silicone Finger Sensor - Paediatric							
0014774	SpO2 Silicone Finger Sen				0			
0014775	SpO2 Silicone Finger Sen	4			5			
0014777	SpO2 Silicone Finger Sensor - Paediatric							
0014778	SpO2 Silicone Finger Sen	4						
0014779	SpO2 Silicone Finger Sensor - Paediatric							
0014780	SpO2 Silicone Finger Sensor - Paediatric							
0014781	SpO2 Silicone Finger Sensor - Paediatric							
0014782	SpO2 Silicone Finger Sensor - Paediatric							
0014785	SpO2 Silicone Finger Sen		45	7	1	1	3	
0014790	SpO2 Silicone Finger Sen							
0014801	Pulse Oximetry Ear Senso	4	6	12	8	16	26	
0014802	Pulse Oximetry Ear Senso	14	11	9	28	71	54	-3
0014803	Pulse Oximetry Ear Sensor							
0014804	Pulse Oximetry Ear Senso	15	15	20	25	6		0
0014805	Pulse Oximetry Ear Senso	2						0
0014806	Pulse Oximetry Ear Senso							

Stock Reference	Description	2010	2011	2012	2013	2014	2015	2016
0014807	Pulse Oximetry Ear Senso			6	3	5	5	
0014808	Pulse Oximetry Ear Sensor							
0014816	Pulse Oximetry Ear Senso			18	17	16	12	1
0014820	Pulse Oximetry Wrap Sens							
0014821	Pulse Oximetry Wrap Sensor			4	15			
0014822	Pulse Oximetry Wrap Sens		3	5			0	
0014823	Pulse Oximetry Wrap Sens			4				
0014824	Pulse Oximetry Wrap Sens		7	7	6	10	1	
0014825	Pulse Oximetry Ear Senso	16	47	34	70	82	104	32
0014826	Pulse Oximetry Wrap Sensor							
0014827	Pulse Oximetry Wrap Sens							
0014828	Pulse Oximetry Wrap Sens	7	4	3				
0014829	Pulse Oximetry Wrap Sens	1			1		1	
0014830	Pulse Oximetry Wrap Sens		1					
0014831	Pulse Oximetry Wrap Sens						13	20
0014832	Pulse Oximetry Wrap Sensor							
0014833	Pulse Oximetry Wrap Sensor							
0014834	Pulse Oximetry Wrap Sensor							
0014835	Pulse Oximetry Wrap Sens	91	37	36	60	26	18	8
0014836	Pulse Oximetry Wrap Sensor							
0014837	Pulse Oximerty Wrap Sensor							
0014838	Pulse Oximetry Wrap Sensor							
0014840	Pulse Oximetry Wrap Sensor							
0014841	Pulse Oximetry Wrap Sensor							
0014842	Pulse Oximetry Wrap Sensor							
0014843	Pulse Oximetry Wrap Sensor							
0014844	Pulse Oximetry Wrap Sensor							
0014845	Pulse Oximetry Multi-Site Y Sensor							
0014848	Pulse Oximetry Wrap Sens			22				
0014895	Pulse Oximetry Extension	11	18	26	35	36	18	91
0014896	Pulse Oximetry Extension	6	4	3	3			
0014901	Disposable SpO2 Sensors	106	91	22	15	10	3	2
0014902	Disposable SpO2 Sensors	15	10	15	15	13	32	12
0014903	Disposable SpO2 Sensors	6	16			3		
0014904	Disposable SpO2 Sensors	84	56	18	15		20	26
0014905	Disposable SpO2 Sensors	30	54	49	42	45	24	8
0014906	Disposable SpO2 Sensors	295	34	49	5	15		
0014909	Disposable SpO2 Sensors	23	29		318	499	658	
0014910	Disposable SpO2 Sensors	24	14	5	13	35		2
0014911	Disposable SpO2 Sensors	31	54	30	232	84	8	5
0014912	Disposable SpO2 Sensors	29	152	206	66	98	47	53
0014913	Disposable SpO2 Sensors - Adult.							
0014914	Disposable SpO2 Sensors - Paediatric.							
0014915	Disposable SpO2 Sensors - Infant.							
0014916	Disposable SpO2 Sensors - Neonatal.							

Stock Reference	Description	2010	2011	2012	2013	2014	2015	2016
0014917	Disposable SpO2 Sensors - Adult.							
0014918	Disposable SpO2 Sensors - Paediatric.							
0014921	Disposable SpO2 Sensors					1		
0014922	Disposable SpO2 Sensors					1		
0014923	Disposable SpO2 Sensors							
0014924	Disposable SpO2 Sensors	25	1	2	3	2	2	
0014925	Disposable SpO2 Sensors	3		6				
0014926	Disposable SpO2 Sensors							
0014927	Disposable SpO2 Sensors							
0014928	Disposable SpO2 Sensors	2						
0014929	Disposable SpO2 Sensors - Adult.							
0014930	Disposable SpO2 Sensors - Paediatric.							
0014933	Disposable SpO2 Sensors - Adult.							
0014934	Disposable SpO2 Sensors - Paediatric.							
0014935	Disposable SpO2 Sensors - Infant.							
0014936	Disposable SpO2 Sensors - Neonatal.							
0014937	Disposable SpO2 Sensors - Adult.							
0014938	Disposable SpO2 Sensors - Paediatric.							
0014939	Disposable SpO2 Sensors - Infant.							
0014940	Disposable SpO2 Sensors			5				
0014941	Disposable SpO2 Sensors - Adult.							
0014942	Disposable SpO2 Sensors - Paediatric.							
0014945	Disposable SpO2 Sensors - Adult.							
0014946	Disposable SpO2 Sensors - Paediatric.							
0014947	Disposable SpO2 Sensors - Infant.							
0014948	Disposable SpO2 Sensors - Neonatal.							
0014949	Disposable SpO2 Sensors - Adult.							
0014950	Disposable SpO2 Sensors - Paediatric.							
0014951	Disposable SpO2 Sensors - Infant.							
0014952	Disposable SpO2 Sensors - Neonatal.							
0014953	Disposable SpO2 Sensors - Adult.							
0014954	Disposable SpO2 Sensors - Paediatric.							
0014957	Disposable SpO2 Sensors	3	8	4	12	51	11	2
0014958	Disposable SpO2 Sensors	10	14	18	24	55	27	8
0014959	Disposable SpO2 Sensors	7	12	28	45	90	48	22
0014960	Disposable SpO2 Sensors		5	24	25	29	13	12
0014961	Disposable SpO2 Sensors - Adult.							
0014962	Disposable SpO2 Sensors - Paediatric.							
0014963	Disposable SpO2 Sensors - Infant.							
0014964	Disposable SpO2 Sensors		9	3	2	4	2	2
0014965	Disposable SpO2 Sensors - Adult.							
0014966	Disposable SpO2 Sensors - Paediatric.							
0014969	Disposable SpO2 Sensors - Adult.							
0014970	Disposable SpO2 Sensors - Paediatric.							
0014971	Disposable SpO2 Sensors - Infant.							

Stock Reference	Description	2010	2011	2012	2013	2014	2015	2016
0014972	Disposable SpO2 Sensors - Neonatal.	2	1					
0014973	Disposable SpO2 Sensors - Adult.							
0014974	Disposable SpO2 Sensors - Paediatric.							
0014975	Disposable SpO2 Sensors - Infant.							
0014976	Disposable SpO2 Sensors - Neonatal.							
0014977	Disposable SpO2 Sensors - Adult.							
0014978	Disposable SpO2 Sensors - Paediatric.							
0014981	Disposable SpO2 Sensors						5	
0014982	Disposable SpO2 Sensors						11	12
0014983	Disposable SpO2 Sensors						30	30
0014984	Disposable SpO2 Sensors						9	
0014985	Disposable SpO2 Sensors - Adult.							
0014986	Disposable SpO2 Sensors - Paediatric.							
0014987	Disposable SpO2 Sensors - Infant.							
0014988	Disposable SpO2 Sensors - Neonatal.							
0014989	Disposable SpO2 Sensors - Adult.		1		7	11	2	
0014990	Disposable SpO2 Sensors		2		1		5	13
0014993	Disposable SpO2 Sensors		2			7	44	4
0014994	Disposable SpO2 Sensors	4	3		10	17	30	7
0014995	Disposable SpO2 Sensors	2	5	5	4	46	5	
0014996	Disposable SpO2 Sensors	8	22	28	22	4	8	
0019001	Nellcor DEC-4 compatible	9	39	18	6	5	8	5
0019002	Nellcor DEC-8 compatible		2	1	2	5		

Comments on Sales Information:

units still selling

Countries Review

Country	2010	2011	2012	2013	2014	2015	2016
AU Austria	[X]	[X]	[X]	[X]	[X]	[X]	[X]
B Belgium	[X]	[X]	[X]	[X]	[X]	[X]	[X]
BA State of Bahrain				[X]			
C China				[X]	[X]		
CAN Canada				[X]			
CH Chile			[X]				
CI Canary Islands		[X]					
CRO Croatia		[X]					
CY Cyprus		[X]		[X]			
CZ Czech Republic	[X]	[X]	[X]	[X]	[X]	[X]	[X]
DE Denmark	[X]	[X]	[X]	[X]	[X]	[X]	[X]
EG Egypt		[X]					
ES Estonia		[X]		[X]	[X]	[X]	
F France	[X]	[X]	[X]	[X]	[X]	[X]	[X]
FI Finland	[X]	[X]	[X]	[X]	[X]	[X]	
G Germany	[X]	[X]	[X]	[X]	[X]		
GR Greece	[X]	[X]	[X]	[X]	[X]	[X]	[X]
HK Hong Kong		[X]	[X]	[X]		[X]	
I Iran	[X]	[X]	[X]	[X]			
IN India	[X]	[X]	[X]				
IND Indonesia				[X]			
IQ IRAQ				[X]	[X]		
IRE Ireland	[X]	[X]	[X]	[X]	[X]	[X]	[X]
IS Israel	[X]		[X]		[X]	[X]	
JO Jordan			[X]				
KE Kenya	[X]		[X]			[X]	
KU Kuwait		[X]			[X]	[X]	[X]
LEB Lebanon	[X]		[X]				
LI Lithuania							[X]
MA Malaysia	[X]	[X]	[X]	[X]	[X]	[X]	
MAU Mauritius				[X]			
MEX Mexico					[X]		[X]
NE Netherlands	[X]		[X]			[X]	[X]
NG Nigera	[X]						
NO Norway	[X]	[X]	[X]	[X]	[X]	[X]	[X]
NZ New Zealand				[X]		[X]	
OM Oman	[X]		[X]	[X]	[X]	[X]	[X]
P Poland	[X]	[X]	[X]	[X]	[X]	[X]	[X]
PE Peru							[X]
PO Portugal	[X]	[X]	[X]	[X]	[X]	[X]	[X]
RO Romania	[X]	[X]	[X]	[X]		[X]	[X]
RU Russia				[X]		[X]	

Country	2010	2011	2012	2013	2014	2015	2016
SA South Africa		[X]	[X]				
SAU Saudi Arabia	[X]	[X]	[X]	[X]	[X]	[X]	[X]
SI Singapore			[X]				
SL Slovenia				[X]	[X]	[X]	
SP Spain	[X]	[X]	[X]	[X]	[X]	[X]	[X]
SRI Sri Lanka				[X]			
SW Sweden	[X]	[X]	[X]	[X]	[X]	[X]	[X]
SWI Switzerland		[X]		[X]	[X]	[X]	
SY Syria	[X]	[X]					
T Turkey	[X]	[X]	[X]	[X]	[X]	[X]	
TH Thailand		[X]	[X]	[X]	[X]	[X]	
UAE U A E	[X]	[X]	[X]	[X]	[X]	[X]	
UG UGANDA				[X]			
UK United Kingdom	[X]	[X]	[X]	[X]	[X]	[X]	[X]
VI Vietnam		[X]	[X]	[X]	[X]	[X]	

Comments on Sales to Countries:

no new countries

Comments on Risks with Sales to Countries:

no risks identified selling to countries

Returns Review

Stock Reference	Fault	2010	2011	2012	2013	2014	2015	2016
0014592	Unchecked - Returned to Stock							
0014601	Faulty probe		1					
0014601	Intermittent		1					
0014601	Low SpO2 readings		1					
0014601	Unchecked - Returned to Stock							
0014602	Faulty probe					9	2	
0014602	Unchecked - Returned to Stock							
0014605	Unchecked - Returned to Stock							
0014606	Connector damaged		1					
0014606	Damaged finger clip		1					
0014606	Faulty Electronics		1					
0014606	Faulty probe		2					
0014606	Intermittent		1	1				
0014606	Unchecked - Returned to Stock							
0014607	Faulty probe		1	2				
0014607	No Fault Return to Stock							
0014607	Unchecked - Returned to Stock							
0014608	Break at the connector					1		
0014608	Faulty probe			2				
0014608	Intermittant	1	1					
0014608	Intermittent	1		2				
0014608	Investigate		1					
0014608	No fault found		1	1				
0014608	Photo-diode intermittent in finger clip						1	
0014609	Unchecked - Returned to Stock							
0014610	Unchecked - Returned to Stock							
0014611	Unchecked - Returned to Stock							
0014612	Unchecked - Returned to Stock							
0014625	Cable replaced connectors rewired							1
0014625	Intermittent		1	1				
0014649	Faulty probe				3	3	1	
0014649	Sensor Fault				3	2		
0014650	Break at the connector					1		
0014650	Chip shorts to emitter				1			
0014650	Chip shorts to photodiode				4			
0014650	Emitter intermittent at the connector				2			
0014650	Faulty probe		3	19	13	5	30	
0014650	Intermittant						1	
0014650	Intermittent		2					
0014650	No Fault Found		1	2				
0014650	Photodiode intermittent at the connector				2			
0014650	Sensor Fault			4	7	12	17	2

Stock Reference	Fault	2010	2011	2012	2013	2014	2015	2016
0014701	Intermittent		1	1				
0014701	IR emitter failed						1	
0014701	No Fault Found			1				
0014701	Photo-diode intermittent in finger clip						1	
0014701	Red emitter faulty.			1				
0014701	Unchecked - Returned to Stock							
0014703	Faulty Electronics	1						
0014703	Faulty probe		1	1				
0014703	Silicone moulding split	1						
0014703	Unchecked - Returned to Stock							
0014704	Break at the connector			4				
0014704	Damaged finger clip							
0014704	Faulty Electronics		1					
0014704	Faulty probe			7				
0014704	Intermittant	1						
0014704	Intermittent		1	2				
0014704	No Fault Found		1					
0014704	No Fault Return to Stock							
0014704	Red emitter faulty.			1				
0014704	Unit dead			1				
0014706	Connector damaged		2					
0014706	Not working on monitor		1					
0014709	No Fault Found						1	
0014715	Damaged finger clip	3						
0014715	Faulty Electronics	1						
0014715	Intermittent	4	2					
0014715	Not working on monitor	2						
0014750	Break in sensor circuit						1	
0014750	Connector damaged	1		1				
0014750	Damaged finger clip		1					
0014750	Declined quote			2				
0014750	Faulty Electronics	1						
0014750	Faulty probe		1	4	8	1	6	
0014750	Intermittant				4			
0014750	Intermittent				1			
0014750	No fault found	1	1	3	3			
0014750	No Fault Return to Stock							
0014750	Not working on monitor	1						
0014750	Photodiode intermittent at the connector				2			
0014750	Sensor Fault			1	1		6	
0014750	Silicone moulding split					1		
0014750	Unable to repair	2						
0014750	Unchecked - Returned to Stock							
0014751	Faulty probe			1				
0014751	No Fault Found	1					1	

Stock Reference	Fault	2010	2011	2012	2013	2014	2015	2016
0014751	Silicone moulding split						2	
0014762	Unchecked - Returned to Stock							
0014774	Unchecked - Returned to Stock							
0014802	No Fault Found							6
0014804	Faulty electronics		1					
0014804	Unchecked - Returned to Stock							
0014822	Unchecked - Returned to Stock							
0014825	Faulty electronics				1			
0014825	No Fault Found		2	1	3	1		
0014825	Sensor fault				1	1	2	1
0014825	Unable to repair			1				
0014835	Intermittent	1			1			
0014835	No Fault Found				1			
0014835	Unchecked - Returned to Stock							
0014905	Unchecked - Returned to Stock							
0014909	No Fault Return to Stock							
0014910	No Fault Return to Stock							
0014911	No Fault Return to Stock							
0014911	Unchecked - Returned to Stock							
0014912	Unchecked - Returned to Stock							
0014924	Unchecked - Returned to Stock							
0014960	Unchecked - Returned to Stock							
0014995	Unchecked - Returned to Stock							

Comments on Returns:

low return rate

Comments on Risks with Returns and Potential Re-work:

no risks identified with returns

Design Changes Review

Showing Documents Filed in Y 14 Design Changes

Comments on Design Changes:

Design files held at OEM

Comments on Risks with Design Changes:

Design files held at OEM

User Instructions Review

Showing Documents Filed in F 5 User Instructions

Document ID	Description	Date Added/Updated
14707	4000 Series SpO2 Sensors Disposable Sensors User Instructions	29/01/15
14477	Pulse Oximeter Sensors BPM Instructions for Use / User Manual ...	24/10/14
14476	Pulse Oximeter Sensors BPM Instructions for Use / User Manual ...	24/10/14

Comments on User Instructions:

no updates to the IFU

Comments on Risks User Instructions:

no risk identified with IFU

Labels Review

Showing Documents Filed in F 7 / F 8 Labels

Document ID	Description	Date Added/Updated
14961	Pulse Oximetry Sensors 4000 Series Labels	29/04/15
11735	Pulse Oximetry Extension Cable 0014896 Dum Bell Label	27/03/12
11734	Pulse Oximetry Extension Cable 0014895 Dum Bell Label	27/03/12
6694	Pulse Oximetry Sensors 4000 Series Labels probe label	13/11/09
6693	Pulse Oximetry Sensors 4000 Series Labels pkg label	13/11/09

Comments on Labels:

Blue man icon added to labels

Comments on Risks Labels:

no risks identified with labeling

Documentation Updates / Changes

Document ID	Description	Date Added/Updated
15418	PPE Equipment	11/08/15
15139	Pulse Oximetry Sensors 4000 Series Competitors Information and P...	09/07/15
14961	Pulse Oximetry Sensors 4000 Series Labels	29/04/15
14948	Pulse Oximetry Sensors 4000 Series Promotional Leaflets 0014890 ...	23/04/15
14886	Capnograph Expected Life of Products Various Products Blue poin...	27/03/15
14857	Pulse Oximetry Sensors 4000 Post Market Surveillance	17/03/15

Comments on Document Changes:

labels updated

Comments on Risks with Document Changes:

no risks identified with documentation

Internal Issues Review

Number of Issues reviewed: 238

Issue ID	Subject
12302	
12301	
50228	0014602 failures from Yorkshire Ambulance
54353	0014623 GBP Distributor Price
57276	0014650 feedback - improved cabling assembly and strain relief
64115	0014802 cable length 2.0m
69433	0014904 lead times
63146	Backorder 00000019 0014802 ORD76220
63893	Backorder 00000019 0014802 ORD76702
64285	Backorder 00000333 0014996 ORD76900
57900	Backorder 00000637 0014625 ORD73447
57718	Backorder 00000900 0014618 ORD73360
64263	Backorder 00000970 0014835 ORD76865
60171	Backorder 00001661 0014606 ORD74889
59890	Backorder 00001661 0018720 ORD74707
58989	Backorder 00002020 0014924 ORD74178
64108	Backorder 00002020 0014924 ORD76819
60612	Backorder 00003007 0014650 ORD75075
67170	Backorder 00003220 0014750 ORD78554
67849	Backorder 00003220 0014750 ORD78870
58014	Backorder 00003310 0014612 ORD73552
57439	Backorder 00003836 0014609 ORD73235
60922	Backorder 00003902 0014606 ORD75172
60953	Backorder 00003902 0014606 ORD75187
61984	Backorder 00004651 0014958 ORD75676
69088	Backorder 00004729 0014990 ORD79439
65455	Backorder 00004729 0014994 ORD77530
68330	Backorder 00005030 0014608 ORD79071
58035	Backorder 00005110 0014725 ORD73574
66811	Backorder 00005110 0014725 ORD78339
59441	Backorder 00005430 0014611 ORD74397
61567	Backorder 00006039 0014750 ORD75476
67794	Backorder 00006089 0014750 ORD78793
62270	Backorder 00006144 0014909 ORD75834
67730	Backorder 00006277 0014601 ORD78814
68201	Backorder 00006277 0014601 ORD79042
68539	Backorder 00006277 0014601 ORD79190
63525	Backorder 00006277 0014904 ORD76435
65172	Backorder 00006277 0014904 ORD77357
66002	Backorder 00006277 0014904 ORD77910
67001	Backorder 00006277 0014904 ORD78487
68561	Backorder 00006277 0014904 ORD79190

Issue ID	Subject
68257	Backorder 00006297 0014590 ORD79055
66630	Backorder 00006326 0014750 ORD78188
60087	Backorder 00006772 0014829 ORD74831
61851	Backorder 00007025 0014750 ORD75604
60606	Backorder 00007025 0014755 ORD75058
65177	Backorder 00007088 0014750 ORD77358
65045	Backorder 00007246 0014751 ORD77247
65282	Backorder 00007327 0014650 ORD77420
62858	Backorder 00007356 0014905 ORD76084
63725	Backorder 00007356 0014905 ORD76573
66984	Backorder 00007356 0014905 ORD78479
65493	Backorder 00007356 0014911 ORD77557
60974	Backorder 00007374 0014608 ORD75141
62053	Backorder 00007374 0014611 ORD75710
59886	Backorder 00007374 0014704 ORD74712
63293	Backorder 00007374 0014704 ORD76263
68503	Backorder 00007374 0014704 ORD79164
60525	Backorder 00007374 0014712 ORD75018
59199	Backorder 00007374 0014715 ORD74264
60127	Backorder 00007374 0014715 ORD74844
64254	Backorder 00007374 0014715 ORD76884
64953	Backorder 00007374 0014715 ORD77218
58594	Backorder 00007374 0014724 ORD73917
63230	Backorder 00007374 0014724 ORD76263
63500	Backorder 00007374 0014724 ORD76402
64994	Backorder 00007374 0014724 ORD77218
65431	Backorder 00007374 0014724 ORD77517
66653	Backorder 00007374 0014724 ORD78194
69500	Backorder 00007374 0014724 ORD79561
58797	Backorder 00007374 0014750 ORD74065
60306	Backorder 00007374 0014750 ORD74933
69479	Backorder 00007374 0014802 ORD79585
63263	Backorder 00007374 0014824 ORD76270
58661	Backorder 00007374 0014957 ORD73970
60307	Backorder 00007374 0014958 ORD74934
61321	Backorder 00007374 0014958 ORD75397
68409	Backorder 00007374 0014958 ORD79138
69795	Backorder 00007374 0014958 ORD79737
58233	Backorder 00007374 0014959 ORD73668
62298	Backorder 00007374 0014959 ORD75873
63708	Backorder 00007374 0014959 ORD76548
66701	Backorder 00007374 0014959 ORD78225
67978	Backorder 00007374 0014959 ORD78931
69694	Backorder 00007374 0014959 ORD79675
57349	Backorder 00007374 0014981 ORD73109

Issue ID	Subject
58388	Backorder 00007374 0014981 ORD73797
57721	Backorder 00007374 0014982 ORD73109
58161	Backorder 00007374 0014982 ORD73568
58647	Backorder 00007374 0014982 ORD73970
66172	Backorder 00007374 0014982 ORD77958
58527	Backorder 00007374 0014983 ORD73797
59958	Backorder 00007374 0014983 ORD74771
65556	Backorder 00007374 0014983 ORD77546
66057	Backorder 00007374 0014983 ORD77849
66308	Backorder 00007374 0014983 ORD78029
66686	Backorder 00007374 0014983 ORD78225
67221	Backorder 00007374 0014983 ORD78587
57272	Backorder 00007374 0014984 ORD73109
58528	Backorder 00007374 0014984 ORD73797
58556	Backorder 00007374 0021013 ORD73864
60057	Backorder 00007415 0014750 ORD74820
63707	Backorder 00007415 0014750 ORD76559
62340	Backorder 00007415 0014995 ORD75892
57524	Backorder 00007515 0014650 ORD73283
58242	Backorder 00007515 0014650 ORD73608
59319	Backorder 00007515 0014650 ORD74271
61215	Backorder 00007515 0014650 ORD75085
61216	Backorder 00007515 0014650 ORD75147
62469	Backorder 00007515 0014650 ORD75781
63027	Backorder 00007515 0014650 ORD75899
63879	Backorder 00007515 0014650 ORD76057
66219	Backorder 00007515 0014650 ORD77234
66220	Backorder 00007515 0014650 ORD77292
66221	Backorder 00007515 0014650 ORD77825
66218	Backorder 00007515 0014650 ORD77826
58053	Backorder 00007515 0014750 ORD73472
66744	Backorder 00007515 0014750 ORD78273
62283	Backorder 00007515 0014825 ORD75863
59165	Backorder 00007515 0014994 ORD74244
66532	Backorder 00007515 0014994 ORD78148
59246	Backorder 00007530 0014704 ORD74301
64708	Backorder 00007530 0014723 ORD77085
56542	Backorder 00007530 0014902 ORD72577
57370	Backorder 00007530 0014902 ORD73168
60997	Backorder 00007530 0014902 ORD75211
65804	Backorder 00007530 0014902 ORD77746
66500	Backorder 00007530 0014902 ORD78118
67599	Backorder 00007530 0014902 ORD78772
68597	Backorder 00007542 0014825 ORD79166
60129	Backorder 00007542 0014989 ORD74867

Issue ID	Subject
66843	Backorder 00007542 0014989 ORD78385
63587	Backorder 00007586 0014650 ORD76351
68174	Backorder 00007592 0014751 ORD78983
58702	Backorder 00007688 0014751 ORD73965
61132	Backorder 00007688 0014988 ORD75284
61594	Backorder 00007688 0014996 ORD75284
60089	Backorder 00007759 0014822 ORD74833
62860	Backorder 00007759 0014831 ORD76093
62861	Backorder 00007759 0014831 ORD76094
66685	Backorder 00007759 0014831 ORD78224
68857	Backorder 00007759 0014831 ORD79312
59853	Backorder 00007828 0014750 ORD74690
60995	Backorder 00007828 0014909 ORD75114
63324	Backorder 00007883 0014723 ORD76334
68071	Backorder 00007885 0014912 ORD78968
68072	Backorder 00007885 0014912 ORD78969
68073	Backorder 00007885 0014912 ORD78970
63372	Backorder 00007954 0014750 ORD76344
69006	Backorder 00007995 0014750 ORD79413
59168	Backorder 00008001 0014602 ORD74173
66173	Backorder 00008016 0014910 ORD77995
60642	Backorder 00009044 0014750 ORD75023
61506	Backorder 00012369 0014592 ORD75451
57716	Backorder 00012580 0014993 ORD73362
66933	Backorder 00012672 0014990 ORD78438
68565	Backorder 00012672 0014990 ORD79204
69286	Backorder 00012672 0014990 ORD79521
69927	Backorder 00012672 0014990 ORD79772
64418	Backorder 00012742 0014994 ORD76974
66983	Backorder 00012906 0014904 ORD78477
69929	Backorder 00012923 0014750 ORD79668
69665	Backorder 00012923 4410520 ORD79668
66687	Backorder 00012949 0014592 ORD78231
69492	BCI compatible disposables
69491	BCI compatible disposables
69490	BCI compatible disposables
69489	BCI compatible disposables
69488	BCI compatible disposables
69487	BCI compatible disposables
69486	BCI compatible disposables
68928	BCI Disposables Sell
68931	Blank BCI Disposables Sell Potential Customers/Distributors
68929	Blank BCI Disposables Sell UK Selling Sites
63353	Bluepoint - Part No: 0014831
66871	Call Log #771 - Ohmeda SpO2 sensors

Issue ID	Subject
69530	Credit for SRS65286, Tunbridge Wells Hospital
66894	Credit required when goods returned
57284	Crown Agents updates
59692	CSI disposable SpO2 sensors for Medivent
59691	CSI disposable SpO2 sensors for Medivent
60865	Discontinue 800 series
66282	Disposable probes for disposal
13578	Distributor prices - VM2160 sensors
65925	DL3000 results for 0014802 HP compatible ear sensor
68491	Document Update Request (16438)
67486	Harley Health Village - returns and credits
68404	Investigation into 0014706 and 0014804
60660	Jamjoom PO: NMOG-636-0
60658	Jamjoom: Weights & Dimensions
60986	New account required
60985	New account required
50787	Ohmeda disposable 4000 series - suggest immediate discontinuation
63139	ORD76220 - Tunbridge Wells Hospital - On hold
63138	ORD76220 - Tunbridge Wells Hospital - On hold
60099	Please reserve 10x 0014909 for IASO (Greece)
57388	PMS Viamed OBL product Pulse Oximetry Sensors 4000 Series Supplier References
66775	POR10495 0014983 Have Been Delivered or Cancelled in Opera
68745	POR10554 0014983 Have Been Delivered or Cancelled in Opera
69605	POR10638 0014904 28/04/2016
69607	POR10638 0014904 29/04/2016
69606	POR10638 0014904 29/04/2016
14261	Provisional prices for Scandanavian company.
70151	Request Scanning Item to Sale or return
70150	Request Scanning Item to Sale or return
66845	Request Scanning Item to Sale or return
57852	Request Scanning Item to Sale or return
57851	Request Scanning Item to Sale or return
58501	Returned Items SOR457
59967	Returned Items SOR622
57766	Returned Items SOR631
58075	Returned Items SOR632
64547	Returned Items SOR652
64264	Returned Items SOR655
57855	Sale or return giving
61622	Sale or Return SOR655 - VM-2160 West Cumberland Hospital
13681	Single stock request
13680	Single stock request
14198	SpO2 probes JK Medical
64113	Viamed SMARTsat SpO2 sensors
64112	Viamed SMARTsat SpO2 sensors

Issue ID	Subject
63735	VM-2160 SmartSat Leaflet
69438	Watch List function to notify changes to shipping date
58190	Web Site Question ID 445
58567	Web Site Question ID 474
68832	WYA Consignment Stock (348)
68513	WYA Consignment Stock (348)
65830	WYA Consignment Stock (349)
64946	WYA Consignment Stock (349)
65788	WYA Consignment Stock reductions - stock returned
65794	Yorkshire Ambulance - invoice to credit
68398	Yorkshire Ambulance Consignment
68294	Yorkshire Ambulance Consignment Stock
65614	Yorkshire Ambulance Order
65613	Yorkshire Ambulance Order
65612	Yorkshire Ambulance Order
65611	Yorkshire Ambulance Order

Comments on Issues:

normal Issues orders / stocking, a limited number returns

Comments on Risks with Issues:

no risks identified with issues

Clinical / FDA Incidents online search

Clinical Investigation online review

Do any of the Results indicate a Risk / Problem : No
Do any of the Results indicate outdated Technology : No
Comments on Clinical Search :

No clinical reports on probes as they as a component of pulse oximetry.

Review of online FDA Incident reports

Do any of the Results indicate a Risk / Problem : No
Do any of the Results indicate outdated Technology : No
Comments on Clinical Search :

No FDA reports on probes as they as a component of pulse oximetry.