

Viamed Limited · 15 Station Road · Cross Hills Keighley · West Yorkshire · BD20 7DT · United Kingdom Tel: +44 (0)1535 634 542 · Fax: +44 (0)1535 635 582 Email: info@viamed.co.uk Website: www.viamed.co.uk

Service Repair Sheet 68942

Contact Name

Company/ Hospital Name

Department Position

Direct Phone General Phone

Opera Account Email

Order Number

Date Received Booked in By

Adam Keeton

Nottingham City Hospital EBME Department Clinical Engineering

01159691169 01159691169 00003910

adam.keeton@nhs.net

202006021 31/Oct/2024 Robert Connor

Repair

Ref S/N

SRN37180 5610152 AL000097

Equipment Type

Blender

Complete - Repaired Time: 0 Hour(s)

Parts Replaced

Qty: 1 0380050 Qty: 1 0332004

Blender AL000097 has been stripped and ultrasonically cleaned. It has been rebuilt, using the service kit, and tested as per the service manual.

Repair

Ref S/N SRN37181 5610152 AL000092 Equipment Type

Blender

Complete - Repaired Time :0 Hour(s)

Parts Replaced

0332004 Qty: 1 0380050

Blender AL000092 has been stripped and ultrasonically cleaned. It has been rebuilt, using the service kit, and tested as per the service manual.

Repair SRN37182 Ref

S/N 151109

Equipment Type Accessory

No Fault Found Time: 0 Hour(s)

Blender accessory not needed for service.