

Keighley · West Yorkshire · BD20 7DT · United Kingdom Viamed Limited to otation Tel: +44 (0)1535 634 542 · Fax: +44 (0)1535 635 582 Email: info@viamed.co.uk Website: www.viamed.co.uk

Service Repair Sheet 69012

Contact Name

Company/ Hospital Name

Direct Phone Opera Account

Email

Order Number **Date Received** Booked in By

Richard Wood

University Hospital Monklands

01698752074 00000070

Richard.Woods@lanarkshire.scot.nhs.uk

LMP15029679 29/Jan/2025 Robert Connor

Repair SRN37593

S/N

1410000 PR02894A14

Equipment Type

Foetal Simulator

Complete - Repaired

Time: 0 Hour(s)

Parts Replaced

1430309 Qty: 1 1480000 Qty: 1

V1000 has had transducer interface cushion replaced and calibration has been checked.