

Service Repair Sheet 69044

Contact Name	Mark Phipps
Company/ Hospital Name	Smit Mobile Equipment
Department	Mobile Medical Field Service
Position	Service Manager
Direct Phone	0127664490
General Phone	0127664490
Opera Account	00020025
Email	mark.phipps@smit.eu
Order Number	MP10606
Date Received	25/Feb/2025
Booked in By	Robert Connor

Repair	Ref	S/N	Equipment Type	
SRN37646	0111235	247391	Oxygen Monitor	Complete - Repaired Time :0.5 Hour(s)
Parts Replaced				
	Qty: 0.5	0140005		
	Qty: 1	0131215		

The MX300 Oxygen monitor monitor has been serviced and the sensor input connector has been replaced.

Repair	Ref	S/N	Equipment Type	
SRN37647	0111235	287917	Oxygen Monitor	Complete - Repaired Time :0.5 Hour(s)
Parts Replaced				
	Qty: 0.5	0140005		
	Qty: 1	0121217		

The MX300 Oxygen monitor monitor has been repaired. The case and sensor input connector have been replaced.

Repair	Ref	S/N	Equipment Type	
SRN37648	0111235	331755	Oxygen Monitor	Complete - Repaired Time :0.25 Hour(s)
Parts Replaced				
	Qty: 1	0121208		
	Qty: 0.25	0140005		

The MX300 Oxygen monitor monitor has been repaired. The keypad has been replaced.

Repair	Ref	S/N	Equipment Type	
SRN37649	0131202	0000003	Misc. Cable	End Of Life Time :0 Hour(s)
Parts Replaced				
	Qty: 1	0131202		

The Oxygen sensor cable was found to be open circuit and has been replaced with SN:001119.

Service Repair Sheet 69044

Contact Name	Mark Phipps
Company/ Hospital Name	Smit Mobile Equipment
Department	Mobile Medical Field Service
Position	Service Manager
Direct Phone	0127664490
General Phone	0127664490
Opera Account	00020025
Email	mark.phipps@smit.eu
Order Number	MP10606
Date Received	25/Feb/2025
Booked in By	Robert Connor

Repair	Ref	S/N	Equipment Type	
SRN37650	0131202	0004767	Misc. Cable	No Fault Found Time :0 Hour(s)

The Oxygen sensor cable has been tested. No fault was found.

Repair	Ref	S/N	Equipment Type	
SRN37651	0131202	0010575	Misc. Cable	No Fault Found Time :0 Hour(s)

The Oxygen sensor cable has been tested. No fault was found.
