














VIA MED

V1000 Foetal Heart Simulator

Class I



VIA MED

-  Stock References Review
-  Supplier Review
-  Sales Review
-  Countries Review
-  Returns Review
-  Design Changes Review
-  User Instructions Review
-  Labels Review
-  Document Updates Review
-  Internal Issues Review
-  Clinical / FDA Incidents online search

Stock References Review

Stock Reference	Description
0013110	Carrying case.
1410000	Foetal Heart Simulator model V1000.
1420000	Wall mount for Foetal Heart Simulator
1420001	Case protector for Foetal Heart
1430270	PCB - Foetal Heart Simulator.
1430284	V1000 Case with keypad cut outs.
1430285	V1000 Case.
1430286	V1000 Case Seal Kit.
1430287	V1000 Battery Compartment 4 x AA
1430288	V1000 Battery Clip 4 x AA
1430295	V1000 Keypad Membrane
1430296	V1000 label set (2) - case rear panel.
1430297	V1000 intensity control potentiometer.
1430300	V1000 PCB Assembly.
1430301	V1000 Transducer Mount.
1430302	V1000 Transducer.
1430303	V1000 Transducer mounting screw.
1430304	V1000 Transducer mounting washer.
1430305	V1000 Transducer mounting nut/spacer.
1430306	V1000 Transducer screw.
1430307	V1000 Transducer washer.
1430308	V1000 Transducer nut.
1430309	V1000 Transducer interface cushion.
1430319	V1000 Transducer interface cushion (v2).
1430320	V1000 Transducer v2.
1440000	Labour charge - Repair of V1000 Foetal
1450000	Packing box - V1000
1480000	Foetal Heart Simulator V1000 - Service
1480001	Foetal Heart Simulator V1000 -
1490000	Leaflet - V1000 Foetal heart simulator.
1490001	Instructions for Use - V1000 Foetal
1490002	V1000 Product Application Sheet

Comments on Stock references review:

Supplier Review

Stock Ref.	Description	Supplier A/C	Supplier P/N	Supplier Name	Rating
0013110	Carrying case.	00009445		Openhouse Products Ltd	B
1430270	PCB - Foetal Heart Simul	00010020		P & M Services	B
1430285	V1000 Case.	00011450		OKW Enclosures	B
1430286	V1000 Case Seal Kit.	00011450		OKW Enclosures	B
1490000	Leaflet - V1000 Foetal h	00009360		Print Centre	C
1430300	V1000 PCB Assembly.	00012189		Northern Hi-Tec - NHT	B
1430301	V1000 Transducer Mount.	00012189	V1000-2	Northern Hi-Tec - NHT	B
1420001	Case protector for Foeta	00011450		OKW Enclosures	B
1420000	Wall mount for Foetal He	00011450	A 92 26 227	OKW Enclosures	B
1430284	V1000 Case with keypad	00010665		Pendle Signs & Plastics	B
1430288	V1000 Battery Clip 4 x A	00011450		OKW Enclosures	B
1430308	V1000 Transducer nut.	00009314		Bradford Bolt & Nut Comp	B
1430306	V1000 Transducer screw.	00009314		Bradford Bolt & Nut Comp	B
1430307	V1000 Transducer washer.	00009314		Bradford Bolt & Nut Comp	B
1430287	V1000 Battery Compartment	00011450		OKW Enclosures	B
1430302	V1000 Transducer.	00012227	WM-R57A	Omni Circuits Internatio	B
1430305	V1000 Transducer mountin	00009061		Farnell	A
1430303	V1000 Transducer mountin	00009061		Farnell	A
1450000	Packing box - V1000	00011818		Macfarlane Group	B
1430296	V1000 label set (2) - ca	00010665		Pendle Signs & Plastics	B
1430297	V1000 intensity control	00009061		Farnell	A
1430297	V1000 intensity control	00009123		RS Components Ltd	A
1430295	V1000 Keypad Membrane	00012225		Danielson - Schurter	B
1430320	V1000 Transducer v2.	00009123		RS Components Ltd	A
1430304	V1000 Transducer mountin	00009061		Farnell	A
1490001	Instructions for Use - V	00009366		Warwick Printing	B
1490000	Leaflet - V1000 Foetal h	00009366		Warwick Printing	B
1430285	V1000 Case.				A
1430286	V1000 Case Seal Kit.				A
1430287	V1000 Battery Compartment				A
1430288	V1000 Battery Clip 4 x A				A
1430284	V1000 Case with keypad				A
1430296	V1000 label set (2) - ca				A
0013110	Carrying case.				A
1430300	V1000 PCB Assembly.				A
1420001	Case protector for Foeta				A
1430306	V1000 Transducer screw.				A
1430307	V1000 Transducer washer.				A
1430308	V1000 Transducer nut.				A
1430320	V1000 Transducer v2.	4248			Needs Grading
1490001	Instructions for Use - V				A
1430295	V1000 Keypad Membrane				A

Comments on Suppliers:

Sales Information

Stock Reference	Description	2019	2020	2021	2022	2023	2024	2025
0013110	Carrying case.	6	7			2		
1410000	Foetal Heart Simulator model V1000.	43	17	18		30		
1420000	Wall mount for Foetal He	1	1					
1420001	Case protector for Foeta	6	11			2		
1430270	PCB - Foetal Heart Simul							
1430284	V1000 Case with keypad				4	2		
1430285	V1000 Case.							
1430286	V1000 Case Seal Kit.				2			
1430287	V1000 Battery Compartment							
1430288	V1000 Battery Clip 4 x A				3			
1430295	V1000 Keypad Membrane	5			2			
1430296	V1000 label set (2) - ca				4			
1430297	V1000 intensity control				1			
1430300	V1000 PCB Assembly.	1				1		
1430301	V1000 Transducer Mount.							
1430302	V1000 Transducer.	1						
1430303	V1000 Transducer mountin							
1430304	V1000 Transducer mountin							
1430305	V1000 Transducer mountin							
1430306	V1000 Transducer screw.							
1430307	V1000 Transducer washer.							
1430308	V1000 Transducer nut.							
1430309	V1000 Transducer interface cushion.	23	12	12	30	27		
1430319	V1000 Transducer interface cushion (v2).							
1430320	V1000 Transducer v2.				1			
1440000	Labour charge - Repair of V1000 Foetal				3	2		
1450000	Packing box - V1000							
1480000	Foetal Heart Simulator V1000 - Service	26	12	11	30	27		
1480001	Foetal Heart Simulator V1000 -	20	1	8	12	4		
1490000	Leaflet - V1000 Foetal h							
1490001	Instructions for Use - V							
1490002	V1000 Product Application Sheet							

Comments on Sales Information:

part numbers upto date and correct

Countries Review

Country	2019	2020	2021	2022	2023	2024	2025
AUS Australia							
Australia							
Belgium	[X]	[X]					
C China							
Chile							
CO Colombia	[X]						
Colombia							
Czech Republic		[X]					
Finland							
France	[X]		[X]				
Germany	[X]				[X]		
Greece							
HK Hong Kong							
Hungary							
IN India							
India	[X]	[X]	[X]	[X]	[X]		
Indonesia	[X]	[X]					
Ireland					[X]		
IT Italy							
Kingdom of Bahrain							
LIB Libya							
MA Malaysia							
Malaysia	[X]	[X]	[X]	[X]			
NO Norway							
Northern Ireland							
Norway			[X]		[X]		
Oman							
PA Pakistan							
Pakistan	[X]						
Peru							
PH Philippines							
Philippines		[X]					
Poland					[X]		
Portugal	[X]						
Romania					[X]		
Saudi Aarabia	[X]						
Saudi Arabia	[X]						
Scotland	[X]	[X]	[X]	[X]			
Singapore					[X]		
South Africa	[X]		[X]				
SP Spain							
Switzerland							

Country	2019	2020	2021	2022	2023	2024	2025
Thailand					[X]		
UAE	[X]	[X]					
UAE U A E							
UK	[X]	[X]	[X]	[X]	[X]		
UK United Kingdom	[X]						
United Kingdom					[X]		
USA	[X]	[X]					
USA USA							

Comments on Sales to Countries:

no sales in last 18 months, just servicing

Comments on Risks with Sales to Countries:

Returns and Q.A. Fails Review

Stock Reference	Fault	2019	2020	2021	2022	2023	2024	2025
0013110	No Fault Found							
0013110	No Fault Return To Stock							
1410000	Needs Evaluating							1
1410000	Battery Contacts Replaced And Calibrated							
1410000	Broken Wire							
1410000	Calibration Checked	21	6	17	24	8	2	
1410000	Clean And Check							
1410000	Dead							
1410000	Faulty Board							
1410000	Membrane Keys Not Working							
1410000	New					2		7
1410000	No Fault Found	3						
1410000	Replaced Interface Cushion							
1410000	Replaced Keyboard Interface Cushion							
1410000	Replaced Keypad							
1410000	Replaced Keypad And Checked Calibration	2	2		1			
1410000	Replaced PCB And Calibrate							
1410000	Replaced PCB and Keypad then Calibrate	1						
1410000	Replaced Potentiometer							
1410000	Replaced Potentiometer and Serviced.						1	
1410000	Replaced Transducer And Serviced	1						
1410000	Replaced Transducer Interface Cushion And Calibrate							
1410000	Return to stock						1	
1410000	Rewired and checked calibration			1	1			
1410000	Serviced And Calibrated	1						
1410000	Serviced And Calibration Checked	21	34	31	40	54	87	40
1410000	Serviced with new battery case and contacts				1			
1410000	Serviced with new battery contacts and calibration checked.			1	3	1	1	3
1410000	Serviced with new battery contacts and pcb.					1		
1410000	Serviced With New Keypad	2	1	1	2		2	
1410000	Serviced with new lower case and calibration checked.				4	6		
1410000	Serviced With New Pot And Keypad							
1410000	Transducer Cushion Interface Replaced.							
1410000	Unchecked - Returned To Stock							
1410000	Unit Leaking							
1420001	No Fault Found							
1420001	No Fault Return To Stock							
1420001	Unchecked - Returned To Stock							
1430300	Needs Evaluating				2			

Comments on Returns:

nothing new in the returns file

Comments on Risks with Returns and Potential Re-work:

Design Changes Review

Showing Documents Filed in Y 14 Design Changes

Comments on Design Changes:

no changes

Comments on Risks with Design Changes:

User Instructions Review

Showing Documents Filed in F 5 User Instructions

Document ID	Description	Date Added/Updated
128797	V1000 Foetal Heart Simulator Instructions for Use User Manual Op...	08/09/23

Comments on User Instructions:

no changes

Comments on Risks User Instructions:

Labels Review

Showing Documents Filed in F 7 / F 8 Labels

Document ID	Description	Date Added/Updated
136745	V1000 Foetal Heart Simulator Labels	08/12/23

Comments on Labels:

no changes

Comments on Risks Labels:

Documentation Updates / Changes

Document ID	Description	Date Added/Updated
181887	V1000 Foetal Heart Simulator Training Information - Competency A...	04/04/25
181883	V1000 Foetal Heart Simulator Training Information - V1000 Foetal...	04/04/25
174487	GBP Price List - Technical Services - V1000 Foetal Heart Simulat...	30/01/25
174485	GBP Price List - V1000 Foetal Heart Simulator	30/01/25
174473	GBP Price List - Spare Parts - V1000 Foetal Heart Simulator	30/01/25
168868	V1000 Foetal Heart Simulator PAQ Questionnaire	26/11/24
150013	V1000 Foetal Heart Simulator Clinical Trials Reports Reviews and...	03/05/24

Comments on Document Changes:

Comments on Risks with Document Changes:

Internal Issues Review

Number of Issues reviewed: 23

Issue ID	Subject
320303	1410000 Price USD Distributor 0006 request Large Quantity
346946	Colsanitas: 1430309, 1430320
332477	Doto Health Pvt. Ltd: 1410000
339251	Evotec: V1000 & Carrying Case
351439	Fernastrqa Technomed Systems: V1000
328623	File 66 V1000 Foetal Heart Simulator Risk File Updated
342347	Juan Sebastian Pardo: V1000 (1410000)
359469	Kamil Dobransky: 1410000
348229	Machining of V1000 cases
348228	Machining of V1000 cases
348227	Machining of V1000 cases
346413	Pits Medical: V1000 1410000
334346	PVM3711 1430295 02/08/2024
347469	Request Scanning Item to Sale or return SOR968
83776	Returned Items SOR745
354148	Returned Items SOR968
343734	Sales forecasts (57)
110807	Sending Service Offers
352206	Ukusiza Abafundi: 1410000
331800	Umer Tariq: 1480000
331949	University Hospitals Plymouth: 1410000
343347	V1000 Cases to be drilled
339843	ZNS GmbH: 1410000

Comments on Issues:

just tracking service, and supplier type issues, nothing of concern

Comments on Risks with Issues:

Additional Sales Orders - Additional Sales Orders

1410000 Price USD Distributor 0006 request Large Quantity

Issue ID #320303

Date Created 09/Feb/2024

Issued To: Other

Completed Status: 08 Jan 2025

Notes:

Price tagged by Catrin Hird:

Add notes as to whom was asking for the price/quantity and Forward Issue to relevant sales person

Issues:320303

Projects - Projects

Colsanitas: 1430309, 1430320

Issue ID #346946

Date Created 14/Nov/2024

Issued To: General

Completed Status: 15 Nov 2024

Notes:

back and forth emails since 5.11.24 trying to decipher which spares he requires for V1000. Had help from Ryan, then when he was in Medica, SH and Phil advised me

Issues:346946

Projects - Projects

Doto Health Pvt. Ltd: 1410000

Issue ID #332477

Date Created 14/Jun/2024

Issued To: General

Completed Status: 14 Oct 2024

Notes:

price of 1410000, end user USD pricing given

Catrin Hird 27 Aug 2024

Neha Chandrakumar

Chased enquiry to see if still of interest.

Catrin Hird 16 Sep 2024

Neha Chandrakumar

Chased enquiry

Catrin Hird 30 Sep 2024

Neha Chandrakumar

Chased V1000 enquiry

Issues:332477

Projects - Projects

Evotec: V1000 & Carrying Case

Issue ID #339251

Date Created 27/Aug/2024

Issued To: General

Completed Status: 14 Oct 2024

Notes:

Confirmed pricing still valid and requested reply if require updated proforma invoice.

Catrin Hird 02 Sep 2024

Mariolize Van Der Mescht

Chased enquiry

Catrin Hird 16 Sep 2024

Mariolize Van Der Mescht

Chased enquiry

Catrin Hird 30 Sep 2024

Mariolize Van Der Mescht

Chased enquiry 1410000 V1000, provided leaflet and IFU.

Issues:339251

Projects - Projects

Fernastrqa Technomed Systems: V1000

Issue ID #351439

Date Created 07/Jan/2025

Issued To: General

Completed Status: Still Outstanding

Notes:

Interested in the V1000 through our website. Based in India, informed can not complete export orders through the website. Customer informed me that they will purchase for delivery within the UK and will arrange delivery to India themselves.

Customer requested cert of calibration with unit, informed of extra charge of 1480001.

Issues:351439

VIAMED Sales And Marketing - Post Market Surveillance

File 66 V1000 Foetal Heart Simulator Risk File Updated

Issue ID #328623

Date Created 03/May/2024

Issued To: Derek Lamb

Completed Status: 07 May 2024

Stock Ref: 0013110

Notes:

Risk file has changed, VM3COP27.12 Needs reviewing for 66 V1000 Foetal Heart Simulator as the Risk assesment has changed

03 May 2024 Derek Lamb
done

Issues:328623

Projects - Projects

Juan Sebastian Pardo: V1000 (1410000)

Issue ID #342347

Date Created 30/Sep/2024

Issued To: General

Completed Status: 17 Dec 2024

Notes:

Abandoned basket, based in Columbia. Emailed leaflet and instruction manual.

Catrin Hird 17 Dec 2024

Juan Sebastian Pardo

Customer requested proforma for V1000, awaiting delivery and invoice addresses.

Issues:342347

Projects - Projects

Kamil Dobransky: 1410000

Issue ID #359469

Date Created 24/Mar/2025

Issued To: General

Completed Status: Still Outstanding

Notes:

price of 1410000

Issues:359469

VIAMED Stock Meeting - Any Other Business

Machining of V1000 cases

Issue ID #348229

Date Created 28/Nov/2024

Issued To: Philip Crossley

Completed Status: 19 Mar 2025

Notes:

Copy Issue Main Issue #348227

Added by Steve Nixon sent to Helen Lamb

Please can a supplier account be set up for:

Maker CNC

Unit 1B

Gledholt Sidings

Huddersfield

HD1 4RR

Contact: Scott Dimelow

makercnchudd@gmail.com

I will then process an order for 50 x 1430284

V1000 cases machined with required holes for components.

Cathy - can you then send a copy of the order with the top half of 50 cases to Scott.

I will also provide you with a sample case (with potentiometer and keypad connector).

13 Feb 2025 Philip Crossley

These have been done and returned.

Issues:348229

VIAMED Stock Meeting - Any Other Business

Machining of V1000 cases

Issue ID #348228

Date Created 28/Nov/2024

Issued To: Catherine Spence

Completed Status: 19 Mar 2025

Notes:

Copy Issue Main Issue #348227

Added by Steve Nixon sent to Helen Lamb

Please can a supplier account be set up for:

Maker CNC

Unit 1B

Gledholt Sidings

Huddersfield

HD1 4RR

Contact: Scott Dimelow

makercnchudd@gmail.com

I will then process an order for 50 x 1430284

V1000 cases machined with required holes for components.

Cathy - can you then send a copy of the order with the top half of 50 cases to Scott.

I will also provide you with a sample case (with potentiometer and keypad connector).

09 Dec 2024 Catherine Spence

Done 1Z9W96386841781547

Issues:348228

VIAMED Stock Meeting - Any Other Business

Machining of V1000 cases

Issue ID #348227

Date Created 28/Nov/2024

Issued To: Helen Lamb

Completed Status: 19 Mar 2025

Notes:

Added by Steve Nixon sent to Helen Lamb
Please can a supplier account be set up for:
Maker CNC
Unit 1B
Gledholt Sidings
Huddersfield
HD1 4RR

Contact: Scott Dimelow
makercnchudd@gmail.com

I will then process an order for 50 x 1430284
V1000 cases machined with required holes for components.

Cathy - can you then send a copy of the order with the top half of 50 cases to Scott.

I will also provide you with a sample case (with potentiometer and keypad connector).

02 Dec 2024 Helen Lamb
account is on CID31323

20 Dec 2024 Helen Lamb

Issues:348227

Projects - Projects

Pits Medical: V1000 1410000

Issue ID #346413

Date Created 08/Nov/2024

Issued To: General

Completed Status: 23 Dec 2024

Notes:

Back in Stock Notification for V1000. Supplying MVM153281.

Catrin Hird 27 Nov 2024

Juan Pardo O`Larte

Chased Proforma MVM153281

Issues:346413

- Back Order Report

PVM3711 1430295 02/08/2024

Issue ID #334346

Date Created 04/Jul/2024

Issued To: Robert Connor

Completed Status: 04 Jul 2024

Notes:

PVM3711 1430295 Estimated ship to us date changed to 02 / 08 / 2024

Issues:334346

VIAMED Stock Meeting - Sale or Returns

Request Scanning Item to Sale or return SOR968

Issue ID #347469

Date Created 20/Nov/2024

Issued To: Catherine Spence

Completed Status: 26 Nov 2024

Notes:

Added by Sophie Lines sent to Catherine Spence

Please scan to Sale or Returns Location ID34132 Stock Reference:

1410000 V1000 x1

SOR968 Wrexham Maelor Hospital

26 Nov 2024 Catherine Spence

S/N:PR03665A14 id 486918

26 Nov 2024 Sophie Lines

Issues:347469

VIAMED Stock Meeting - Sale or Returns

Returned Items SOR745

Issue ID #83776

Date Created 17/Nov/2016

Issued To: Katie Evans

Completed Status: 02 Jun 2025

Notes:

Added by Jonathan Connor sent to Katie Evans

SOR ID745

Originally added by :Katie Evans

Royal Blackburn Hospital

Tahani Suliman

E.B.M.E. (Electronics Dept.)

Haslingden Road

Blackburn

Lancashire

BB2 3HH

Returned : 1410000 PR02194A10

18 Nov 2016 Katie Evans

Completed

24 Nov 2016 Jonathan Connor

Thanks

Issues:83776

VIAMED Stock Meeting - Sale or Returns

Returned Items SOR968

Issue ID #354148

Date Created 30/Jan/2025

Issued To: Catherine Spence

Completed Status: 17 Mar 2025

Notes:

Added by Catherine Spence sent to Catherine Spence

SOR ID968

Originally added by :Sophie Lines

Wrexham Maelor Hospital

Huw Bevan

Medical Engineering

Croesnewydd Road

Wrexham

LL13 7TD

Returned : 1410000 PR03665A14

Sold On CVM154707-1

13 Feb 2025 Catherine Spence

17 Mar 2025 Catherine Spence

done

Issues:354148

VIAMED Sales And Marketing - Post Market Surveillance

Sales forecasts (57)

Issue ID #343734

Date Created 14/Oct/2024

Issued To: Derek Lamb

Completed Status: 14 Oct 2024

Notes:

System Generated

Please fill in any outstanding sales forecasts,
Intrastats - Sales - Sales - Stock Sales Forecasts

Take a Look at the Sales Warnings.
Intrastats - Sales - Sales - Sales Warnings

14 Oct 2024 Derek Lamb

Parts with Continuous Sales That Recently Dropped or Stopped, No More Than 1 Month of No Sales:
These parts had consistent sales but have now seen a significant drop or stopped completely, without
having 2 or more months with no sales in the past:

0110057 - Oxygen sensor T-7:

Average: Â£462

Last month: Â£71 (85% drop)

0110058 - Oxygen sensor T-7V:

Average: Â£160

Last month: Â£32 (80% drop)

0110435 - Oxygen sensor MAX-250TM:

Average: Â£1,232

Last month: Â£0 (100% drop)

1410000 - Foetal Heart Simulator model V1000:

Average: Â£3,858

Last month: Â£659 (83% drop)

14 Oct 2024 Derek Lamb

Done

Issues:343734

Office Meeting - Office Procedures

Sending Service Offers

Issue ID #110807

Date Created 21/Dec/2017

Issued To: Derek Lamb

Completed Status: 02 Jun 2025

Linked Files

21 Dec 2017 V1000 Email Template

21 Dec 2017 V1000 Letter Template

Notes:

Added by Katie Evans sent to Katie Evans
Derek,

A brief overview of what I have been doing.
From lookup tab>Worklists?Chase Possible Service Items.

Clicking on 1410000 as TT letter not yet signed off, see issue 89789.

Click on get info, click on barcode, if the invoice number starts with CRE or it has an (r) then I haven't been sending as they are credits or returns (to my understanding) sent the letter as standard through GoldMine with template.

The GBP end user letter template is in Q>goldmine>Template>Katie>V1000 Service AMENDED.dot For GBP distributors I have simply amended the price.

The tricky part is separating the end users from the distributors.

For EUR/USD service offers, I made a very similar email template but used a table (particularly useful for distributors who may have bought loads) please find attached

To help me with the distributor purchases, you built me a table on this page:
<https://www.vmsecure.me.uk/intranet/operainfo/showstockaccount.php?zz=1&vui=57&user=Katie%20Evans&idp=raztfgzrBcdOg&submenu=&sessionid=717df487131d62116fdcdb61812981af&catagory=141%20%20Foetal%20Heart%20-%20Prod&stk=1410000&account=I0007909&name=METLOG%20Deutschland%20GmbH>

I hope this is useful. When I went through this with Emily, the parts she found to be too difficult were figuring out who is a distributor and figuring out which were credits or returns.

21 Dec 2017 Katie Evans

Next Action Changed From Katie Evans To Derek Lamb

27 Dec 2017 Katie Evans

Issues:110807

Projects - Projects

Ukusiza Abafundi: 1410000

Issue ID #352206

Date Created 10/Jan/2025

Issued To: General

Completed Status: Still Outstanding

Notes:

price of v1000 and carry pouch, checking lead time with SN for another customer too

Issues:352206

Projects - Projects

Umer Tariq: 1480000

Issue ID #331800

Date Created 07/Jun/2024

Issued To: General

Completed Status: 16 Sep 2024

Notes:

price of V1000 calibration. No srs given yet, he will email when one is required

Issues:331800

Projects - Projects

University Hospitals Plymouth: 1410000

Issue ID #331949

Date Created 10/Jun/2024

Issued To: General

Completed Status: 16 Sep 2024

Notes:

price of v1000

Issues:331949

VIAMED General Issue - VIAMED General Issue

V1000 Cases to be drilled

Issue ID #343347

Date Created 08/Oct/2024

Issued To: Steve Nixon

Completed Status: 08 Jan 2025

Notes:

Added by Catherine Spence sent to Steve Nixon

50 x 1430285 V1000 cased to be drilled, A new supplier was being sourced. Any luck

27 Oct 2024 Steve Nixon

Awaiting quotation from OKW.

Issues:343347

Projects - Projects

ZNS GmbH: 1410000

Issue ID #339843

Date Created 02/Sep/2024

Issued To: General

Completed Status: 18 Sep 2024

Notes:

price of 1410000

Catrin Hird 16 Sep 2024

Edwin Lian

Chased enquiry

Catrin Hird 18 Sep 2024

Edwin Lian

V1000 Information request: Project is temporarily discontinued due to the project termination.

Issues:339843

Clinical / FDA Incidents online search

Clinical Investigation online review

Do any of the Results indicate a Risk / Problem : No

Do any of the Results indicate outdated Technology : No

Comments on Clinical Search :

Review of online FDA Incident reports

Do any of the Results indicate a Risk / Problem : No

Do any of the Results indicate outdated Technology : No

Comments on Clinical Search :