



VIA MED

Tom Thumb

Hand held, gas powered resuscitation unit.

Class IIb
Via Rule Special Rule 9
Assesment Route Annex II
NBOG MD 1100

Carried out by Derek Lamb
02 / 06 / 25



VIA MED

-  Stock References Review
-  Supplier Review
-  Sales Review
-  Countries Review
-  Returns Review
-  Design Changes Review
-  User Instructions Review
-  Labels Review
-  Document Updates Review
-  Internal Issues Review
-  Clinical / FDA Incidents online search

Stock References Review

Stock Reference	Description
0310027	Tom Thumb Infant Resuscitator (Air).
0310029	Tom Thumb Infant Resuscitator (Air).
0310030	Tom Thumb Infant Resuscitator.
0310031	Tom Thumb Infant Resuscitator (Oxygen).
0310032	Tom Thumb Infant Resuscitator (Oxygen).
0310033	Tom Thumb Infant Resuscitator (Oxygen).
0310034	Tom Thumb Infant Resuscitator (Oxygen).
0310072	Tom Thumb Infant Resuscitator.
0310074	Tom Thumb Infant Resuscitator.
0310075	Tom Thumb Infant Resuscitator.
0310076	Tom Thumb Infant Resuscitator.
0310080	Tom Thumb Infant Resuscitator (RO).
0310082	Tom Thumb Infant Resuscitator (RO).
0310092	Tom Thumb Infant Resuscitator (TG).
0310093	Tom Thumb Infant Resuscitator.
0310094	Tom Thumb Infant Resuscitator.
0310096	Tom Thumb Oxygen Vent with Pole Mount
0310097	Tom Thumb Infant Resuscitator (TG).
0310098	Tom Thumb Oxygen Flow.
0310099	Tom Thumb Oxygen Vent.
0310250	Tom Thumb Conversion Kit.
0320102	Tom Thumb Rectangular Mounting Block.
0320103	Tom Thumb Rectangular Mounting Block.
0330185	Wing Nut Screw Black 25mm Thread
0330190	Universal clamp - 44mm horizontal
0330191	Clamp Spacer
0330193	Pressure gauge: 63mm dial -5 to +60 cmWG
0330200	TT460 MAIN Body Block
0330201	TT480/490 MAIN Body Block
0330202	TT495 MAIN Body Block
0330203	Pressure gauge for `Tom Thumb` Infant
0330204	Flowmeter (air) 0-15 Lpm
0330205	Flowmeter Meter Bolt
0330206	Blocking Bolt And Nut
0330207	Blanking Bolt
0330208	Inlet 10 mm
0330209	Outlet 15mm
0330210	Precision Valve
0330211	Adjustable Valve
0330212	Gauge Seal
0330213	O Ring Thick
0330214	O Ring Thin

Stock Reference	Description
0330215	Clamp Spacer
0330216	Screws For Clamp M5
0330217	Oxygen hose 1m. C/W indirect probe &
0330218	Oxygen Hose 3m. C/W indirect probe &
0330219	Right Angle For Hose
0330220	Fomblin OT20N Grease 100g tube
0330221	Lock Tight
0330222	Precision Valve Spring
0330223	Adjustable Valve Spring
0330224	Precision Valve Body
0330225	Precision/Adjustable Valve Seating
0330227	Precision/Adjustable Valve Seat
0330228	Precision Valve Adjustment Screw
0330229	Adjustable Valve Body
0330230	Adjustable Valve Screw
0330231	Adjustable Valve Internal Stop
0330232	Adjustable Valve Collar
0330233	Adjustable Valve Centre Screw
0330234	Adjustable Valve Side Screw
0330235	Inlet Restrictor
0330236	Replacement Gauge Face.
0330241	TT480/490 MAIN Body Block
0330242	Medirail Mounting Bracket
0330243	Adjustable Valve Side Screw. Bag 100
0330246	Blocking Bolt And Nut
0330247	Blanking Bolt
0330248	Inlet 10 mm
0330249	Outlet 15mm
0330250	Oxygen hose 1.5m. C/W DIN fitting:
0330251	Oxygen hose 3.0m. C/W DIN fitting:
0330263	Oxygen Hose 2m. C/W indirect probe &
0330268	Air Hose 3m. C/W indirect probe &
0330291	Inlet 10mm Thumbelina Inlet
0330292	Outlet 15mm Thumbelina
0330295	Outlet 10mm do not use
0330296	Blanking Bolt
0330297	Clamp Spacer - dove tailed.
0330300	Key Ring chain 1 inch ball
0330800	Label - pressure control for Tom Thumb
0340000	Labour - Tom Thumb infant resuscitator.
0340001	Training - Tom Thumb infant resuscitator
0340002	Training - Tom Thumb infant resuscitator
0350120	Packing box for Tom Thumb TT480.
0350121	Packing box for Tom Thumb TT490 Series.
0380000	Annual service - Tom Thumb

Stock Reference	Description
0380030	Tom Thumb Infant Resuscitator.
0390003	Tom Thumb Instruction Manual
0390004	Leaflet - Tom Thumb Infant Resuscitator.
0390005	Tom Thumb Training Manual
0390006	Tom Thumb Instruction Manual
0390007	Tom Thumb Instruction Manual
0390008	Tom Thumb Serial Number Label
0390009	`Tom Thumb` Label
0390011	Tom Thumb CE Label
0390015	Tom Thumb label
0390016	Address label for Tom Thumb gauge.
0390017	Gauge label - Tom Thumb
0390020	Tom Thumb Service Manual
0390021	Tom Thumb Service Manual TT460
0390022	Tom Thumb Service Manual TT480
0390023	Tom Thumb Service Manual TT490-15
0390024	Tom Thumb Service Manual TT496
0390060	Tom Thumb Instruction Manual - TT495

Comments on Stock references review:

Supplier Review

Stock Ref.	Description	Supplier A/C	Supplier P/N	Supplier Name	Rating
0330210	Precision Valve	00009117	9901-40	JNK Tooling	C
0330211	Adjustable Valve	00009117	9901-18	JNK Tooling	C
0330214	O Ring Thin	00009123		RS Components Ltd	A
0390004	Leaflet - Tom Thumb Infa	00009360		Print Centre	C
0330201	TT480/490 MAIN Body Bloc	00010160	0330201	KSJ Engineering Limited	B
0330216	Screws For Clamp M5	00009123	483-9616	RS Components Ltd	A
0330030	Dovetail (male)	00009150	4505	Therapy Equipment Ltd	B
0330203	Pressure gauge for `Tom	00010395		Wika Instruments Ltd	B
0330205	Flowmeter Meter Bolt	00010160	0330205	KSJ Engineering Limited	B
0330209	Outlet 15mm	00010160		KSJ Engineering Limited	B
0330215	Clamp Spacer	00010160	9901-22	KSJ Engineering Limited	B
0330212	Gauge Seal	00009150		Therapy Equipment Ltd	B
0330218	Oxygen Hose 3m. C/W indi	00009150		Therapy Equipment Ltd	B
0320060	Flowmeter (Oxygen) 0-15	00009150	6501	Therapy Equipment Ltd	B
0330220	Fomblin OT20N Grease 10	00009123		RS Components Ltd	A
0330221	Lock Tight	00009123	514-515	RS Components Ltd	A
0330219	Right Angle For Hose	00010770	2020	Chelmer Pneumatics Ltd	B
0330201	TT480/490 MAIN Body Bloc	00010725	0330201	MAJ Hi-Spec Ltd	B
0330205	Flowmeter Meter Bolt	00010725	0330205	MAJ Hi-Spec Ltd	B
0330209	Outlet 15mm	00010725		MAJ Hi-Spec Ltd	B
0330234	Adjustable Valve Side Sc	00009117	0330234	JNK Tooling	C
0330235	Inlet Restrictor	00009117	0330235	JNK Tooling	C
0330200	TT460 MAIN Body Block	00010160	0330200	KSJ Engineering Limited	B
0330202	TT495 MAIN Body Block	00010160	0330202	KSJ Engineering Limited	B
0330204	Flowmeter (air) 0-15 Lpm	00009150	6502	Therapy Equipment Ltd	B
0330200	TT460 MAIN Body Block	00010725	0330200	MAJ Hi-Spec Ltd	B
0330202	TT495 MAIN Body Block	00010725	0330202	MAJ Hi-Spec Ltd	B
0330055	Universal clamp - 44mm h	00009150	2002	Therapy Equipment Ltd	B
0330208	Inlet 10 mm	00010160		KSJ Engineering Limited	B
0330236	Replacement Gauge Face.	00010395		Wika Instruments Ltd	B
0330207	Blanking Bolt	00009123	0330207	RS Components Ltd	A
0330208	Inlet 10 mm	00010725		MAJ Hi-Spec Ltd	B
0330800	Label - pressure control	00009930		Renlim Labels	B
0330222	Precision Valve Spring	00012028	C6409	Harlow Springs Limited	C
0330223	Adjustable Valve Spring	00012028	9968	Harlow Springs Limited	C
0330217	Oxygen hose 1m. C/W indi	00009150		Therapy Equipment Ltd	B
0330213	O Ring Thick	00009123		RS Components Ltd	A
0330207	Blanking Bolt	00010725	0330207	MAJ Hi-Spec Ltd	B
0330206	Blocking Bolt And Nut	00010160		KSJ Engineering Limited	B
0390015	Tom Thumb label	00010665	TT IR	Pendle Signs & Plastics	B
0390016	Address label for Tom Th	00010665	TT ADR	Pendle Signs & Plastics	B
0390017	Gauge label - Tom Thumb	00010665	TT VIA	Pendle Signs & Plastics	B

Stock Ref.	Description	Supplier A/C	Supplier P/N	Supplier Name	Rating
0330216	Screws For Clamp M5	00009061	483-9616	Farnell	A
0330224	Precision Valve Body	00009117		JNK Tooling	C
0330228	Precision Valve Adjustme	00009117		JNK Tooling	C
0330227	Precision/Adjustable Val	00009117		JNK Tooling	C
0330295	Outlet 10mm do not use	00010160	0330295	KSJ Engineering Limited	B
0330300	Key Ring chain 1 inch ba	00012395	BC6	Seton	B
3210071	Single use silicone roun	00009721		Marshall Products Limite	B
3210011	NeoPEEP neonatal single	00009721		Marshall Products Limite	B
0120140	Adapter 15mm I.D. - 15m	00011844	PF 1515	Instrumentation Industri	B
0330250	Oxygen hose 1.5m. C/W DI	00009150		Therapy Equipment Ltd	B
0330251	Oxygen hose 3.0m. C/W DI	00009150		Therapy Equipment Ltd	B
0320103	Tom Thumb Rectangular Mo	00009022	PP8694	Tanda Engineering Ltd	A
0330190	Universal clamp - 44mm h	00009032		R & L Enterprises Ltd	B
0330212	Gauge Seal				A
3210011	NeoPEEP neonatal single	FL39B01			Needs Grading
0330268	Air Hose 3m. C/W indirec	00009150	ZZ HOSE 3M	Therapy Equipment Ltd	B
0330296	Blanking Bolt	00009123	367-6098	RS Components Ltd	A
0390004	Leaflet - Tom Thumb Infa	00009363		ACI Print Group Impress	B
0320221	Pole Clamp without slide	00007317	IHCPC1-SPECIAL	Inspiration Healthcare L	B
3210011	NeoPEEP neonatal single	00009602		Bunzl Healthcare	B
3210071	Single use silicone roun	00009602		Bunzl Healthcare	B
0330206	Blocking Bolt And Nut	00010725		MAJ Hi-Spec Ltd	B
0330243	Adjustable Valve Side Sc	00009073	020865070206020	Modern Screws	B
0390006	Tom Thumb Instruction Ma	00009345		Duffield Printers	A
0390007	Tom Thumb Instruction Ma	00009345		Duffield Printers	A
0350120	Packing box for Tom Thum	00009353	RPT10B	Rajapack Ltd	B
0350121	Packing box for Tom Thum	00009353	RPT17B	Rajapack Ltd	B
0330241	TT480/490 MAIN Body Bloc	00010160	9901-01	KSJ Engineering Limited	B
0330246	Blocking Bolt And Nut	00010160	9901-06	KSJ Engineering Limited	B
0330248	Inlet 10 mm	00010160	4701-61	KSJ Engineering Limited	B
0330249	Outlet 15mm	00010160	9901-34	KSJ Engineering Limited	B
0330234	Adjustable Valve Side Sc	00009073	0330234	Modern Screws	B
0330191	Clamp Spacer	00009032	A4917	R & L Enterprises Ltd	B
0330242	Medirail Mounting Bracke	00009032	A5090	R & L Enterprises Ltd	B
0330185	Wing Nut Screw Black 2	00009032	3335	R & L Enterprises Ltd	B
0330263	Oxygen Hose 2m. C/W indi	00009150	ZZ HOSE 2M	Therapy Equipment Ltd	B
0330230	Adjustable Valve Screw	00009117		JNK Tooling	C
0330193	Pressure gauge: 63mm dia	00009097	74656	GNW - Gauges North West	D
0390006	Tom Thumb Instruction Ma	00009366		Warwick Printing	B
0390003	Tom Thumb Instruction Ma	00009366		Warwick Printing	B
0330209	Outlet 15mm	00009032		R & L Enterprises Ltd	B
0330208	Inlet 10 mm	00009032		R & L Enterprises Ltd	B
0330206	Blocking Bolt And Nut	00009032		R & L Enterprises Ltd	B
0330224	Precision Valve Body	00009032		R & L Enterprises Ltd	B
0330228	Precision Valve Adjustme	00009032		R & L Enterprises Ltd	B

Stock Ref.	Description	Supplier A/C	Supplier P/N	Supplier Name	Rating
0330225	Precision/Adjustable Val	00009032		R & L Enterprises Ltd	B
0330227	Precision/Adjustable Val	00009032		R & L Enterprises Ltd	B
0330229	Adjustable Valve Body	00009032		R & L Enterprises Ltd	B
0330234	Adjustable Valve Side Sc	00009032	0330234	R & L Enterprises Ltd	B
0330233	Adjustable Valve Centre	00009032		R & L Enterprises Ltd	B
0330205	Flowmeter Meter Bolt	00009032	0330205	R & L Enterprises Ltd	B
0330230	Adjustable Valve Screw	00009032		R & L Enterprises Ltd	B
0330232	Adjustable Valve Collar	00009032		R & L Enterprises Ltd	B
0330231	Adjustable Valve Interna	00009032		R & L Enterprises Ltd	B
0330201	TT480/490 MAIN Body Bloc	00009032	0330201	R & L Enterprises Ltd	B
0330217	Oxygen hose 1m. C/W indi				A
0330203	Pressure gauge for `Tom				A
3210071	Single use silicone roun	FL39B01			Needs Grading
0330218	Oxygen Hose 3m. C/W indi				A
0330220	Fomblin OT20N Grease 10				A
0330218	Oxygen Hose 3m. C/W indi	00009211		MEC Medical Limited	B
0330236	Replacement Gauge Face.				A

Comments on Suppliers:

Sales Information

Stock Reference	Description	2019	2020	2021	2022	2023	2024	2025
0310027	Tom Thumb Infant Resuscitator (Air).							
0310029	Tom Thumb Infant Resuscitator (Air).							
0310030	Tom Thumb Infant Resuscitator.	109	6					
0310031	Tom Thumb Infant Resuscitator (Oxygen).							
0310032	Tom Thumb Infant Resuscitator (Oxygen).							
0310033	Tom Thumb Infant Resuscitator (Oxygen).							
0310034	Tom Thumb Infant Resuscitator (Oxygen).							
0310072	Tom Thumb Infant Resuscitator.							
0310074	Tom Thumb Infant Resuscitator.							
0310075	Tom Thumb Infant Resuscitator.							
0310076	Tom Thumb Infant Resuscitator.							
0310080	Tom Thumb Infant Resuscitator (RO).							
0310082	Tom Thumb Infant Resuscitator (RO).							
0310092	Tom Thumb Infant Resuscitator (TG).							
0310093	Tom Thumb Infant Resuscitator.							
0310094	Tom Thumb Infant Resuscitator.							
0310096	Tom Thumb Oxygen Vent with Pole Mount							
0310097	Tom Thumb Infant Resuscitator (TG).							
0310098	Tom Thumb Oxygen Flow.							
0310099	Tom Thumb Oxygen Vent.							
0310250	Tom Thumb Conversion Kit.	1		14				
0320102	Tom Thumb Rectangular Mounting Block.							
0320103	Tom Thumb Rectangular Mo							
0330185	Wing Nut Screw Black 2							
0330190	Universal clamp - 44mm h	100			56	50		
0330191	Clamp Spacer				6			
0330193	Pressure gauge: 63mm dia	2						
0330200	TT460 MAIN Body Block							
0330201	TT480/490 MAIN Body Bloc				6			
0330202	TT495 MAIN Body Block							
0330203	Pressure gauge for `Tom	24	14	30	11	9		
0330204	Flowmeter (air) 0-15 Lpm							
0330205	Flowmeter Meter Bolt							
0330206	Blocking Bolt And Nut				2			
0330207	Blanking Bolt							
0330208	Inlet 10 mm		1					
0330209	Outlet 15mm	15			2			
0330210	Precision Valve	4	1	2	5	4		
0330211	Adjustable Valve	13	1	11	2	3		
0330212	Gauge Seal	69	39	89	13	18		
0330213	O Ring Thick	77	44	75	125	40		
0330214	O Ring Thin	979	182	651	170	384		

Stock Reference	Description	2019	2020	2021	2022	2023	2024	2025
0330215	Clamp Spacer							
0330216	Screws For Clamp M5				2			
0330217	Oxygen hose 1m. C/W indi	5			9			
0330218	Oxygen Hose 3m. C/W indi	7						
0330219	Right Angle For Hose	10						
0330220	Fomblin OT20N Grease 10	1			2			
0330221	Lock Tight							
0330222	Precision Valve Spring							
0330223	Adjustable Valve Spring		6					
0330224	Precision Valve Body							
0330225	Precision/Adjustable Val							
0330227	Precision/Adjustable Val			2				
0330228	Precision Valve Adjustme							
0330229	Adjustable Valve Body							
0330230	Adjustable Valve Screw							
0330231	Adjustable Valve Interna							
0330232	Adjustable Valve Collar							
0330233	Adjustable Valve Centre							
0330234	Adjustable Valve Side Sc							
0330235	Inlet Restrictor							
0330236	Replacement Gauge Face.	6	1	6	4	3		
0330241	TT480/490 MAIN Body Bloc							
0330242	Medirail Mounting Bracke							
0330243	Adjustable Valve Side Sc							
0330246	Blocking Bolt And Nut							
0330247	Blanking Bolt							
0330248	Inlet 10 mm							
0330249	Outlet 15mm							
0330250	Oxygen hose 1.5m. C/W DI							
0330251	Oxygen hose 3.0m. C/W DI							
0330263	Oxygen Hose 2m. C/W indi	2						
0330268	Air Hose 3m. C/W indirec							
0330291	Inlet 10mm Thumbelina Inlet							
0330292	Outlet 15mm Thumbelina							
0330295	Outlet 10mm do not use							
0330296	Blanking Bolt				2			
0330297	Clamp Spacer - dove tailed.							
0330300	Key Ring chain 1 inch ba							
0330800	Label - pressure control							
0340000	Labour - Tom Thumb infant resuscitator.		1		1			
0340001	Training - Tom Thumb infant resuscitator							
0340002	Training - Tom Thumb infant resuscitator							
0350120	Packing box for Tom Thum							
0350121	Packing box for Tom Thum							
0380000	Annual service - Tom Thumb		8	24	15	9		

Stock Reference	Description	2019	2020	2021	2022	2023	2024	2025
0380030	Tom Thumb Infant Resuscitator.							
0390003	Tom Thumb Instruction Ma							
0390004	Leaflet - Tom Thumb Infa							
0390005	Tom Thumb Training Manual							
0390006	Tom Thumb Instruction Ma							
0390007	Tom Thumb Instruction Ma							
0390008	Tom Thumb Serial Number Label							
0390009	`Tom Thumb` Label							
0390011	Tom Thumb CE Label							
0390015	Tom Thumb label							
0390016	Address label for Tom Th							
0390017	Gauge label - Tom Thumb							
0390020	Tom Thumb Service Manual							
0390021	Tom Thumb Service Manual TT460							
0390022	Tom Thumb Service Manual TT480							
0390023	Tom Thumb Service Manual TT490-15							
0390024	Tom Thumb Service Manual TT496							
0390060	Tom Thumb Instruction Manual - TT495							

Comments on Sales Information:

no new sales, just spare parts

Countries Review

Country	2019	2020	2021	2022	2023	2024	2025
Belgium			[X]				
Canada							
China	[X]						
Czech Republic							
Finland							
France	[X]	[X]	[X]	[X]	[X]		
Germany	[X]						
IN India							
Italy		[X]					
Nigeria							
Portugal							
Russia		[X]					
Saudi Arabia							
Scotland	[X]	[X]	[X]	[X]			
Sri Lanka							
Switzerland							
Turkey							
UK	[X]	[X]	[X]	[X]	[X]		
UK United Kingdom	[X]			[X]			
United Kingdom	[X]						

Comments on Sales to Countries:

no new sales, just spare parts

Comments on Risks with Sales to Countries:

Returns and Q.A. Fails Review

Stock Reference	Fault	2019	2020	2021	2022	2023	2024	2025
0310030	Needs Evaluating							
0310030	Adjusted Factory Set Valves	2						
0310030	Calibrated Gauge							
0310030	Converted To TT480	2						
0310030	Converted to TT480 and service					4	5	
0310030	Customer Declined Quotation							
0310030	Damaged Pressure Gauge	1						
0310030	N/a							
0310030	New					2		2
0310030	No Fault Found	2						
0310030	No Fault Return To Stock							
0310030	Renovated			3	13	3		
0310030	Repaired Clamp And Serviced							
0310030	Replace Gauge And Service	2	11	2	4	6	2	
0310030	Replaced Adjustable Valve	1	1					
0310030	Replaced Gauge	2	1		1	1		1
0310030	Replaced Gauge and Adjustable Valve	2						
0310030	Replaced gauge face	1						
0310030	Replaced Gauge Face And Service	1	2	1	3	1	3	
0310030	Replaced precision valve and serviced				1			
0310030	Serviced	3	4		15	23	5	
0310030	Serviced / Adjusted Factory Set Valves	15	5	10	2	2	58	17
0310030	Serviced and Upgraded		1					
0310030	Stuck Gauge							
0310030	Tested					3		1
0310031	New					1		
0310031	Serviced							
0310032	Adjusted Factory Set Valves							
0310032	Calibrated Gauge							
0310032	Converted To TT480							
0310032	Damaged Pressure Gauge							
0310032	New							
0310032	Replaced Gauge Face And Service	1						
0310032	Serviced	1						
0310032	Serviced / Adjusted Factory Set Valves							
0310032	Serviced and Upgraded							
0310032	Tested							
0310033	Replace Gauge And Service							
0310033	Serviced							
0310034	Adjusted Factory Set Valves							
0310034	Converted To TT480							
0310034	Converted to TT480 and service						1	

Stock Reference	Fault	2019	2020	2021	2022	2023	2024	2025
0310034	Damaged Pressure Gauge							
0310034	Replace Gauge And Service	1						
0310034	Replaced hose and service.						3	1
0310034	Serviced		2	1				
0310034	Serviced / Adjusted Factory Set Valves	2		2				3
0310034	Unchecked - Returned To Stock							
0310076	New							
0310080	Replace Gauge And Service							
0310094	Converted To TT480							
0310094	Not Delivering More Than 10cm							
0310094	Replace Gauge And Service							
0310094	Replaced Gauge Face And Service				1			
0310097	Serviced		1		1	1		
0310097	Serviced / Adjusted Factory Set Valves	1					1	2
0310098	Serviced / Adjusted Factory Set Valves							
0310099	Serviced / Adjusted Factory Set Valves							
0330193	No Fault Found							
0330193	Unrepaired							
0330201	Needs Evaluating							
0330217	Unchecked - Returned To Stock							
0330218	Tested for leaks			4				
0380030	Needs Evaluating	1						
0380030	New			3	3			
0380030	Replaced gauge face					1		
0380030	Service exchange unit							
0380030	Serviced / Adjusted Factory Set Valves						3	3
0380030	Tested	3				3		5
0390023	Serviced							

Comments on Returns:

no new repair types, no issues with service

Comments on Risks with Returns and Potential Re-work:

Design Changes Review

Showing Documents Filed in Y 14 Design Changes

Comments on Design Changes:

no changes

Comments on Risks with Design Changes:

User Instructions Review

Showing Documents Filed in F 5 User Instructions

Document ID	Description	Date Added/Updated
15572	Tom Thumb Instructions for Use / User Manual TT490-15	01/09/15
15571	Tom Thumb Instructions for Use User Manual TT480	01/09/15
15570	Tom Thumb Instructions for Use User Manual 0310027	01/09/15
15511	Tom Thumb Instructions for Use User Manual TT480 Instructions D...	24/08/15
15510	Tom Thumb Instructions for Use User Manual TT480 Instructions D...	24/08/15
15509	Tom Thumb Instructions for Use User Manual TT480 Instructions D...	24/08/15
13047	Tom Thumb Instructions for Use User Manual French	14/10/13
12211	Tom Thumb Instructions for Use User Manual TT495-15	30/11/12
11732	Tom Thumb Instructions for Use User Manual	23/03/12
8512	Tom Thumb Instructions for Use User Manual NeoPEEP	06/07/11
8506	Instructions For Use Tom Thumb TT495	06/07/11
5107	Instructions For Use Tom Thumb TT480	23/02/09

Comments on User Instructions:

no changes

Comments on Risks User Instructions:

Labels Review

Showing Documents Filed in F 7 / F 8 Labels

Document ID	Description	Date Added/Updated
28357	Tom Thumb Labels Original labels as at Dec 2018	07/12/18
20970	Tom Thumb Labels Top Box Sticker example	17/07/17
17722	0330193 Pressure gauge: 63mm dial -5 to +60 cmWG (1416 / 1437)	19/10/16
11748	Tom Thumb Labels 0390016 Address Sticker Round outside Gauge	27/03/12
11747	Tom Thumb Labels Caution Label 0390018 Top Block	27/03/12
11746	Tom Thumb Labels 0390017 Front of Gauge Sticker	27/03/12
11745	Tom Thumb Labels 0390015 Front Block Logo Sticker	27/03/12
11744	Tom Thumb Labels 0390010 Adjustable Valve Sticker	27/03/12
11743	Tom Thumb Labels Part number to Label	27/03/12
11742	Tom Thumb Labels TT491 Serial	27/03/12
11741	Tom Thumb Labels TT490 Serial	27/03/12
11740	Tom Thumb Labels Serial TT480	27/03/12
7658	0310030 TT480 labelling modification 140211	14/02/11
5922	Tom Thumb Labels mediprema labels	14/05/09
2254	Tom Thumb Parts List Adj label 04	19/06/07
2204	Tom Thumb Labels Tomb Thumb Infant Resuscitator	19/06/07

Comments on Labels:

no changes

Comments on Risks Labels:

Documentation Updates / Changes

Document ID	Description	Date Added/Updated
174691	GBP Price List - 0310030 Tom Thumb TT480 Spare Parts	02/02/25
174689	GBP Price List - 0310034 Tom Thumb TT490-15-3 Spare Parts	02/02/25
174687	GBP Price List - 0310032 Tom Thumb TT490-15 Spare Parts	02/02/25
167103	VM3COP02.02 Viamed Company Responsibility's organisation chart s...	06/11/24
158002	WIKA ISO 9001:2015 Certificate	29/07/24
150015	Tom Thumb Post Market Surveillance Periodic Safety Update Repor...	03/05/24

Comments on Document Changes:

no issues with documents

Comments on Risks with Document Changes:

Internal Issues Review

Number of Issues reviewed: 27

Issue ID	Subject
107694	0310030 Tom Thumb Infant Resuscitator.
359589	Blank Intrastats 2025 Barcode Showing as Vandagraph
360286	Intelligent Control Co Ltd: 0330220
108577	Leaflet Tom Thumb Leaflet Amend 3rd Checking Stage
114403	Leaflet Update Infant Resus Cabinet Leaflet Leaflet Designed
125584	Photos/Images Tom Thumb for Cabinet Obtain Samples
329856	Product Failures New Codes Review (750)
155166	Proposed order from STL for 0310030
343695	PVM4043 0330190 31/10/2024
359068	PVM4352 0330218 19/03/2025
359763	Request Scanning Item to Sale or return
348580	Request Scanning Item to Sale or return
328809	Returned Items SOR965
341727	Returned Items SOR967
340674	Returned Items SOR967
364853	Returned Items SOR969
361443	Returned Items SOR971
364717	Returned Items SOR974
343515	Review The Tom Thumb Grease Date (1003)
359935	Sale Return SOR974 Leicester Royal Infirmary
107601	Tom Thumb Review 8th November 2017 - Discontinued Variants
352196	Tom Thumbs scanned to wrong location
330881	TomThumb needing stock ref change
329089	TomThumbs need a service form
113113	TT480 with pole mount
62275	User Manual Instructions and Labels Tom thumb Artwork to Intrastats
359229	Watch List Update 20 Mar 2025

Comments on Issues:

normal servicing type issues, nothing of concern

Comments on Risks with Issues:

VIAMED Stock Meeting - Stock Queries
0310030 Tom Thumb Infant Resuscitator.
Issue ID #107694
Date Created 09/Nov/2017
Issued To: Catrin Hird
Completed Status: 02 Jun 2025

Notes:

Copy Issue Main Issue #107693
Added by Derek Lamb sent to Steve Nixon
Design input specifications with mbar as test numbers.

worst case i.e. we set the values to the maximums as per the QA

Adjustable Value @ 10lpm 42.1 +/- 1 mbar
Setting on Test unit 43.0

Precision Value @ 10lpm Maximum 44.1 mbar +/- 1
Setting on the Test unit 45.0

At 5 Lpm

Adjustable 40.0
Precision Valve 42.5

At 10 Lpm

Adjustable 43.0
Precision Valve 45.0

At 15 Lpm

Adjustable 45.7
Precision Valve 47.2

20 Nov 2017 Catrin Hird
Do we have the recommended body weight range also?

20 Nov 2017 Catrin Hird

Urgent Flag Changed To On Also Gas Inlet flow range and operating time (400 litre cylinder) plus any additional notes to add to the specification screen. Thanks

20 Nov 2017 Catrin Hird
Added to marketing job 829 and issue reference 108587.

If you know the answer to the above 2 questions let me know otherwise I have included in issue to SN.

21 Nov 2017 Catrin Hird
Artwork in job 829, Sn to review and signed off for printing when suitable.

Issues:107694

Marketing Job Logger - Marketing Job Logger

Blank Intrastats 2025 Barcode Showing as Vandagraph

Issue ID #359589

Date Created 25/Mar/2025

Issued To: Derek Lamb

Completed Status: 02 Jun 2025

Notes:

Catrin Hird 25/03/2025 Issue Generated
Blank Intrastats 2025 Barcode Showing as Vandagraph

25 Mar 2025 Catrin Hird

Next Action Changed From To Derek Lamb Barcode 38100
Stock Ref 0310030

Showing as Vandagraph.

Showing on a few different ones, may be the stock ref which is doing it.

25 Mar 2025 Derek Lamb
seems its a historic mistag at some point

26 Mar 2025 Catrin Hird
OK
Done

Issues:359589

Projects - Projects

Intelligent Control Co Ltd: 0330220

Issue ID #360286

Date Created 01/Apr/2025

Issued To: General

Completed Status: Still Outstanding

Notes:

price of 0330220

Issues:360286

Marketing Job Logger - Marketing Job Logger

Leaflet Tom Thumb Leaflet Amend 3rd Checking Stage

Issue ID #108577

Date Created 20/Nov/2017

Issued To: Steve Nixon

Completed Status: 02 Jun 2025

Linked Files

20 Nov 2017 Tom Thumb Leaflet 2017-11-20

20 Nov 2017 Tom Thumb Leaflet 2017-11-20

20 Nov 2017 Tom Thumb Leaflet 2017-11-20 AMENDED

21 Nov 2017 Tom Thumb Leaflet 2017-11-21

Notes:

Catrin Hird 20/11/2017 Issue Generated
Leaflet Tom Thumb Leaflet Amend 3rd Checking Stage

20 Nov 2017 Catrin Hird

Next Action Changed From To Steve Nixon Please find attached revised artwork for Tom Thumb Leaflet removing relevant products.

Hard copy with proof sheet on SN desk.

20 Nov 2017 Catrin Hird

Updated reference to Configuration One on reverse of leaflet. Reattached.

21 Nov 2017 Catrin Hird

Revisions received, amended and provided back with fresh proof sheet.

Price list on issue 107601.

12 Mar 2018 Catrin Hird

0310094 to be added
Cabinet Tom Thumb (new part)

Cabinet Variant requires new image

Notes require clarification from 25/01/18 meeting above

26 Mar 2019 Steve Nixon

Product discontinued.

Issues:108577

Marketing Job Logger - Marketing Job Logger

Leaflet Update Infant Resus Cabinet Leaflet Leaflet Designed

Issue ID #114403

Date Created 13/Feb/2018

Issued To: Steve Hardaker

Completed Status: 02 Jun 2025

Notes:

Catrin Hird 13/02/2018 Issue Generated
Leaflet Update Infant Resus Cabinet Leaflet Leaflet Designed

13 Feb 2018 Catrin Hird

Next Action Changed From To Steve Hardaker
Priority Changed From 5 To 0 Hi,

I have been asked to amend the current infant resus cabinet leaflet so that it lists the individual parts that can be purchased...

As standard I've got:

Cabinet 0330105
Tom Thumb 0310030
Low suction controllers; oxytitre 0310140?
Therapy 7725 0310048 or 0310050
Cereatherm 0310302
Blender 0310200

But if we are offering a pick and mix sort of option then are there a large number of parts we can include?

Blenders - including precision or you happy to stick with Tom Thumb?

I guess we should also be including Apgar, AX300, Liners?

I think we have gone to look at this years ago but not sure where we got to.

15 Feb 2018 Steve Hardaker

Low suction controllers; oxytitre 0310140?
- This is not a standard unit, was used as a problem solver for 1 customer with no piped vacuum, offer 0310035 S714R with 3m hose instead.

Therapy 7725 0310048 or 0310050 - go with 0310050 which has a 3m hose but state shorter hoses available upon request.

Blenders - stick with Inspiration for now, I will discuss Precision Medical ones if the customer wants an alternative.

You can include the part numbers of the other main items, such as Apgar timer, Vacsax jars, liners, bed, mattress etc but leave off the AX300: we will offer it as a custom-requirement if it comes up but only 1 customer ever had them and took them all out again immediately.

The complete parts list can be found in the user instructions, you don't want to add absolutely everything as it goes down to spare parts, but the part numbers for the main items are all in there. Note: the version on L is out of date, the latest (v11.0) is attached to Issue 70173 for SN to upload.

16 Feb 2018 Catrin Hird
Fab thanks

19 Feb 2018 Steve Hardaker

Issues:114403

Marketing Job Logger - Marketing Job Logger

Photos/Images Tom Thumb for Cabinet Obtain Samples

Issue ID #125584

Date Created 06/Aug/2018

Issued To: Catherine Spence

Completed Status: 02 Jun 2025

Notes:

Catrin Hird 06/08/2018 Issue Generated
Photos/Images Tom Thumb for Cabinet Obtain Samples

06 Aug 2018 Catrin Hird

Next Action Changed From To Catherine Spence
Priority Changed From 5 To 0 Hi

Trying to tidy up my issues, I still require a Tom Thumb for Cabinet with a WIKA gauge 0310095 for photography. Thanks

I also require the pole clamp version 0310094 also with Wika gauge.

07 Aug 2018 Catherine Spence
i think we have a tom thumb in stock with the wika gauge, but i dont have the clamp

07 Aug 2018 Catrin Hird
I believe the 2 types of tom thumb as different, the cabinet version is different from the pole clamp version.
Do we have the clamps in stock or does SN need to order?

11 Dec 2018 Catrin Hird
Cabinet has been discontinued therefore no images of the tom thumb for cabinet to be taken.

11 Dec 2018 Catrin Hird
Cabinet has been discontinued therefore no images of the tom thumb for cabinet to be taken.

Issues:125584

VIAMED Management Meeting - VIAMED Management Meeting

Product Failures New Codes Review (750)

Issue ID #329856

Date Created 17/May/2024

Issued To: Derek Lamb

Completed Status: 21 May 2024

Notes:

System Generated

Task To be Completed

Review Intrastats Page:

Returns-> Repair codes Clean up

If New Faults have entered the system,

Review the Stock Repairs and apply a clean code and tag Fault / user or No Fault.

Complete all Stock References until the list is empty

17 May 2024 Derek Lamb

0110217 requires Cleaning! - SERVICE REPLACEMENT

0380030 requires Cleaning! - SERVICED / ADJUSTED FACTORY SET VALVES

2810004 requires Cleaning! - DISPLAY FAULT

7910110 requires Cleaning! - SENSOR LEAD FAILED

8000004 requires Cleaning! - NO FAULT FOUND

no capa required

Issues:329856

VIAMED Sales And Marketing - Sales Viamed Medical Export

Proposed order from STL for 0310030

Issue ID #155166

Date Created 25/Sep/2019

Issued To: Derek Lamb

Completed Status: 02 Jun 2025

Notes:

Copy Issue Main Issue #155165

Added by Steve Nixon sent to Ryan Swaine

Proposed order from STL for 0310030

Can you please confirm quantity (35) and timescale required.

Derek: If we have sufficient notes on Intrastats/Opera we could use the same part number. If so, we can use the existing labels, the CE mark is in an ideal position to be cut off.

or

We can generate a new part number to make it absolutely clear in terms of traceability.

Ryan have you registered the existing part number or would we be free to change it.

27 Sep 2019 Derek Lamb

see 155165

10 Feb 2020 Steve Nixon

I believe that they have ordered 6.

Suggest that you add link to sales projects. The main part to keep an eye on is the gauge, which we order from WIKA. They have been accepting orders for 50 rather than 100.

04 Mar 2020 Derek Lamb

job 1096

Issues:155166

- Back Order Report

PVM4043 0330190 31/10/2024

Issue ID #343695

Date Created 11/Oct/2024

Issued To: Robert Connor

Completed Status: 14 Oct 2024

Notes:

PVM4043 0330190 Estimated ship to us date changed to 31 / 10 / 2024

Issues:343695

- Back Order Report

PVM4352 0330218 19/03/2025

Issue ID #359068

Date Created 19/Mar/2025

Issued To: Robert Connor

Completed Status: 19 Mar 2025

Notes:

PVM4352 0330218 Estimated ship to us date changed to 19 / 03 / 2025

Issues:359068

VIAMED Stock Meeting - Sale or Returns

Request Scanning Item to Sale or return

Issue ID #359763

Date Created 26/Mar/2025

Issued To: Catherine Spence

Completed Status: 28 Mar 2025

Notes:

Added by Steve Hardaker sent to Catherine Spence

Please scan to Sale or Returns Location ID34132 Stock Reference:

3x 0310030 or 0380030 service loan Tom Thumbs for Leicester Royal. SOR974

27 Mar 2025 Catherine Spence

ID: 836909,836910,836911

SE004-SE006

28 Mar 2025 Steve Hardaker

28 Mar 2025 Steve Hardaker

Issues:359763

VIAMED Stock Meeting - Sale or Returns

Request Scanning Item to Sale or return

Issue ID #348580

Date Created 02/Dec/2024

Issued To: General

Completed Status: 08 Jan 2025

Linked Files

04 Dec 2024 Director Authorisation - SOR969 North Manchester 20241204

04 Dec 2024 SOR T&Cs Signed SOR969

04 Dec 2024 MIA Delivery Note - SOR969 20241204

04 Dec 2024 Letter - NMGH SOR969 Service Loan 20241204

04 Dec 2024 Email - Loan Tom Thumbs request SOR969

Notes:

Added by Steve Hardaker sent to Catherine Spence

Please scan to Sale or Returns Location ID34132 Stock Reference: 2x 0310030 service loan units.

I have asked Phil to check the calibration on the 2 in Main Office Stock Room 14, Row 1 Block 2 Shelf 2 Sale and Return / Service Exchange Stock.

03 Dec 2024 Catherine Spence

id 836912 and 836912 SE007, SE008

04 Dec 2024 Steve Hardaker

I think I have attached everything you need, please let me know if I've missed anything.

04 Dec 2024 Catherine Spence

1Z9W96386842060323

06 Dec 2024 Steve Hardaker

Next Action Changed From Catherine Spence To General Delivered On Friday, December 06 at 9:20 A.M. at Office

19 Dec 2024 Steve Hardaker

2x Tom Thumbs received from customer on SRS68974

Issues:348580

VIAMED Stock Meeting - Sale or Returns

Returned Items SOR965

Issue ID #328809

Date Created 07/May/2024

Issued To: Steve Hardaker

Completed Status: 07 May 2024

Notes:

Added by Robert Connor sent to Steve Hardaker
INFORMATION ONLY ISSUE DO NOT ADD NOTES!
SOR ID965

Originally added by :Steve Hardaker

Leicester Royal Infirmary

Roger Smith

Neonatal Unit

Infirmary Square

Leicester

LE1 5WW

Returned : 0310030 SE020

Returned : 0310030 SE022

Returned : 0380030 SE006

Issues:328809

VIAMED Stock Meeting - Sale or Returns

Returned Items SOR967

Issue ID #341727

Date Created 23/Sep/2024

Issued To: Steve Hardaker

Completed Status: 23 Sep 2024

Notes:

Added by Robert Connor sent to Steve Hardaker

INFORMATION ONLY ISSUE DO NOT ADD NOTES!

SOR ID967

Originally added by :Steve Hardaker

Royal Preston Hospital

Dave Edmondson

Medical Engineering

Sharoe Green Lane

Fulwood

Preston

PR2 9HT

Returned : 0310200 BX1651896

Returned : 0380030 SE008

Issues:341727

VIAMED Stock Meeting - Sale or Returns

Returned Items SOR967

Issue ID #340674

Date Created 11/Sep/2024

Issued To: Steve Hardaker

Completed Status: 11 Sep 2024

Notes:

Added by Robert Connor sent to Steve Hardaker

INFORMATION ONLY ISSUE DO NOT ADD NOTES!

SOR ID967

Originally added by :Steve Hardaker

Royal Preston Hospital

Dave Edmondson

Medical Engineering

Sharoe Green Lane

Fulwood

Preston

PR2 9HT

Returned : 0310200 BX1653066

Returned : 0380030 SE007

Issues:340674

VIAMED Stock Meeting - Sale or Returns

Returned Items SOR969

Issue ID #364853

Date Created 20/May/2025

Issued To: Steve Hardaker

Completed Status: 22 May 2025

Notes:

Added by Robert Connor sent to Steve Hardaker
INFORMATION ONLY ISSUE DO NOT ADD NOTES!
SOR ID969

Originally added by :Steve Hardaker
North Manchester General Hospital
Tim Holland
Clinical Engineering Dept
Delaunays Road

Manchester
M8 5RB

Returned : 0380030 SE007
Returned : 0380030 SE008

Issues:364853

VIAMED Stock Meeting - Sale or Returns

Returned Items SOR971

Issue ID #361443

Date Created 14/Apr/2025

Issued To: Steve Hardaker

Completed Status: 14 Apr 2025

Notes:

Added by Robert Connor sent to Steve Hardaker
INFORMATION ONLY ISSUE DO NOT ADD NOTES!
SOR ID971

Originally added by :Steve Hardaker
Darlington Memorial Hospital
David Watson
Clinical Engineering Department
Hollyhurst Road

Darlington
DL3 6HX

Returned : 0310030 SE020
Returned : 0310030 SE022

Issues:361443

VIAMED Stock Meeting - Sale or Returns

Returned Items SOR974

Issue ID #364717

Date Created 19/May/2025

Issued To: Steve Hardaker

Completed Status: 19 May 2025

Notes:

Added by Robert Connor sent to Steve Hardaker

INFORMATION ONLY ISSUE DO NOT ADD NOTES!

SOR ID974

Originally added by :Steve Hardaker

Leicester Royal Infirmary

Roger Smith

Neonatal Unit

Infirmary Square

Leicester

LE1 5WW

Returned : 0380030 SE004

Returned : 0380030 SE005

Returned : 0380030 SE006

Issues:364717

VIAMED Stock Meeting - VIAMED Stock Meeting
Review The Tom Thumb Grease Date (1003)

Issue ID #343515

Date Created 10/Oct/2024

Issued To: Helen Lamb

Completed Status: 25 Oct 2024

Notes:

System Generated

Task To be Completed

Check the grease used in the production and service of the Tom Thumb is with in it life span

Grease needs to be disposed of 4 years from the date of purchase.

10 Oct 2024 Helen Lamb

The Fomblin Expires 30/8/25 0330220 all good
Done

25 Oct 2024 Helen Lamb

Issues:343515

VIAMED Stock Meeting - Sale or Returns
Sale Return SOR974 Leicester Royal Infirmary

Issue ID #359935

Date Created 28/Mar/2025

Issued To: Steve Hardaker

Completed Status: Still Outstanding

Linked Files

28 Mar 2025 Director authorisation - Leicester Royal service loans

28 Mar 2025 SOR974 - Leicester Royal T&Cs 20250327

08 Apr 2025 322022.13872 Covering Letter - Leicester Royal SOR974 20250408

08 Apr 2025 MIA Agreement Delivery Note SOR974 20250408

08 Apr 2025 Checklist SOR974

Notes:

28 Mar 2025 Steve Hardaker Issue Generated from SOR

28 Mar 2025 Steve Hardaker

Next Action Changed From AUTO_SOR To Steve Hardaker Scanned SE004, 5 & 6 to SOR974.
Attached Director authorisation and signed T&Cs.

08 Apr 2025 Steve Hardaker

Next Action Changed From Steve Hardaker To Catherine Spence All paperwork attached, also printed out.
Dated for shipping on Wed 9th April, please add tracking details to this issue and close off when sent, thanks.

10 Apr 2025 Catherine Spence
1Z9W96386840215213

07 May 2025 Steve Hardaker

The services are done, just waiting for the revised PO to be invoiced.

07 May 2025 Steve Hardaker

Advised Silvia the tracking details, have issued SRS69096 For the return of 3x Tom Thumb TT480 service loan units.

09 May 2025 Steve Hardaker

Next Action Changed From Catherine Spence To Steve Hardaker Delivered 8/5/25. Awaiting return of loan units.

15 May 2025 Steve Hardaker

Followed up to ensure that they have received the 3 serviced units back and whether they have yet been able to ship the Service Loan units back to us on SRS69096. Asked for tracking details.

16 May 2025 Steve Hardaker

He took it to the post room yesterday, should be back next week.

19 May 2025 Steve Hardaker

Returned : 0380030 SE004

Returned : 0380030 SE005

Returned : 0380030 SE006

Closed off.

VIAMED Sales And Marketing - Sales

Tom Thumb Review 8th November 2017 - Discontinued Variants

Issue ID #107601

Date Created 08/Nov/2017

Issued To: Catrin Hird

Completed Status: 02 Jun 2025

Linked Files

20 Nov 2017 Infant Resuscitation Catalogue 2017-11-20

20 Nov 2017 Tom Thumb Price List 2017-11-20 DRAFT WEIGHTINGS ISSUE OUTSTANDING WITH DL

21 Nov 2017 Tom Thumb Price List 2017-11-21 DRAFT WEIGHTINGS ISSUE OUTSTANDING WITH DL

21 Nov 2017 Tom Thumb GBP Price List 2017-11-21

21 Nov 2017 Tom Thumb GBP Price List 2017-11-21 b

21 Nov 2017 Infant Resuscitation Catalogue 2017-11-21

21 Nov 2017 Infant Resuscitation Catalogue 2017-11-21 b

Notes:

Copy Issue Main Issue #107599

Added by Steve Nixon sent to Derek Lamb

As of today we will discontinuing non relevant Tom Thumb (Infant T-piece Resuscitator) devices, either due to using non blended gas, specials or low/no sales. See attached list for further details.

Variants retained & continued:

0310030 TT-480

Required for general sales and for use with resuscitation cabinets

0310094 TT-480 with pole mount

Continue, but monitor - low sales to date.

Catrin, please edit price lists, leaflet, catalogues & web site accordingly, i.e. just retain 0310030 TT-480

Once reviewed and appraised, may add 0310094.

20 Nov 2017 Catrin Hird

All variants of Tom thumb apart from 0310030, removed from website.

Infant Resus Catalogue - removed tom thumb variants and updated image - still needs work regarding acrylics being removed - attached current artwork for review by SN - i.e. page 9 has some white space left - would you like further information regarding specifications to be included (even though we are trying to steer away from that in catalogues)?!?

Amended pages 14/15 regarding "also available" pages.

Hard copy provided to SN with proofing sheet.

Price List

9901-22 shown on price list - is 0330191 correct replacement part? Have updated price list to generate from Intrastats - issue with DL regarding weightings not working (108584) but attached for reference.

The following prices as required to update other price list (if required)... Hard copy provided to SN with proofing sheet.

0310030 0005 International End User USD

0330191 0002 UK Distributor

0330191 0004 International Distributor Euros

0330191 0006 International Distributor USD

3210011 0003 International End User Euros

3210011 0005 International End User USD

3210013 0003 International End User Euros

3210013 0005 International End User USD
3210067 0002 UK Distributor
3210067 0003 International End User Euros
3210067 0005 International End User USD
3210067 0006 International Distributor USD
3210068 0002 UK Distributor
3210068 0003 International End User Euros
3210068 0004 International Distributor Euros
3210068 0005 International End User USD
3210068 0006 International Distributor USD
3210069 0002 UK Distributor
3210069 0003 International End User Euros
3210069 0004 International Distributor Euros
3210069 0005 International End User USD
3210069 0006 International Distributor USD
3210070 0003 International End User Euros
3210070 0005 International End User USD
3210071 0003 International End User Euros
3210071 0005 International End User USD
3210073 0003 International End User Euros
3210073 0005 International End User USD
3210074 0003 International End User Euros
3210074 0005 International End User USD
3210076 0003 International End User Euros
3210076 0005 International End User USD
3210077 0003 International End User Euros
3210077 0005 International End User USD

21 Nov 2017 Catrin Hird

9901-22/0330191 removed from price list as requested.

Incorrect part numbers given for silicone adapters, should have been 0120150 (instead of 0120149),
0120152 (instead of 0120159).

Updated prices required list...

0120151 0003 International End User Euros
0120151 0005 International End User USD
0120153 0002 UK Distributor
0120153 0003 International End User Euros
0120153 0004 International Distributor Euros
0120153 0005 International End User USD
0310030 0005 International End User USD
3210013 0003 International End User Euros
3210013 0005 International End User USD
3210023 0003 International End User Euros
3210023 0005 International End User USD
3210067 0002 UK Distributor
3210067 0003 International End User Euros
3210067 0005 International End User USD
3210067 0006 International Distributor USD
3210068 0002 UK Distributor
3210068 0003 International End User Euros
3210068 0004 International Distributor Euros
3210068 0005 International End User USD
3210068 0006 International Distributor USD
3210069 0002 UK Distributor
3210069 0003 International End User Euros
3210069 0004 International Distributor Euros
3210069 0005 International End User USD
3210069 0006 International Distributor USD
3210070 0003 International End User Euros
3210070 0005 International End User USD
3210071 0003 International End User Euros

3210071 0005 International End User USD
3210073 0003 International End User Euros
3210073 0005 International End User USD
3210074 0003 International End User Euros
3210074 0005 International End User USD
3210076 0003 International End User Euros
3210076 0005 International End User USD
3210077 0003 International End User Euros
3210077 0005 International End User USD

21 Nov 2017 Catrin Hird

Next Action Changed From Catrin Hird To Katie Evans Please review price list for correct descriptions.

Hard copy provided and saved in Z Checking, proof sheet provided.

21 Nov 2017 Catrin Hird
Updated and reattached.

21 Nov 2017 Catrin Hird
SN - Infant Resus catalogue provided to you for checking with proof sheet.

27 Nov 2017 Catrin Hird
Price list updated (still waiting for product weighing coding to be updated by DL #108584. Hard copy with proof sheet with KE.

28 Nov 2017 Catrin Hird
Created Related Issue #109225
Added by Catrin Hird sent to Steve Nixon
Price list updated and sent to SN for review with proof sheet. Pricing missing as follows:

0120151 0003 International End User Euros
0120151 0005 International End User USD
0120153 0002 UK Distributor
0120153 0003 International End User Euros
0120153 0004 International Distributor Euros
0120153 0005 International End User USD
0310030 0005 International End User USD
3210013 0003 International End User Euros
3210013 0005 International End User USD
3210023 0003 International End User Euros
3210023 0005 International End User USD
3210067 0002 UK Distributor
3210067 0003 International End User Euros
3210067 0005 International End User USD
3210067 0006 International Distributor USD
3210068 0002 UK Distributor
3210068 0003 International End User Euros
3210068 0004 International Distributor Euros
3210068 0005 International End User USD
3210068 0006 International Distributor USD
3210069 0002 UK Distributor
3210069 0003 International End User Euros
3210069 0004 International Distributor Euros
3210069 0005 International End User USD
3210069 0006 International Distributor USD
3210070 0003 International End User Euros
3210070 0005 International End User USD
3210071 0003 International End User Euros
3210071 0005 International End User USD
3210073 0003 International End User Euros
3210073 0005 International End User USD
3210074 0003 International End User Euros
3210074 0005 International End User USD

3210076 0003 International End User Euros
3210076 0005 International End User USD
3210077 0003 International End User Euros
3210077 0005 International End User USD

Also saved in: Z:Information for LibraryIn Checking with SN

28 Nov 2017 Catrin Hird

Next Action Changed From Katie Evans To Catrin Hird See related issue.

28 Nov 2017 Catrin Hird
see related issue

Issues:107601

VIAMED Stock Meeting - Missing Stock or Adjustments

Tom Thumbs scanned to wrong location

Issue ID #352196

Date Created 10/Jan/2025

Issued To: Catherine Spence

Completed Status: 13 Feb 2025

Notes:

Added by Steve Hardaker sent to Catherine Spence

Please can you look at part no 0310030, there are 6 units showing in a Viamed stock location but these are all customer units. It looks like they were incorrectly scanned:

- 6 Location ID 33278 Main Office Goods In corridor - Stock Ready for QA Row 3, 4 and 5 Shelves 2, 3 and 4

- VMDF01,2,5,&7 are South Tyneside District Hospital.

- VMDI25 is Royal Preston (note: it is the letter I but it was booked in for repair as VMD125).

- VMCJ48 is Hull Royal.

13 Feb 2025 Catherine Spence
derek is sorting for me

13 Feb 2025 Steve Hardaker

Issues:352196

VIAMED General Issue - VIAMED General Issue

TomThumb needing stock ref change

Issue ID #330881

Date Created 29/May/2024

Issued To: Derek Lamb

Completed Status: 30 May 2024

Notes:

Added by Philip Crossley sent to Derek Lamb

As discussed, I have a service exchange TomThumb on SRN36457. It is Serial number SE003, barcode 108869, part number 0310032.

Can you edit the part number to 0310030 as it is now a TT480.

Thanks.

29 May 2024 Derek Lamb
done

30 May 2024 Philip Crossley

Issues:330881

VIAMED General Issue - VIAMED General Issue

TomThumbs need a service form

Issue ID #329089

Date Created 09/May/2024

Issued To: Derek Lamb

Completed Status: 21 May 2024

Notes:

Added by Philip Crossley sent to Derek Lamb

I'm working on SRS68765 - 3 TomThumbs back from service exchange. I intend to service all 3 with new O rings ready for their next service exchange.

These TomThumbs have 2 different part numbers 0380030 and 0310030 and we have no service forms for any of them.

Could you clone the QA forms over to the service form please

20 May 2024 Derek Lamb
done

21 May 2024 Philip Crossley
Thanks

Issues:329089

VIAMED Stock Meeting - Photography

TT480 with pole mount

Issue ID #113113

Date Created 25/Jan/2018

Issued To: Jonathan Connor

Completed Status: 02 Jun 2025

Notes:

Added by Catrin Hird sent to Catrin Hird

25 Jan 2018 Catrin Hird

Next Action Changed From Catrin Hird To Jonathan Connor Please can I have a TT480 with pole clamp 0310094 for photography.

Thanky

19 Feb 2018 Catrin Hird

Header Changed From 332 VIAMED Stock Stock Queries To VIAMED Sales And Marketing Photography

Issues:113113

Marketing Job Logger - Marketing Job Logger

User Manual Instructions and Labels Tom thumb Artwork to Intrastats

Issue ID #62275

Date Created 01/Sep/2015

Issued To: Derek Lamb

Completed Status: 02 Jun 2025

Linked Files

01 Sep 2015 Tom Thumb TT490 Instructions for Use (high)

01 Sep 2015 Tom Thumb TT480 Instructions for Use (high)

01 Sep 2015 Tom Thumb TT480TG Instructions for Use (high)

Notes:

Catrin Hird 01/09/2015 Issue Generated

User Manual Instructions and Labels Tom thumb Artwork to Intrastats

01 Sep 2015 Catrin Hird

Next Action Changed From To Derek Lamb Please add attached to intrastats, print and email versions in Z for L

0310027

01 Sep 2015 Derek Lamb

Done

Issues:62275

Auto Heading Watch Lists - Auto Heading Watch Lists

Watch List Update 20 Mar 2025

Issue ID #359229

Date Created 20/Mar/2025

Issued To: Philip Crossley

Completed Status: 20 Mar 2025

Notes:

SupplierOrder PVM4352 0330218 Qty 5 Booked in on 20 Mar 2025

Issues:359229

Clinical / FDA Incidents online search

Clinical Investigation online review

Do any of the Results indicate a Risk / Problem : No

Do any of the Results indicate outdated Technology : No

Comments on Clinical Search :

Review of online FDA Incident reports

Do any of the Results indicate a Risk / Problem : No

Do any of the Results indicate outdated Technology : No

Comments on Clinical Search :