

MAXTEC

**REF** R125P02-011 MAX-250+  
**INTERNAL MEDICAL  
OXYGEN SENSOR**

Handle with care - Dropping can adversely affect performance.



Optimal Storage Temperature Range: 5 to 25° C  
Maximum Storage Temperature Range: -15 to 50° C

Operating Temperature Range: 5 to 40° C

90% Response Time @ 23° C: 15 seconds

Operation Life: 24 Months

Remove Protective Tab at time of use

MADE IN THE USA



(01)00853061006470(11)230530(21)JC68699052

\*JC68699052\*



**EC REP**

EMERGO EUROPE  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands



2023-05-30

**R<sub>x</sub> ONLY**

**CE** 0123

**MD**



TP-0057 (01)

**CH REP**

Medicel AG  
Dornierstrasse NR.11  
923 Altenrhein  
Switzerland

**UK RP**

Emergo Consulting (UK) Limited  
c/o Cr 360 - UL International,  
Compass House, Vision Park  
Histon, Cambridge CB24 9BZ, United Kingdom

**UK  
CA**

TP-0148 Rev 01