

Service Repair Sheet 69167

Contact Name	Matthew Jameson
Company/ Hospital Name	Sunderland Royal Hospital
Department	Clinical Engineering
Position	Senior Clinical Engineering Technician
Direct Phone	01915699933
General Phone	01915656256
Opera Account	00004992
Email	matthew.jameson@nhs.net
Order Number	RLNP400126294
Date Received	02/Jul/2025
Booked in By	Robert Connor

Repair	Ref	S/N	Equipment Type	
SRN37970	0310030	040026	Tom Thumb	Complete - Repaired Time :0 Hour(s)

Parts Replaced

Qty: 1	0380000
--------	---------

Tom Thumb TT480 040026 has been stripped and cleaned then reassembled with new O rings. The two valves have been checked and reset before QA.

Repair	Ref	S/N	Equipment Type	
SRN37971	0310030	040090	Tom Thumb	Complete - Repaired Time :0 Hour(s)

Parts Replaced

Qty: 1	0380000
Qty: 1	0330203
Qty: 1	0330212

Tom Thumb TT480 040090 has been stripped and cleaned then reassembled with new O rings and a new gauge. The two valves have been checked and reset before QA.

Repair	Ref	S/N	Equipment Type	
SRN37972	0310030	VMA113	Tom Thumb	Complete - Repaired Time :0 Hour(s)

Parts Replaced

Qty: 1	0380000
--------	---------

Tom Thumb TT480 VMA113 has been stripped and cleaned then reassembled with new O rings. The two valves have been checked and reset before QA.

Repair	Ref	S/N	Equipment Type	
SRN37973	0310030	VMA0084	Tom Thumb	Complete - Repaired Time :0 Hour(s)

Parts Replaced

Qty: 1	0380000
--------	---------

Tom Thumb TT480 VMA0084 has been stripped and cleaned then reassembled with new O rings. The two valves have been checked and reset before QA.

Service Repair Sheet 69167

Contact Name	Matthew Jameson
Company/ Hospital Name	Sunderland Royal Hospital
Department	Clinical Engineering
Position	Senior Clinical Engineering Technician
Direct Phone	01915699933
General Phone	01915656256
Opera Account	00004992
Email	matthew.jameson@nhs.net
Order Number	RLNP400126294
Date Received	02/Jul/2025
Booked in By	Robert Connor

Repair	Ref	S/N	Equipment Type	
SRN37974	0310030	VMA115	Tom Thumb	Complete - Repaired Time :0 Hour(s)
Parts Replaced				
	Qty: 1	0380000		

Tom Thumb TT480 VMA115 has been stripped and cleaned then reassembled with new O rings. The two valves have been checked and reset before QA.

Repair	Ref	S/N	Equipment Type	
SRN37975	0310030	VMA9001	Tom Thumb	Complete - Repaired Time :0 Hour(s)
Parts Replaced				
	Qty: 1	0380000		

Tom Thumb TT480 VMA9001 has been stripped and cleaned then reassembled with new O rings. The two valves have been checked and reset before QA.

Repair	Ref	S/N	Equipment Type	
SRN37976	0310030	VMC9072	Tom Thumb	Complete - Repaired Time :0 Hour(s)
Parts Replaced				
	Qty: 1	0380000		

Tom Thumb TT480 VMC9072 has been stripped and cleaned then reassembled with new O rings. The two valves have been checked and reset before QA.

Repair	Ref	S/N	Equipment Type	
SRN37977	0310030	VMA0103	Tom Thumb	Complete - Repaired Time :0 Hour(s)
Parts Replaced				
	Qty: 1	0380000		
	Qty: 1	0330236		

Tom Thumb TT480 VMA0103 has been stripped and cleaned then reassembled with new O rings and a new gauge face. The two valves have been checked and reset before QA.

Service Repair Sheet 69167

Contact Name	Matthew Jameson
Company/ Hospital Name	Sunderland Royal Hospital
Department	Clinical Engineering
Position	Senior Clinical Engineering Technician
Direct Phone	01915699933
General Phone	01915656256
Opera Account	00004992
Email	matthew.jameson@nhs.net
Order Number	RLNP400126294
Date Received	02/Jul/2025
Booked in By	Robert Connor

Repair	Ref	S/N	Equipment Type	
SRN37978	0310030	VMA0107	Tom Thumb	Complete - Repaired Time :0 Hour(s)
Parts Replaced				
	Qty: 1	0380000		

Tom Thumb TT480 VMA0107 has been stripped and cleaned then reassembled with new O rings. The two valves have been checked and reset before QA.

Repair	Ref	S/N	Equipment Type	
SRN37979	0310030	VMA0111	Tom Thumb	Complete - Repaired Time :0 Hour(s)
Parts Replaced				
	Qty: 1	0380000		
	Qty: 1	0330236		

Tom Thumb TT480 VMA0111 has been stripped and cleaned then reassembled with new O rings and a new gauge face. The two valves have been checked and reset before QA.

Repair	Ref	S/N	Equipment Type	
SRN37980	0310030	VMA0112	Tom Thumb	Complete - Repaired Time :0 Hour(s)
Parts Replaced				
	Qty: 1	0380000		

Tom Thumb TT480 VMA0112 has been stripped and cleaned then reassembled with new O rings. The two valves have been checked and reset before QA.
