

Photo-therapy Light Shield

COMPANY OPERATING PROCEDURE

Q.A. - Packing Checklist Procedure VM3/COP/52.03

Date: 29-Sep-03 Revision Date: 24-Jan-05	ssue: 1
--	---------

- 1. Unpack the Light shield from its original packing, and check that:
 - a. The product is the correct item
 - b. The product is the correct configuration:
 - i) PP Full Amber Hood
 - ii) PP Semi Hood
 - iii) PP Open Light Shield
 - iv) Special



- 2. Check the product to ensure that there are no sharp edges.
- 3. Check that there are no scratches on the product.
- 4. Check the neck aperture to ensure that there are no sharp edges.
- 5. Add the appropriate labels.
- 6. Add the insert leaflet
- 7. Clean the product and repack.
- 8. Ensure that there is an adequate amount of packaging material around the product prior to taping up the box.