

# How to choose the ideal Ventilator Tube Holder

1) **Choose the thickness**

Viamed supply two standard thickness's, 5mm and 10mm, as standard . Other thickness' are available by request.

5mm is lightweight, sits under the mattress easily but can be broken if treated harshly 10mm is heavy and is thick under the mattress but is hard to break

2) **Choose overall size**

Stipulate a height and length of base plate if different from the diagram measurements

3) **Hole positions and sizes**

It is possible to have any number of holes in any position in a variety of sizes.

The neck joining the hole to the outside edge is 5mm

4) **To choose a hole**

Choose a position and required size. Mark the table accordingly

e.g. E Yes                      **22mm**

Only holes marked on the Hole selection can will be drilled

5) **Centre holes**

If a hole is required in the centre of the vertical plate please draw a circle in position on the diagram with measurements in millimetres from each edge and state the hole size.

6) **Non standard Ventilator Tube Holders**

If you are not happy with the standard versions offered please submit a diagram of your requirements and we will give you a quotation

Please specify in millimetres

Thickness

Height of vertical plate

Length of base plate

Position of all holes

Size of Holes

