



Main Account <viamedinbox@gmail.com>

Reference Guide - Maxell KE Series Lead-Free Oxygen Sensors

1 message

Steve Nixon <steve.nixon@viamed.co.uk>
To: Cameron.Billingham@spectrummedical.com
Bcc: kate.griffiths@viamed.co.uk

22 April 2026 at 15:12

Dear Cameron

Thank you for your order for 100 x 0110545 KE-25LF sensors (measurement range 0 - 30% oxygen, output 10.0 mV - 15.5 mV).

I have attached a quick reference guide in case you require any other sensor variants in the future. There is the KE-25F3LFM which has the same range (0 - 100% oxygen) of the original KE25 (leaded) sensor, but the output is much lower at 5.0 - 8.0 mV, and it has a threaded M16 input connection.


Kind regards

--
Steve

Steve Nixon
Director - Viamed Ltd.

Office: +44 (0)1535 634542
Mobile: +44 (0)7850 252267

GDPR / Privacy Policy available at www.viamed.co.uk/policy

 **Reference Guide - Maxell KE Series Lead-Free Oxygen Sensors ID153212 (2).pdf**
661K