




VIAMED



TOM THUMB INFANT RESUSCITATOR

TOM THUMB INFANT RESUSCITATOR DATA SHEET



Please specify MAX Pressure for Safety blow-off valve PT. NO. 9901-17

50 [] 40 [] 30 [] 20 [] Other []

TO ORDER

Choose layout of components and number the circle with the hole number.

Limitation:

Flowmeter PT. NO. 6575 must be vertical in port '6'.

Adjustable Blow-off valve must be vertical in the position shown above.

Safety Blow-off valve PT. NO. 9901-17 must be in a vertical position on the bottom face of the block.

Contact	Order Number:
Title	Invoice Number:
Department	Serial Number:
Hospital	Q/A
Tel. No.	Date:

Additional Information:

Viamed Ltd. 15 Station Road, Cross Hills, Keighley, West Yorkshire BD20 7DT

Telephone: 0535 634542 / 636757 Fax: 0535 635582

Registration No. 1291765 in England