

Packaging Trial & Validation

The following should be checked:

- The unit being sent is to the specification as detailed on the official order.
- The disposable tubing INT-6181 is included.
- The instructions for use and service sheet are included.

Packaging:

- Use one cardboard box: 14 x 35 x 23 cm (H x L x W).
- Wrap instrument in bubble-pack.
- Place in box with airbags surrounding – applied with use of the Specific airbag machine – (Cell-O-Ez – Air Cushion Machine).

Postage:

- For single and multiple orders – Despatch is via an applicable courier with adequate insurance.
- Alternative Delivery is via a Hospital visit by a Sales Representative.