


TOM THUMB INFANT RESUSCITATOR

TOM THUMB INFANT RESUSCITATOR DATA SHEET



Please specify MAX Pressure for Safety blow-off valve PT. NO. 9901-17

50 [] 40 [] 30 [] 20 [] Other []

TO ORDER

Choose layout of components and number the circle with the hole number.

Limitation:

Flowmeter PT. NO. 6575 must be vertical in port '6'.

Adjustable Blow-off valve must be vertical in the position shown above.

Safety Blow-off valve PT. NO. 9901-17 must be in a vertical position on the bottom face of the block.

Contact

Order Number:

Title

Invoice Number:

Department

Serial Number:

Hospital

Q/A

Tel. No.

Date:

Additional Information:

FLOW METER
FIXED TO BLOCK
PT NOS 6575

Hose RT. Angle

Hose Assembly
PT NOS 9901-16

Flowmeter to
Block Adapter

RESUSCITATION
BLOCK
9901-01

ADJUSTABLE
BLOW-OFF
VALVE
PTNO. 9901-18

MODURA
CLAW
PT NOS 2002
TMS 2002

BLANKING
PLATE
PT NO 4850-12

O₂ TUBE
OUTLET
PT NO 4701-61

15mm TAPER
MALE

SAFETY
BLOW-OFF
VALVE
PRESIST
20cm H₂O - 50cm H₂O
PT NO 9901-17

PRESSURE
GAUGE
0-100 CM H₂O
50 → 50cm H₂O

PT NO 9901-15




VIAMED




1. All Ports interconnected
2. Inlet & Outlet Ports 1,2,3,4, provide choice of fitting for Tubing connectors
3. Two Blanking plates required
4. Dimensions not critical

(0-15 L/min ST0)





Flow meter
0-10 L/min
Fixed to Block




Adjustable
Blow-off
Valve

Modura
Claw



Pressure Gauge 0-100 cm H2O



Safety
Blow-off
Valve

Newcastle upon Tyne Resuscitation Block

Separate Flow meter & Block for

Emergency ventilation of vent fresh cdy blow off can be set to read vent Pres.

26/04/93

V 1.1

J.S.L.


Viamed Ltd. 15 Station Road, Cross Hills, Keighley, West Yorkshire BD20 7DT

Telephone: 0535 634542 / 636757 Fax: 0535 635582

Registration No. 1291765 in England

TT86

1. All Ports interconnected
2. Inlet & Outlet Ports 1,2,3,4, provide choice of fitting for Tubing connectors
3. Two Blanking plates required
4. Dimensions not critical



Newcastle upon Tyne Resuscitation Block


26/04/93
V 1.1
JJESLL.



VIAMED



1. All Ports interconnected
2. Inlet & Outlet Ports 1,2,3,4, provide choice of fitting for Tubing connectors
3. Two Blanking plates required
4. Dimensions not critical



Pressure Gauge 0-100 cm H₂O

Newcastle upon Tyne Resuscitation Block

26/04/93
V 1.1
J.S.L.

Viamed Ltd. 15 Station Road, Cross Hills, Keighley, West Yorkshire BD20 7DT
Telephone: 0535 634542 / 636757 Fax: 0535 635582
Registration No. 1291765 in England


TT 84



VIAMED



1. All Ports interconnected
2. Inlet & Outlet Ports 1,2,3,4, provide choice of fitting for Tubing connectors
3. Two Blanking plates required
4. Dimensions not critical



Pressure Gauge 0-100 cm H₂O

Newcastle upon Tyne Resuscitation Block

26/04/93

V 1.1

J.S.L.

Viamed Ltd. 15 Station Road, Cross Hills, Keighley, West Yorkshire BD20 7DT

Telephone: 0535 634542 / 636757 Fax: 0535 635582

Registration No. 1291765 in England

TT88