




VIAMED



1. All Ports interconnected
2. Inlet & Outlet Ports 1,2,3,4, provide choice of fitting for Tubing connectors
3. Two Blanking plates required
4. Dimensions not critical



Pressure Gauge 0-100 cm H₂O

Newcastle upon Tyne Resuscitation Block

26/04/93

V 1.1

J.S.L.

Viamed Ltd. 15 Station Road, Cross Hills, Keighley, West Yorkshire BD20 7DT

Telephone: 0535 634542 / 636757 Fax: 0535 635582

Registration No. 1291765 in England

TT 84

FLOW METER
FIXED TO BLOCK
P7 Nos 6575

ADJUSTABLE
BLOW-OFF
VALVE
P7 No. 9901-18

HOSA RT. ANGLE


HOSA ASSEMBLY
P7 Nos 9901-16

Flowmeter to
Block Adapter

MODURA
CLAW
P7 Nos 2002
TMS 2002

RESUSCITATION
BLOCK
9901-01

BLANKING
PLATE
P7 No 4850-12



P7 No 9901-15

SAFETY
BLOW-OFF
VALVE

PRESSURE
20 cm H₂O - 50 cm H₂O
P7 No 9901-17

O₂ TUBE
No. OUTLET
P7 No 4701-61


15mm TAPER
MALE



VIAMED



1. All Ports interconnected
2. Inlet & Outlet Ports 1,2,3,4, provide choice of fitting for Tubing connectors
3. Two Blanking plates required
4. Dimensions not critical



Pressure Gauge 0-100 cm H₂O

Newcastle upon Tyne Resuscitation Block

26/04/93

V 1.1

J.S.L.

Viamed Ltd. 15 Station Road, Cross Hills, Keighley, West Yorkshire BD20 7DT

Telephone: 0535 634542 / 636757 Fax: 0535 635582

Registration No. 1291765 in England

TT88