



Main Account <viamedinbox@gmail.com>

FW: Bodyjet [Scanned by Freecom.net]

1 message

Tom Ladds <tom@humanmeduk.com>
To: Humanmeduk <logistic@humanmeduk.com>
Cc: Daniel Lear <daniel@humanmeduk.com>

14 December 2017 at 11:01

Hi

Please see below order from Stratford Dermotherapy.

The infiltration kits mentioned is one box of 500001-5.

If we no longer have 670200-5 in stock then he should have 671200-5 as an alternative which will be fine.

Thanks

Tom

From: Stratford Dermatherapy Clinic [mailto:info@skincareclinic.co.uk]
Sent: 14 December 2017 10:31
To: Tom Ladds <tom@humanmeduk.com>
Subject: RE: Bodyjet [Scanned by Freecom.net]

Tom

Hope you are well and all is good with fat sucking...In the New Year can you send me 1 x box of 5 Bodyjet infiltration kits and also a box of 670200 O-rings for the lipo-collector 3 kit please..

Have a great Christmas matey

Hugo

On 04/05/2016 14:29, Tom Ladds wrote:

Yes, it can be difficult to manage. I have seen some inventive use of towel clips to keep it out of the way and the weight at your hand down!

From: Clinic [<mailto:info@skincareclinic.co.uk>]
Sent: 04 May 2016 11:02
To: Tom Ladds <tom@humanmeduk.com>
Subject: Re: Bodyjet [Scanned by Freecom.net]

Good

It's just that the infiltration kit introducer is so long, if you're doing the face?

H

Sent from my iPhone

On 4 May 2016, at 10:26, Tom Ladds <tom@humanmeduk.com> wrote:

Yes. Or you can infiltrate with the Bodyjet in the normal way with the pump pressure set quite low.

Sent from my Samsung device

----- Original message -----

From: Clinic <info@skincareclinic.co.uk>
Date: 04/05/2016 10:19 (GMT+00:00)
To: Tom Ladds <tom@humanmeduk.com>
Subject: Re: Bodyjet [Scanned by Freecom.net]

So

Can I infiltrate manually and suck with the cannula?

H

Sent from my iPhone

On 4 May 2016, at 09:43, Tom Ladds <tom@humanmeduk.com> wrote:

It's a yes to the cannula but a no to the application set. The latter is not compatible.

T

Sent from my Samsung device

----- Original message -----

From: Clinic <info@skincareclinic.co.uk>
Date: 04/05/2016 09:31 (GMT+00:00)
To: Tom Ladds <tom@humanmeduk.com>
Subject: Re: Bodyjet [Scanned by Freecom.net]

So is that a yes, we can adapt?

Sent from my iPhone

> On 4 May 2016, at 09:13, Tom Ladds <tom@humanmeduk.com> wrote:

>

> Good morning Hugo

>

> I'm afraid the standard Bodyjet is not compatible with the cheaper 'Biofill'

> WAL applicator set. Either the Evo or the desk-top Eco can use this.

>

> We have a full range of 15cm cannula that fit, via a spacer, to the end of

> the standard applicator. The smallest of these are a 2.5mm infiltration

> cannula and a 3.5mm for aspiration; most of the guys who do face work use

> these. Happy for me or one of the team to come in if needed.

>

> Kind regards

>

> Tom

>

>

>

>

> -----Original Message-----

> From: Clinic [<mailto:info@skincareclinic.co.uk>]

> Sent: 03 May 2016 22:03

> To: Tom Ladds <tom@humanmeduk.com>

> Subject: Re: Bodyjet [Scanned by Freecom.net]

>

> Tom

>

> Can you tell me what small infusion kit / hand pieces are available to do

> Bodyjet Lipo on a neck using my existing machine? I do not have the Evo and

> have no need for it?

>

> Hugo

>

> Sent from my iPhone

>

>> On 4 Mar 2016, at 18:00, Tom Ladds <tom@humanmeduk.com> wrote:

>>

>> Hi

>>

>> That's no problem. We have a good experienced engineer on Bodyjet so

>> feel free to get in touch if your chap needs any advice.

>>

>> Have a good weekend.

>>

>> Tom

>>

>> -----Original Message-----

>> From: Clinic [<mailto:info@skincareclinic.co.uk>]

>> Sent: 03 March 2016 18:52

>> To: Tom Ladds <tom@humanmeduk.com>

>> Subject: Re: Bodyjet [Scanned by Freecom.net]

>>

>> Hi Tom

>>

>> I'm organising an engineer to service all my machines under one

>> umbrella and he will want to contact you guys if necessary ref parts

>> or consumables as and when. Nice to have a point of contact.

>>

>> Hope you're well
>>
>> Hugo
>>
>> Sent from my iPhone
>>
>>> On 3 Mar 2016, at 18:14, Tom Ladds <tom@humanmeduk.com> wrote:
>>>
>>> Good evening Hugo
>>>
>>> We are indeed. Your local sales consultant is Daniel Lear, e-mail
>>> above or on 07496 479083.
>>>
>>> Colin Pyne has now retired from the business so I am now covering his
>>> old role. Would you like me to ask Daniel to get in touch?
>>>
>>> Kind regards
>>>
>>> Tom
>>>
>>> -----Original Message-----
>>> From: Clinic [<mailto:info@skincareclinic.co.uk>]
>>> Sent: 03 March 2016 13:47
>>> To: tom.ladds@cnpmmedical.com
>>> Subject: Bodyjet
>>>
>>>
>>> Tom
>>>
>>> Are you guys ie Humanmed UK, still responsible for Bodyjet stuff and
>>> spares/disposables ?
>>>
>>> Hugo
>>>
>>>
>>> Sent from my iPhone=
>
>