


versastream luer lock male selection guide



Is the patient intubated?

Yes

Long-Term
≤ 72 hours

Airway Adapter

Adult/Paediatric
4420873

Infant
4420874



Airway adapter
connector sizes:

Adult/Paediatric 15 mm O.D.
to 15 mm I.D. / 22 mm O.D.

Infant 15 mm O.D.
to 4 mm I.D. / 15 mm I.D.



Luer Lock to Luer Lock
Sampling Line
for use with 4420531
4420876



Airway Adapter with Luer
Lock Female Connection
e.g. Adult/Paediatric
4420531

No

Long-Term
≤ 72 hours

Does the patient require oxygen?

Yes

Nasal

Adult
4420860

Paediatric
4420861

Infant
4420862



Nasal/Oral

Adult
4420869

Paediatric
4420870



For use with patients
that exceed 10 kg



For use with patients
that are less than 10 kg

No

Nasal

Adult
4420857

Paediatric
4420858

Infant
4420859



Nasal/Oral

Adult
4420867

Paediatric
4420868



Part Number: 4490117
Date: 10/18