

Product Labelling

Product: TT-480
 Part Number: 0310030
 Date: 6th December 2018

Device Labels

Tom Thumb Infant Resuscitator

CAUTION
 Must be used with a flow
 meter. **DO NOT** exceed
 15 litres/minute.
 Read instructions
 prior to use.

Model: Tom Thumb
 TT480
 PN: 0310030
 SN: 0401757



Viamed Ltd.
 15 Station Road
 Cross Hills
 Keighley
 West Yorkshire BD20 7DT
 United Kingdom
 Tel: +44 (0)1535 634542
 Fax: +44 (0)1535 635582
 e-mail: info@viamed.co.uk
 Web site: www.viamed.co.uk

Box Label



P/N 0310030

LOT PR02680

S/N 0401513

2018-11



Viamed Limited, 15 Station Road, Cross Hills, Keighley
 West Yorkshire BD20 7DT United kingdom



Ryan Swaine - 181206



ISO 13485:2003
 ISO 13485:2003 CMDCAS Recognised
 Company registered in England, No. 01291765