


Viamed R-22Vi Oxygen Sensor Packaging & Sensor Labelling

0110020 R-22Vi

Medical Oxygen Sensor R- 22Vi


REF 0110020 **SN** V123456 2019-03-01



(01)04036616000352(11)190301(21)V123456


VIAMED

 0123

Tel: +44 (0)1535 634542
 Email: info@viamed.co.uk Website: www.viamed.co.uk


 EnviteC-Wismar GmbH
 Alter Holzhafen 18
 D-23966 Wismar
 www.envitec.com

Depending on the ambient conditions the sensor may require up to 30 minutes to reach signal stability after installation. The sensor contains caustic liquid. In case of leakage avoid contact with eyes and skin! Recommended storage temperature: 5°C - 15°C!

Packaging label: 63.5 x 46.6 mm


VIAMED

Medical Oxygen Sensor R- 22Vi

REF 0110020 **SN:** V123456




 EnviteC-Wismar GmbH
CE 0123



 2019-03

Sensor label: 90 x 15 mm