

Viamed R-22MEDV Oxygen Sensor Packaging & Sensor Labelling

0110021 R-22MEDV

Medical Oxygen Sensor R- 22MEDV








REF 0110021

SN V123456

2019-03-01



(01)04036616000376(11)190301(21)V123456



CE 0123

EnviteC-Wismar GmbH
 Alter Holzhafen 18
 D-23966 Wismar
 www.envitec.com

Tel: +44 (0)1535 634542
 Email: info@viamed.co.uk Website: www.viamed.co.uk

Depending on the ambient conditions the sensor may require up to 30 minutes to reach signal stability after installation. The sensor contains caustic liquid. In case of leakage avoid contact with eyes and skin! Recommended storage temperature: 5°C - 15°C!

Packaging label: 63.5 x 46.6 mm



Medical Oxygen Sensor R- 22MEDV

REF 0110021

SN: V123456






2019-03

EnviteC-Wismar GmbH
CE 0123

Tel: +44 (0)1535 634542
 Email: info@viamed.co.uk
 Website: www.viamed.co.uk

Sensor label: 90 x 15 mm