


Viamed R-22VA Oxygen Sensor Packaging & Sensor Labelling

0110039 R-22VA



Medical Oxygen Sensor R- 22VA








REF 0110039 SN V123456 2019-03-01

(01)04036616010740(11)190301(21)V123456



VIAMED

Tel: +44 (0)1535 634542
 Email: info@viamed.co.uk Website: www.viamed.co.uk

CE 0123
 EnviteC-Wismar GmbH
 Alter Holzhafen 18
 D-23966 Wismar
 www.envitec.com

Depending on the ambient conditions the sensor may require up to 30 minutes to reach signal stability after installation. The sensor contains caustic liquid. In case of leakage avoid contact with eyes and skin! Recommended storage temperature: 5°C - 15°C!

Packaging label: 63.5 x 46.6 mm




VIAMED

Tel: +44 (0)1535 634542
 Email: info@viamed.co.uk
 Website: www.viamed.co.uk

Medical Oxygen Sensor R- 22VA
 REF 0110039 SN: V123456





 EnviteC-Wismar GmbH
 CE 0123




 2019-03

Sensor label: 90 x 15 mm