

## Viamed R-41V Oxygen Sensor Packaging & Sensor Labelling

**0110041 R-41V**

**Medical Oxygen Sensor R- 41V**


**REF 0110041**
**SN V123456**
2019-03-01



(01)04036616000512(11)190301(21)V123456

---


**VIAMED**

**CE 0123**  
 EnviteC-Wismar GmbH  
 Alter Holzhafen 18  
 D-23966 Wismar  
 www.envitec.com

Tel: +44 (0)1535 634542  
 Email: info@viamed.co.uk Website: www.viamed.co.uk

Depending on the ambient conditions the sensor may require up to 30 minutes to reach signal stability after installation. The sensor contains caustic liquid. In case of leakage avoid contact with eyes and skin! Recommended storage temperature: 5°C - 15°C!

Packaging label: 63.5 x 46.6 mm


**VIAMED**

**Medical Oxygen Sensor R- 41V**  
**REF 0110041**
**SN: V123456**



2019-03

---


 EnviteC-Wismar GmbH  
**CE 0123**

Tel: +44 (0)1535 634542  
 Email: info@viamed.co.uk  
 Website: www.viamed.co.uk

Sensor label: 90 x 15 mm