







## Viamed R-43V Oxygen Sensor Packaging & Sensor Labelling


**0110043 R-43V (Pack of 2 sensors)**

**Medical Oxygen Sensor R- 43V**


  
  


REF 0110043 SN V123456 2019-03-01




(01)04036616004596(11)190301(21)V123456

---



Tel: +44 (0)1535 634542  
 Email: info@viamed.co.uk Website: www.viamed.co.uk


CE 0123



EnviteC-Wismar GmbH  
 Alter Holzhafen 18  
 D-23966 Wismar  
 www.envitec.com

Depending on the ambient conditions the sensor may require up to 30 minutes to reach signal stability after installation. The sensor contains caustic liquid. In case of leakage avoid contact with eyes and skin! Recommended storage temperature: 5°C - 15°C!


Packaging label: 63.5 x 46.6 mm






Tel: +44 (0)1535 634542  
 Email: info@viamed.co.uk  
 Website: www.viamed.co.uk

**Medical Oxygen Sensor R- 43V**

REF 0110043 SN: V123456



EnviteC-Wismar GmbH  
 CE 0123

2019-03

Sensor label: 90 x 15 mm