pulse oximetry



MD300-C63

Finger Pulse Oximeter

Ultra-compact, lightweight pulse oximeter with colour OLED display



- Adjustable brightness: 10 levels
- Pleth waveform, displayed in 2 formats
- Signal strength indicator
- Easy to use; 1-button operation
- Automatic power-off when finger removed
- Special battery cover for more durability
- Suitable for paediatric and adult patients

The MD300-C63 can be used in a variety of applications, for use in GP surgeries, clinics and dental surgeries, through to spot-checking before, during and after exercise.

The bright organic LED display offers superb visibility in all lighting conditions, including direct sunlight.



Actual Size



MD300-C63 Finger Pulse Oximeter

Measuring Range

SpO₂ 70 - 100% Pulse Rate 30 - 235 bpm

Accuracy

Display

Features Colour OLED, adjustable brightness (10 levels)

Data displayed Oxygen saturation, pulse rate, plethysmogram, pulse indicator bar graph

Indicators

Finger out, signal strength, battery status

Environmental Conditions

Operating temperature 5°C to 40°C Storage temperature -20°C to 55°C

Regulatory Approvals

Standards CE Class Ila

Assembly

Dimensions 58 mm (L) x 34 mm (W) x 36 mm (H)
Weight Approximately 55 g (including batteries)
Power Supply 2 x AAA/LR03/MN2400 1.5V alkaline batteries
Battery Life Approximately 30 hours of continuous operation

Automatically powers off when finger is removed for longer than 8 seconds

Part Number

MD300-C63 2810038

Standard accessories included:

Lanyard

2 x AAA Batteries Operator Manual

Optional Accessories

Soft carry case, with belt loop, colour - black
Soft carry case, with belt loop, colour - orange
Soft carry case, with belt loop, colour - orange
Soft carry case, with belt loop, colour - green
Soft carry case, with belt loop, colour - blue
0022192

Clinical trials carried out in accordance with BS EN ISO 9919:2005 standards. Materials tested in accordance with ISO 10993-, -5, -10.



2nd Rotation



3rd Rotation



4th Rotation



Pleth Waveform 1



Pleth Waveform 2

C € 0123
Specifications Subject to Change

Distributed by:



NIAMED

Viamed Limited · 15 Station Road · Cross Hills Keighley · West Yorkshire · BD20 7DT · United Kingdom Tel: +44 (0)1535 634 542 Fax: +44 (0)1535 635 582 Email: info@viamed.co.uk Website: www.viamed.co.uk

