pulse oximetry



MD300-C15D

Finger Pulse Oximeter

Ultra-compact, lightweight pulse oximeter with LCD display

- Low cost
- Easy to use; 1-button operation
- Automatic power-off when finger is removed

Inexpensive pulse oximeter for general applications



Actual Size

The MD300-C15D can be used in a variety of applications, from use in hospitals, paramedic ambulances, clinics and dental surgeries, through to spot-checking before, during and after exercise.



MD300-C15D Finger Pulse Oximeter

Measuring Range	
SpO ₂	70 - 100%
Pulse Rate	30 - 250 bpm
Accuracy	
SpO ₂ (70 - 100%)	± 2 digits
SpO ₂ (0 - 69%)	No definition
SpO ₂ Resolution	± 1%
Pulse Rate; beats per minute (30 - 99 bpm)	± 2 bpm
Pulse Rate; beats per minute (100 - 250 bpm)	± 2%
Pulse Rate Resolution	1 bpm
Display	

Features LCD display

Data displayed Oxygen saturation, pulse rate, pulse indicator, bar graph

Indicators Signal strength, battery status

Environmental Conditions

5°C to 40°C Operating temperature Storage temperature -20°C to 55°C

Regulatory Approvals

Standards CE Class IIa

Assembly

58.3 mm (L) x 33.2 mm (W) x 37.0 mm (H) **Dimensions** Weight Approx. 50 g (including batteries) Power Supply 2 x AAA/LR03/MN2400 1.5V alkaline batteries

> Automatically powers off when finger is removed for longer than 8 seconds

Part Numbers

MD300-C15D 2810046

Standard accessories included:

Lanyard

2 x AAA Batteries

Operator Maridar	
Optional Accessories	
Soft Carry Case, with belt loop, colour - black	0022190
Soft Carry Case, with belt loop, colour - orange	0022191
Soft Carry Case, with belt loop, colour - green	0022192
Soft Carry Case, with belt loop, colour - blue	0022193



Clinical trials carried out in accordance with BS EN ISO 9919:2005 standards. Materials tested in accordance with ISO 10993-, -5, -10.



Optional carrying cases with belt loop and D-ring

C€0123

Specifications and Case Colouring Subject to Change

Distributed by:



W VIAMED

Viamed Limited · 15 Station Road · Cross Hills Keighley · West Yorkshire · BD20 7DT · United Kingdom Tel: +44 (0)1535 634 542 Fax: +44 (0)1535 635 582 Email: info@viamed.co.uk Website: www.viamed.co.uk

