



Bottom Layer

Brief Summary for the Changes

- 1, Pads coloured in purple are enlarged to 1.75mm or above. (Diameter of circles, Major Diameter of ellipses, and longer side of rectangles)
- 2, Ground plate coloured in green is deleted.
- 3, Trace coloured in orange is added.

			COMPANY:		HEALTHCARE TECHNOLOGY INTERNATIONAL LTD.		
			DRAWN:		TITLE:		
			DATE:		MICROSTIM DB		
			CHECKED:		PAGES:	DOC .NO:	REV: 0
			DATE:		8	ED-68801	
REV	DESCRIPTION	DATE	APPROVED:		FILE:	PCB P/N:	
REVISIONS			DATE:		688 MD 000.PCB	688-66748	

FORM-0055A