

Packaging:

'Pack'

The **NeoPEEPTM** Single Use Resuscitation Circuit will be individually wrapped in clear plastic heat sealed bags. They do not need to be sterile, but must be 'clinically clean'. A Product label is to be stuck on to the outside the plastic bag with the 'pack' section ticked (see figure 1 (a)).

The dimensions of the **NeoPEEPTM** Single Use Resuscitation Circuit in the sealed plastic bags will be:

30cms X 20cms




Figure 1 Product Label



'Inner Box'

The single pack of the **NeoPEEPTM** Single Use Resuscitation Circuit then has to be placed into a box in 2 rows of 10 side by side. The approximate size of the box should be:

36cms X 27cms X 15cms

Inside of this box there should be 1 X Instructions for Use. (Figure 2)




Figure 2
Instructions for Use

On the outside of the box there needs to be a product label stuck to the longer side of the box with the 'inner box' section ticked. (Figure 1(b))




'Outer Case'

The 'Outer Case' should contain 20 boxes of 'inner boxes' totalling 400 NeoPEEPTM Single Use Resuscitation Circuits per case.

A Product Label is to be stuck on the outside of this box with the 'Outer Case' box

ticked. (Figure 1 (c))

There should be the standard symbols displayed on the outer case as in all Marshall Products cases have at present (Figure 3)



N.W

XX Kgs

G.W

XX Kgs

PRODUCT CODE: 8001XX/XXXX

C/No:

XX

(Including Product Label)

Figure 3 **Outer Case Symbols**