



**Revised design of Microstim rear panel label, DB35 ; 19<sup>th</sup> May 03.**

**Microstim DB.**  
**Summary instructions for use.**

Press and hold switch. Release to stop.

**T.O.F :** 4 pulses at 2Hz.

**D.B.S :** 3 pulses at 50Hz, followed by a 0.75 second pause, then 2 pulses at 50Hz.

**P.T.C :** 50Hz for 5 seconds, followed by a 3 second pause, then continuous 1Hz.

**1 Hz :** Continuous 1 Hz.

Set output using intensity control. Audio indicates patient current flowing; 100mA maximum, 200µs pulses, constant current.

Replace battery when flash is deep red.  
Uses 1 x MN1604, 9V alkaline battery.



**Refer to User Manual.**  
**Use with care with**  
**cardiac pacemakers.**



Viamed Ltd, 15 Station Road, Cross Hills,  
West Yorkshire, BD20 7DT, UK.

Tel : +44 (0)1535 634542.

Fax : +44 (0)1535 635582.

Email : info@viamed.co.uk

Overall dimensions 64mm x 45mm (h x w).

*Quote reference:-  
PJS 1828/0503*

*Total costs*

*200 OFF  
@ 66p EACH  
+ VAT.*

*Thanks  
Peter*