

VITAL STATISTICS OF CURRENT EQUIPMENT

VICKERS MOBILE RESUSCITAIRE

Height - 1040 mm
Width of baby tray - 370 mm

FIXED RESUSCITATION TABLES AT NGH

Height - 1100 mm
Width of baby tray - 260 mm
Shelf width - 750 mm
Shelf depth - 530 mm

NEW RESUSCITATION CABINETS IN RVI

We have agreed the following essential modifications need to be made before this equipment becomes safe to use:

1. The shelf supports need to be adjusted so that the shelf cannot come down beyond 90° (85° is better than 95°)
2. The existing baby tray should be scrapped and redesigned. It could take the form of a 1 piece, flat but concave cylinder section made of moulded plastic. Consideration needs to be given to the dimensions and to the positioning on the shelf.
3. The heater should be positioned directly over the baby tray in such a way that it can be swung flat against the wall and locked in that position.
- ④ The entire cabinet and heater should be moved up the wall so that the height of the shelf when horizontal from the floor is 1000 mm.

Bygge Cateh

Heater Bumber