

Service Repair Sheet SRS67570

Date Promised to Return: 01 Jan 1970 -3600

Contact Name

Joseph Bracey

Company/ Hospital Name

Medloc Trade Receipts (MOD)

Department

Position

Direct Phone

0117 332 8732

General Phone

01225 884159

Opera Account

00012481

Email

LEIDOS-TVS-Med-RepairSection@mod.gov.uk

Order Number

Date Received

09/Sep/2020

Booked in By

Robert Connor

Main Company

Viamed

Type Return

Quote

Notes 25/Aug/2020 Zoey Teal

25/Aug/2020 Zoey Teal

Sending in a AX300. Serial number 280730. 3.5mm Jack Connection On The Side Of The Equipment Has Come Away From The PCB Inside Due To Locking Nut Being Missing.

09/Sep/2020 Robert Connor

Received 1 x AX300-i s/n 280730, 0131202 s/n 0010505, R-17MED s/n 158142, with flow diverter, 3 x AA battery and base stand.

Ready For quote

CGreen 22-9-2020

Repair Complete Signed

SRN	Equipment	Stock Ref	Serial Number	Warranty
SRN33103	Oxygen Monitor	0111230	280730	2
SRN33104	Oxygen Monitor Cable	0131202	0010505	1
SRN33105	Oxygen Sensor	0110017	158142	1

0140004 x 1 @ £0.00

SN, SRS, SRN

NO fault found.

0040001 x 1 @ £0.00

SN, SRS, SRN

0140012 x 1 @ £0.00

Old SN, SRS, SRN.

0110017 x 1 @ £42.00

New SN, SRS

UPS x 1 @ £10.00

101 Instruction - Supplier Site

Your demand has been received and assigned the job number below. Please check the details below and inform us of any errors:	
Job Number:	2273055
Purchase Order Number:	SR567570
Instruction Sent:	25/08/2020
Demand Reference:	0034D2D6
Equipment NSN:	6515997437258
Equipment:	Oxygen Analyser AX300
Serial Number:	280730
Asset Tag (if known):	TMW0442
Work Requested by:	Cpl Welsh 01993 842551 x 4148
Dispatch UIN:	F1510A
Dispatch Location:	MDSS TMW, RAF Brize Norton, Carterton, OX18 3LX
Please send this equipment to the following supplier:	
Supplier:	Viamed Ltd
Supplier address:	15 Station Road, Cross Hills, Keighley, West Yorkshire, BD20 7DT
Include a copy of this instruction with this equipment. Keep a copy for your records and notify LEIDOS-TVMS-MED-REPAIRSECTION@MOD.UK of the despatch date:	
Despatch Date:	01/09/2020
On return of the equipment complete the below and return this to LEIDOS-TVMS-MED-REPAIRSECTION@MOD.UK:	
Received date:	
By:	
Condition:	Satisfactory/Unsatisfactory/BER
Comment:	