

Service Repair Sheet SRS67654

Contact Name

Company/ Hospital Name

Department

Position

Direct Phone

General Phone

Opera Account

Email

Order Number

Date Received

Booked in By

Main Company

Type Return

Date Promised to Return: 01 Jan 1970 -3600

Andrew Greaves

Glenfield General Hospital

Medical Physics

Senior Maintenance Assistant

0116 258 3511

00001184

andrew.greaves@uhl-tr.nhs.uk

07/Dec/2020

Catherine Spence

Viamed

Quote

O/N: G4596516

Notes Sending in a DB3. The plastic control knob has broken on the inside and now just spins on top of the device. Also an issue with one of the switches. When you set it to D.B.S. the switch sticks and doesn't unlatch. Serial number M0004823.

26/Nov/2020 Zoey Teal

Sending in a DB3. The plastic control knob has broken on the inside and now just spins on top of the device. Also an issue with one of the switches. When you set it to D.B.S. the switch sticks and does not unlatch. Serial number M0004823.

07/Dec/2020 Catherine Spence

received unit in, no leads or battiers

Ready For quote

CGNEEN 8-12-2020

Repair Complete Signed

CGNEEN 4-1-2021

SRN	Equipment	Stock Ref	Serial Number	Warranty
SRN33266	Microstim DB3	2510000	M0004823	N

2540010 x1 @ E45
S/N, SRS, SRN

UPS x 1 E10