

# **Neuromuscular Monitoring in Clinical Practice**

Possible Solutions and practical realization



## **Neuromuscular Blocking Agents (NMBA)**

- NMBA belong to the most frequently used agents in general anesthesia
- Support Intubation
- Maintain surgical conditions



## Neuromuscular Blocking Agents (NMBA)

- Onset as well as duration depend on the condition of the respective patient.
- Duration of relaxation as well as muscular recovery depend on several external factors
  - Accompanying disease, Interaction with other drugs, pharmacogenetic factors\*



## Neuromuscular Blocking Agents (NMBA)

- Residual paralysis
  - Risk for pulmonary events
  - Impaired ability to swallow
  - Post operative malaise of patients



- "No Monitoring" (Duration of effect)
  - "Clinical Signs"
- Subjective Neuromuscular Monitoring (tactile / visual evaluation)
- Objective (quantitative) Neuromuscular Monitoring



## "No Monitoring"

In many cases NMBA are administered without proper monitoring.

Survey (2018) Use of NMT monitoring

Europe: 20%

**US: 10%** 

Australia: 10%

**Never use NMT Monitoring!** 

Subjective monitoring in < 40% of all patients
Objective monitoring in < 20% of all patients



## "No Monitoring"

- In 12% of hospitals and 67% of private centers NMT monitoring is not available.
- In 44% of hospitals and 63% of private centers decisions about additional relaxation are based on clinical signs! (Coughing, Bear Down) or "fixed time intervals".
- Even if reversal agents are administered the incidence of residual blockade is between 20-40%



## "No Monitoring"

- In 54% of hospitals and 70% of private centers clinical signs are the most frequently used criteria for the decision of neuromuscular reversal
- Evaluation of residual paralysis based on clinical tests
  - Tongue spatula test
  - Head Lifting >5 Sec. (possible at TOF = 0,5!)

The tests can only be applied in the very late phase of anesthesia. Only AFTER extubation!



#### Mechanomyography - MMG

- Measures the force of the muscle (Pressure sensor)
- Preload of muscle is necessary
- Sensitive technical setup
  - Small changes in the positioning of the patient may influence the measurement result
- Mostly used for research and studies
- Reference method for NMT Monitoring



#### **Electromyography - EMG**

- Measures the electric action potentials of the muscle
- Digital signal processing
- Provided a low noise level, the signal quality and measurement results are close to MMG.
- Susceptible to electric interferences, temperature and movement



#### **Acceleromyography - AMG**

- Measures acceleration of the muscle
- 2. Newtons Law: F = m x a (force = mass times acceleration)
- Mass constant → Force proportional to acceleration
- Piezo crystal changes potential during acceleration. Changes are measured and amplified.
- Free movement of the respective muscle is necessary!



- Various peripheral Nerve stimulators
  - Evaluation is done visually / tactile
- Stimulation modes
  - Single Twitch (Single Stimulation)
  - Train of Four (TOF)
  - Double Burst (DBS)
  - Tetanic Stimulation
  - Post Tetanic Stimulation (PTC)



- Studies have shown that when using TOF stimulation even experienced physicians cannot tell the difference between muscular response if TOF is > 0,5.
  - Unexperienced users if TOF = 0,4
- Using DBS the result is "improved" to TOF = 0,6.
- DBS is better than TOF if subjective monitoring is used.




A study recently conducted in the USA showed that the assessment of TOF = 0,9 has a 47% hit rate.



= 50:50



#### Residual paralysis in the PACU



Residual Paralysis in the PACU after a Single Intubating Dose of Non-depolarizing Muscle Relaxant with an Intermediate Duration of Action (Debaene et al.)



# Vital functions during recovery

	TOF 0,5	TOF 0,8	TOF 1
Tidal volume	Normal	Normal	Normal
Forced Vital capacity	Often compromised	Often normal	Normal
Pharyngeal Function (Swallowing)	Certainly compromised	Often compromised	In general normal
Integrity of upper Airway	Certainly compromised	Often compromised	In general normal
Hypoxic breathing response	Often compromised	Often normal	Normal



AMG has shown to be a practical method for NMT monitoring in clinical practice.

Measurement of muscle acceleration using an acceleration transducer (Piezo Crystal) – Quantifiable! Strength of the muscle response can be measured by transducer.


Can be used on different muscles – Adductor pollicis is most suitable.



Calculation of TOF ratio: 4 Stimulation impulses
 0-4 Responses.

**TOF Ratio = T4 / T1** 

- No block
- Partly NM block
- Strong NM block TOF Count
- Complete NM block






- Different NMT modules on the market
  - (GE, Philips)
- Handheld Devices
  - TOF 3D
  - TOFScan
  - StimPod
  - TetraGraph
  - TwitchView



## Alternative electrode placement

Facial nerve: place negative electrode (black) by ear lobe and the positive (red) 2cms from the eyebrow (the along facial nerve inferior and lateral to eye)

Response: Orbicularis occuli muscle – eyelid twitching





## Alternative electrode placement

Sural (posterior tibial) nerve: place the negative electrode (black) over inferolateral aspect of medial malleolus (palpate posterior tibial pulse and place electrode there) and positive electrode (red) 2-3cm proximal to the negative electrode

Response: Fexor hallucis brevis muscle – planter flexion (curl) of big toe





#### References

- Neuromuskuläres Monitoring in Klinik und Forschung (Fuchs-Buder)
- Einsatz des Neuromuskulären Monitorings in Deutschland (Der Anästhesist – 09/2008 pp. 908-914)
- Residual Paralysis in the PACU after a Single Intubating Dose of Nondepolarizing Muscle Relaxant with an Intermediate Duration of Action (Debaene et al. Anesthesiology 2003; 98:1042–8)
- A Survey of Current Management of Neuromuscular Block in the United States and Europe (Naguib et al. International Anesthesia Research Society, 2009)