



Zoey Teal &lt;viamed.zoey.teal@gmail.com&gt;

---

**Fwd: Alternative oxygen sensor**

1 message

---

**Main Account** <office@viamed.co.uk>  
To: Zoey Teal <zoey.teal@viamed.co.uk>

Wed, Apr 7, 2021 at 10:11 AM

----- Forwarded message -----

From: **Steve Nixon** <office@viamed.co.uk>  
Date: Wed, 7 Apr 2021 at 10:09  
Subject: Fwd: Alternative oxygen sensor  
To: Main Account <office@viamed.co.uk>

Please send 1 x 0110435 today and send by UPS Express Saver to Denzil Searle at Inspiration Healthcare.

Supply FOC for engineering evaluation purposes.

Steve

----- Forwarded message -----

From: **Denzil Searle** <Denzil.Searle@coopersurgical.com>  
Date: Wed, 7 Apr 2021 at 09:10  
Subject: RE: Alternative oxygen sensor  
To: [steve.nixon@viamed.co.uk](mailto:steve.nixon@viamed.co.uk) <[steve.nixon@viamed.co.uk](mailto:steve.nixon@viamed.co.uk)>

Hi, Steve,

Ah, I see. OK, thanks for clarifying that.

Denzil

**Denzil Searle**, *Product Development Manager*

Bickland Industrial Park, Falmouth, Cornwall UK TR11 4TA

CooperSurgical

w: +44 1326 370615 | m: +44 7583 099929

[denzil.searle@coopersurgical.com](mailto:denzil.searle@coopersurgical.com)

---

**From:** Main Account <[viamedinbox@gmail.com](mailto:viamedinbox@gmail.com)> **On Behalf Of** Steve Nixon  
**Sent:** 06 April 2021 18:48  
**To:** Denzil Searle <[Denzil.Searle@coopersurgical.com](mailto:Denzil.Searle@coopersurgical.com)>  
**Subject:** Re: Alternative oxygen sensor

**CAUTION:** This email was sent from an EXTERNAL source. Use caution when clicking links or opening attachments.

Hi Denzil

That is the serial number which we can also use to trace details, to the right of the serial number there should be the two digit date code.

Sold to you on 10/11/2020 Invoice RVM127021-1

Sensor date code: I0 = August 2020



On Tue, 6 Apr 2021 at 15:59, Denzil Searle <[Denzil.Searle@coopersurgical.com](mailto:Denzil.Searle@coopersurgical.com)> wrote:

Hi, Steve,

I have been given some 6-digit numbers, though. For example, 531696.

Denzil



**Denzil Searle**, *Product Development Manager*

Bickland Industrial Park, Falmouth, Cornwall UK TR11 4TA

CooperSurgical w: +44 1326 370615 | m: +44 7583 099929

[denzil.searle@coopersurgical.com](mailto:denzil.searle@coopersurgical.com)

---

**From:** Main Account <[viamedinbox@gmail.com](mailto:viamedinbox@gmail.com)> **On Behalf Of** Steve Nixon

**Sent:** 06 April 2021 15:44

**To:** Denzil Searle <[Denzil.Searle@coopersurgical.com](mailto:Denzil.Searle@coopersurgical.com)>

**Subject:** Re: Alternative oxygen sensor

**CAUTION: This email was sent from an EXTERNAL source. Use caution when clicking links or opening attachments.**

Hi Denzil

There is a two digit alphanumeric date code consisting of a letter and a number:

The first is a letter, A = January, B = February, C = March....

The second is a number: 9 = 2019, 0 = 2020, 1 = 2021

So the code for April 2021 will be D1

Steve

On Tue, 6 Apr 2021 at 15:21, Denzil Searle <[Denzil.Searle@coopersurgical.com](mailto:Denzil.Searle@coopersurgical.com)> wrote:

I have been asked if we can tell from the serial number of the Teledyne sensors when they were made. The AI ones have the month coded into to the start of the serial number. Is there anything similar in these?

Denzil



**Denzil Searle**, *Product Development Manager*

Bickland Industrial Park, Falmouth, Cornwall UK TR11 4TA

CooperSurgical w: +44 1326 370615 | m: +44 7583 099929

[denzil.searle@coopersurgical.com](mailto:denzil.searle@coopersurgical.com)

---

**From:** Main Account <[viamedinbox@gmail.com](mailto:viamedinbox@gmail.com)> **On Behalf Of** Steve Nixon

**Sent:** 06 April 2021 12:11

**To:** Denzil Searle <[Denzil.Searle@coopersurgical.com](mailto:Denzil.Searle@coopersurgical.com)>

**Subject:** Re: Alternative oxygen sensor

**CAUTION: This email was sent from an EXTERNAL source. Use caution when clicking links or opening attachments.**

Hi Denzil

Can you please let me know the pin configuration on your equipment, the connector should be a three pin Molex with either 2 or 3 pins loaded. Can you please confirm which pins are loaded and the assignment of function: -ve, +ve.

Steve

On Tue, 6 Apr 2021 at 09:53, Steve Nixon <[steve.nixon@viamed.co.uk](mailto:steve.nixon@viamed.co.uk)> wrote:

Hi Denzil

We may have a suitable oxygen sensor for you, can you please check to see if there is sufficient room within your equipment? The diameter will be OK, but the length may be an issue

Diameter: 0.90" = 22.86mm

Total length including connection thread: 1.85" = 46.99 mm

Regards

--

Steve

Steve Nixon

Director - Viamed Ltd.

Office: +44 (0)1535 634542

Mobile: +44 (0)7850 252267

GDPR / Privacy Policy available at [www.viamed.co.uk/policy](http://www.viamed.co.uk/policy)

.

--

Steve

Steve Nixon

Director - Viamed Ltd.

Office: +44 (0)1535 634542

Mobile: +44 (0)7850 252267

GDPR / Privacy Policy available at [www.viamed.co.uk/policy](http://www.viamed.co.uk/policy)

.

--

Steve

Steve Nixon

Director - Viamed Ltd.

Office: +44 (0)1535 634542

Mobile: +44 (0)7850 252267

GDPR / Privacy Policy available at [www.viamed.co.uk/policy](http://www.viamed.co.uk/policy)

.

--

Steve

Steve Nixon

Director - Viamed Ltd.

Office: +44 (0)1535 634542

Mobile: +44 (0)7850 252267

GDPR / Privacy Policy available at [www.viamed.co.uk/policy](http://www.viamed.co.uk/policy)

.

--

Steve

Steve Nixon  
Director - Viamed Ltd.

Office: +44 (0)1535 634542

Mobile: +44 (0)7850 252267

GDPR / Privacy Policy available at [www.viamed.co.uk/policy](http://www.viamed.co.uk/policy)

.

---

**3 attachments**



**image001.jpg**  
10K

CooperSurgical



**image004.jpg**  
4K

CooperSurgical



**image006.jpg**  
4K

CooperSurgical