

COMPANY OPERATING PROCEDURES			
Microstim Testing and QA			
Created:	22/01/03	VM3/COP/99.00	3
Revised:		Last printed 24 June 2009	Page 1 of 1

Microstim Testing and QA

Microstim DB's are supplied to Viamed fully assembled and Tested,

- Ensure there is no physical damage to the unit or marks to labels.
-
- Follow the Operating Instructions Microstim DB3 Tester as Found in the Document Index.
-
- Disconnect the Microstim, Remove the battery and refit the battery cover.