


Dimension Box Label: 150 mm x 127 mm

Ref.: 7090121001V, Type: 4420800
Single Port Nasal Prong, Adult, Short Term



Ref.: 7090121002V, Type: 4420801
Nasal Sampling Line, Adult, Short Term




Pouch Label



Box Label

Ref.: 7090121003V, Type: 4420802
Nasal Sampling Line, Paediatric, Short Term



Pouch Label



Box Label


VIAMED
 Viamed Limited, 15 Station Road, Cross Hills, Keighley
 West Yorkshire, BD20 7DT, United Kingdom
 www.viamed.co.uk, info@viamed.co.uk, +44(0)1535 634542

VersaStream®

Nasal Sampling Line, Infant, 2m/7ft

Short Term Use ≤ 24h

REF 7090121004V

LOT 22621.04.14

4420803
















UDI

MD

(01) 1 42516796 014 19

(11) 210414

(17) 240414

(10) 22621.04.14



< 10kg
 (< 22 lb)




bluepoint MEDICAL GmbH & Co.KG
 An der Trave 15
 23823 Glindeorf
 Germany




2021-04-14
 2024-04-14



A228v.3

Pouch Label

25 Nasal Sampling Line, disposable, Infant, 2m/7ft

REF 

7090121004V

LOT 

22621.04.14


cs Odběrová nosní hadička, jednorázová, kojenecká da Nasaludsamlingslange, engangsbrug, spædbarn de Einweg-Probennahmung, Nasal, Kleinkinder el Νίνα προβιτών, (μικροδότης) καταστή, βρέφους en Грудной провитель для новорожденных, носовик, новорожден es Línea de muestreo nasal, desechable, infantil fr Tubulure d'échantillonnage nasale, usage unique, nourisson ko 비강 검체 채취관, 일회용, 유아용 hr Linija za uzorkovanje nosnog sekreta, jednokratna, djeca it Linea di campionamento nasale, usa e getta, neonatale Iv Deguna parauga ņemšanas līnija, vienaizmantojama, zīdaijiem it Nosinė mėginų ėmimo linija, vienkartinė, kūdikiams hu Názális mintavételi vezeték, egyszer használatos, csecsemő nl Nasale monsteringslijn, voor eenmalig gebruik, zuigeling ja 鼻サンプライン管、使い捨て、幼児用	no Nasal prøvetakingsledning, til engangsbruk, spedbarn pl Nosowa linia próbkująca, jednorazowa, niemowlę pt Linha de amostragem nasal, descartável, infantil qj-av Linia de amostrogang nasal, descartável, bebê ro Linie de eşantionare nazală, de unică folosință, bebeluși ru Насальная пробоотборная линия, одноразовая, для младенцев sk Nazálna hadička na ober vzoriek, jednorázová, pre digitálnu sr Voz za uzorkovanje za nazalnu primenu, za jednokratnu upotrebu, za novorođenčad fi Nenäilyntälinja, kertakäyttöinen, pikkulapsel sv Nasal provtagningsledning, för engångsbruk, spädbarn tr Burun Numune Alma Hattı, tek kullanımlık, Bebek zh 鼻腔采样管，一次性，婴幼儿
---	--

4420803



10% 100% 100% 140°F 60°F

 ≤ 10 kg / 22 lbs



MD 

UDI 

21 02516790 0114 G

(1) 210414

(2) 240414

(10) 22621.04.14

Short Term Use ≤ 24h





2021-04-14

2024-04-14

blueprint medical GmbH & Co. KG

An der Trade 15

23923 Selmendorf, Germany

www.blueprintmedical.com

Box Label



VIAMED

Viamed Limited, 15 Station Road, Cross Hills, Keighley
West Yorkshire, BD20 7DT, United Kingdom
www.viamed.co.uk, info@viamed.co.uk, +44(0)1535 634542

VersaStream®

Nasal Sampling Line with O2 Supply, Adult, 2m/7ft

Short Term Use ≤ 24h

REF 7090121005V	LOT 22121.04.14	# 4420804	
-----------------	-----------------	-----------	--



10 %
35 %



1060 hPa
680 hPa



140 °F
60 °C
-13 °F
-23 °C

















Diagram illustrating the VersaStream device setup for nasal sampling. The device includes a MONITOR, a flowmeter (O₂ ≤ 5 l/min), a CO₂ sensor (CO₂ ≤ 75 ml/min), and a H₂O sensor. The patient is shown wearing the nasal cannula and the CO₂ sensor.

UDI

MD

(01) 1 42516796 015 16

(11) 210414

(17) 240414

(10) 22121.04.14

 ≥ 20kg
(≥ 44 lb)

CE

0482



bluepoint MEDICAL GmbH & Co. KG
An der Trave 15
23823 Glindeborf
Germany



2021-04-14

2024-04-14



Pouch Label

10 Nasal Sampling Line with O2 Supply, disposable, Adult, 2m/ftt

REF	7090121005V	LOT	22121.04.14
cs	Odběrná roura nasální s O ₂ , jednorázová, pro dospělé	no	Nasal pravtávlingslinje med O ₂ , 1t engangsbruk, voksen
da	Nasaludtagingslange med O ₂ , engangsbrug, voksen	pl	Nosowa linia próbkująca z O ₂ , jednorazowa, dorosły
de	Einweg-Probenabnehmer mit O ₂ , nasal, Erwachsene	pt	Linha de amostragem nasal com O ₂ , descartável, adulto
et	Nosa proovitor O ₂ -ga, ühekordne kasutus, täiskasvanu	ru	Линия от amostragem nasal com O ₂ , descartável, adulto
el	Γραμμή πινετης ρινικής συσκευής με O ₂ , ανατολική, ενήλικας	ro	Linie de esantionare nazală cu O ₂ , de unică folosință, adulți
es	Línea de muestreo nasal con O ₂ , desechable, adulto	ru	Назальная пробоотборная линия для O ₂ , одоразового, для взрослых
fr	Tubulaire d'échantillonnage nasale avec O ₂ , usage unique, adulte	sk	Násalná linka na odber vzorky s O ₂ , jednorázová, pre dospelých
ko	인소 비강 측정 장치, 일회용, 성인용	sr	Vod za uzorkovanje za nazalnu primenu sa O ₂ , za jednokratnu upotrebu, za odrasle
hr	Linija za uzorkovanje nosnog sekreta s O ₂ , jednokratna, odrasli	fi	Nenälyntälinja O ₂ n kanssa, kertakäyttöinen, aikuiset
it	Linea di campionamento nasale con O ₂ , usa e getta, pazienti adulti	sv	Nasal provtagningslinje med O ₂ , engangsbruk, vuxen
Iv	Duguna paragrau përmësimja linja e O ₂ , vlerëzimetqija, përdorqijem	tr	Burun Numune Alma Hattı, O ₂ li tek kullanımlık, Yetişkin
lt	Nosinė mėginų ėmimo linija su O ₂ , vienkartinė, suaugusiems	zh	带 O ₂ 的鼻采样管 (管), 一次性, 成人
hu	Násalis mintavételi vezeték oxigénnel, egyszer használatú, felnőtt	MD	(01) 2: 42516706 011 53
nl	Nasale bemonsteringslijn met O ₂ , voor eenmalig gebruik, volwassen	UDI	(11) 210414 (17) 240414 (10) 22121.04.14
ja	O ₂ 鼻サンプリングライン、使い捨て、成人用		

4420804  ≥ 20 kg / 44 lbs Short Term Use ≤ 24h

 10%  650 HPa  1060 HPa  140 °F  13 °F -25 °C           



 0482


Viamed Limited, 15 Station Road, Cross Hills, Keighley
 West Yorkshire, BD20 7DT, United Kingdom
 www.viamed.co.uk info@viamed.co.uk 01535 535452

2021-04-14
 2024-04-14

bluepoint medical GmbH & Co. KG
 An der Trave 15
 23923 Selmsdorf, Germany
 www.bluepointmedical.com

Box Label

Ref.: 7090121006V, Type: 4420805
Nasal Sampling Line with O2 Supply, Pediatric, Short Term




Pouch Label




Box Label

Ref.: 7090121007V, Type: 4420806
Nasal Sampling Line with O2 Supply, Infant, Short Term



Pouch Label



Box Label

Ref.: 7090121008V, Type: 4420807
Nasal Sampling Line, Adult, Long Term

Is identical with Ref.: 7090121002V, Type: 4420801
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121009V, Type: 4420808
Nasal Sampling Line, Pediatric, Long Term

Is identical with Ref.: 7090121003V, Type: 4420802
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121010V, Type: 4420809
Nasal Sampling Line, Infant, Long Term

Is identical with Ref.: 7090121004V, Type: 4420803
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121011V, Type: 4420810
Nasal Sampling Line with O2 Supply, Adult, Long Term

Is identical with Ref.: 7090121005V, Type: 4420804
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h


Ref.: 7090121012V, Type: 4420811
Nasal Sampling Line with O2 Supply, Pediatric, Long Term

Is identical with Ref.: 7090121006V, Type: 4420805
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121013V Type: 4420812
Nasal Sampling Line with O2 Supply, Infant, Long Term

Is identical with Ref.: 7090121007V, Type: 4420806
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121014V, Type: 4420813
Oral/Nasal Sampling Line, Adult, Short Term




Pouch Label



Box Label

Ref.: 7090121015V, Type: 4420814
Oral/Nasal Sampling Line, Pediatric, Short Term



Pouch Label



Box Label



VIAMED
 Viamed Limited, 15 Station Road, Cross Hills, Keighley
 West Yorkshire, BD20 7DT, United Kingdom
 www.viamed.co.uk, info@viamed.co.uk, +44(0)1535 634542

VersaStream®

Oral/Nasal Sampling Line with O₂ Supply, Adult, 2m/7ft

Short Term Use ≤ 24h

REF 7090121016V
 LOT 24121.04.14
 # 4420815


















O₂ ≤ 5 l/min
 MONITOR
 rO₂ ≤ 75 ml/min
 H₂O

UDI
MD



(01) 1 42516796 0126 2
 (11) 210414
 (17) 240414
 (10) 24121.04.14


**≥ 20kg
(≥ 44 lb)**

CE
0482



bluepoint MEDICAL GmbH & Co.KG
 An der Trave 15
 23023 Glindeorf
 Germany


2021-04-14

2024-04-14

REF



7090121016V

CS

Odslovná /istinnifunnið haddiða z O2, jednorázová, pro dospělosti

dra

Oral-haasutakutingslange med O2, engangsborg, voksen

de

Eintweg-Probenahmeleitung mit O2, oralnasal, Erwachsene

et

Suunina püüvituuri O2-ga, ühekordne kasutatav, täiskasvanu

el

Γραμμή συσπειρωμένης διαπνοής με O2, ανατολική, ενήλικες

es

Línea de muestreo nasal/oral con O2, desechable, adulto

fr

Tubulaire d'extraction nasale avec O2, usage unique, adulte

ko

신호 구멍에서 구멍 자른 것, 일회용, 성인용

hr

Linija za uzorkovanje slinešnog nosnog sekreta z O2, jednokratna, odrasli

it

Linea di campionamento orale/nasale con O2, usa e getta, pazienti adulti

lv

Muves/segura paraugu pērmāšana šūnās ar O2, vienreizējai lietojama, pieaugušajiem

it

Burnos / rosinės mėginėjimo linija su O2, vienkartinė, suaugusiems

hu

Orális/nazális mintavételi vezeték oxigénnel, egyszer használatos, felnőtt


nl

Orale/nasale bemonsteringslijn met O2, voor eenmalig gebruik, volwassen

ja

口・鼻サンプングライン、使い捨て、成人用

LOT



24121.04.14

no

Oralnasal prevaktingslinje med O2, én gangsborg, voksen

pl

Ustna/hosowa linia próbkująca z O2, jednorazowa, dorosły

pt

Linha de amostragem oral/nasal com O2, descartável, adulto

pt-br

Linha de amostragem oral/nasal com O2, descartável, adulto

ro

Linie de examinare oral/nazală cu O2, de unic folosință, adult

ru

Орально-назальный пробоотборный линия для O2, одоразового, для взрослых

sk

Orálna/nazálna hadička na odber vzorky z O2, jednorázová, pre dospelých

sr

Vod za uzorkovanje za oral/nazalnu primenu sa O2, za jednokratnu upotrebu, za odrasle

fi

Suu-/Nenäytäytinlinja O2:n kanssa, kertakäyttöinen, aikuiset

sv

Oral/nasalprovtagningledning med O2, för engångsborg, vuxen

tr

Ağız/Burun Numune Alma Hattı, O2'li, tek kullanımlık, Yetişkin

zh

带 O2 的口腔/鼻腔采样管，一次性，成人

MD



011) 242516796 0129 9
11) 210414
17) 240414
10) 2412.04.14

UDI



0482

CE

0482

4420815



≥ 20 kg / 44 lbs

Short Term Use ≤ 24h



100% RH

140 °F

50 °C

13 °F

-35 °C







VIAMED

Viamed Limited, 15 Station Road, Cross Hills, Kirkcaldy
West Yorkshire, B20 7DT, United Kingdom
www.viamed.co.uk, info@viamed.co.uk, 01937 534542



2021-04-14



2024-04-14



blueprint medical GmbH & Co. KG

An der Trave 15

23933 Selsdorf, Germany

www.blueprintmedical.com



VIAMED

Viamed Limited, 15 Station Road, Cross Hills, Keighley
 West Yorkshire, BD20 7DT, United Kingdom
www.viamed.co.uk, info@viamed.co.uk +44(0)1535 634542

VersaStream®

Oral/Nasal Sampling Line with O2 Supply, Paediatric, 2m/ftt

Short Term Use ≤ 24h

REF 7090121017V

LOT 24221.04.14

4420816



10 L/min 35 L/min



600 hPa 1050 hPa



-13 °C +25 °C 140 °F 68 °F




















UDI

MD



(01) 1 42516796 012 9
 (11) 210414
 (17) 240414
 (10) 24221.04.14



≥ 10kg
(≥ 22 lb)



0482



bluepoint MEDICAL GmbH & Co.KG
 An der Trave 15
 23823 Glindeorf
 Germany

2021-04-14
 2021-04-14



[illegible]

Company registered in England, No. 01291765

Ref.: 7090121018V, Type: 4420817
Nasal Sampling Line, Adult, Long Term

Is identical with Ref.: 7090121014V, Type: 4420813
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121019V, Type: 4420818
Nasal Sampling Line, Pediatric, Long Term

Is identical with Ref.: 7090121015, Type: 4420814
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h


Ref.: 7090121020, Type: 4420819
Nasal Sampling Line with O2 Supply, Adult, Long Term

Is identical with Ref.: 7090121016V, Type: 4420815
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121021V, Type: 4420820
Nasal Sampling Line with O2 Supply, Pediatric, Long Term

Is identical with Ref.: 7090121017V, Type: 4420816
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121022V, Type: 4420821
Airway Adapter Sampling Line, Adult/Pediatric, Short Term




Pouch Label



Box Label

Ref.: 7090121023V, Type: 4420822
Airway Adapter Sampling Line, Infant, Short Term



Pouch Label



Box Label

Ref.: 7090121024, Type: 4420823

Airway Adapter Sampling Line, Adult/Pediatric, Long Term

Is identical with Ref.: 7090121022V, Type: 4420821

Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121025V, Type: 4420824

Airway Adapter Sampling Line, Infant, Long Term

Is identical with Ref.: 7090121023V, Type: 4420822

Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

2. Label Layout: LLM Connector Sampling Lines

Ref.: 7090121033V, Type: 4420857
Nasal Sampling Line, Adult, Long Term

Is identical with Ref.: 7090121002V, Type: 4420801
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121034V, Type: 4420858
Nasal Sampling Line, Pediatric, Long Term

Is identical with Ref.: 7090121003V, Type: 4420802
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121035V, Type: 4420859
Nasal Sampling Line, Infant, Long Term

Is identical with Ref.: 7090121004V, Type: 4420803
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121036V, Type: 4420860
Nasal Sampling Line with O2 Supply, Adult, Long Term

Is identical with Ref.: 7090121005V, Type: 4420804
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121037V, Type: 4420861
Nasal Sampling Line with O2 Supply, Pediatric, Long Term

Is identical with Ref.: 7090121006V, Type: 4420805
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121038V, Type: 4420862
Nasal Sampling Line with O2 Supply, Infant, Long Term

Is identical with Ref.: 7090121007V, Type: 4420806
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121043V, Type: 4420867
Nasal Sampling Line, Adult, Long Term

Is identical with Ref.: 7090121014V, Type: 4420813
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121044V, Type: 4420868
Nasal Sampling Line, Pediatric, Long Term

Is identical with Ref.: 7090121015V, Type: 4420814
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121045V, Type: 4420869
Nasal Sampling Line with O2 Supply, Adult, Long Term

Is identical with Ref.: 7090121016V, Type: 4420815
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121046V, Type: 4420870
Nasal Sampling Line with O2 Supply, Pediatric, Long Term

Is identical with Ref.: 7090121017V, Type: 4420816
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121049V, Type: 4420873
Airway Adapter Sampling Line, Adult/Pediatric, Long Term

Is identical with Ref.: 7090121022V, Type: 4420821
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121050V, Type: 4420874
Airway Adapter Sampling Line, Infant, Long Term

Is identical with Ref.: 7090121023V, Type: 4420822
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

3. Label Layout: ORI Connector Sampling Lines

Ref.: 7090121102V, Type: 4420900
Single Port Nasal Prong, Adult, Short Term

Is identical with Ref.: 7090121001V, Type: 4420800

Ref.: 7090121103V, Type: 4420901
Nasal Sampling Line, Adult, Short Term

Is identical with Ref.: 7090121002V, Type: 4420801

Ref.: 7090121104V, Type: 4420902
Nasal Sampling Line, Pediatric, Short Term

Is identical with Ref.: 7090121003V, Type: 4420802

Ref.: 7090121105V, Type: 4420903
Nasal Sampling Line, Infant, Short Term

Is identical with Ref.: 7090121004V, Type: 4420803

Ref.: 7090121106V, Type: 4420904
Nasal Sampling Line with O2 Supply, Adult, Short Term

Is identical with Ref.: 7090121005V, Type: 4420804

Ref.: 7090121107V, Type: 4420905
Nasal Sampling Line with O2 Supply, Pediatric, Short Term

Is identical with Ref.: 7090121006V, Type: 4420805

Ref.: 7090121108V, Type: 4420906
Nasal Sampling Line with O2 Supply, Infant, Short Term

Is identical with Ref.: 7090121007V, Type: 4420806

Ref.: 7090121109V, Type: 4420907
Nasal Sampling Line, Adult, Long Term

Is identical with Ref.: 7090121002V, Type: 4420807
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121110V, Type: 4420908
Nasal Sampling Line, Pediatric, Long Term

Is identical with Ref.: 7090121003V, Type: 4420808
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121111V, Type: 4420909
Nasal Sampling Line, Infant, Long Term

Is identical with Ref.: 7090121004V, Type: 4420809
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121112V, Type: 4420910
Nasal Sampling Line with O2 Supply, Adult, Long Term

Is identical with Ref.: 7090121005V, Type: 4420810
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121113V, Type: 4420911
Nasal Sampling Line with O2 Supply, Pediatric, Long Term

Is identical with Ref.: 7090121006V, Type: 4420811
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121114V, Type: 4420912
Nasal Sampling Line with O2 Supply, Infant, Long Term

Is identical with Ref.: 7090121007V, Type: 4420812
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121115V, Type: 4420913
Oral/Nasal Sampling Line, Adult, Short Term

Ref.: 7090121116V, Type: 4420914
Oral/Nasal Sampling Line, Pediatric, Short Term

Is identical with Ref.: 7090121015V, Type: 4420814

Ref.: 7090121117V, Type: 4420915
Oral/Nasal Sampling Line with O2 Supply, Adult, Short Term

Is identical with Ref.: 7090121016V, Type: 4420815

Ref.: 7090121118V, Type: 4420916
Oral/Nasal Sampling Line with O2 Supply, Pediatric, Short Term

Is identical with Ref.: 7090121017V, Type: 4420816

Ref.: 7090121119V, Type: 4420917
Nasal Sampling Line, Adult, Long Term

Is identical with Ref.: 7090121014V, Type: 4420817
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121120V, Type: 4420918
Nasal Sampling Line, Pediatric, Long Term

Is identical with Ref.: 7090121015V, Type: 4420818
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 70901210121V, Type: 4420919
Nasal Sampling Line with O2 Supply, Adult, Long Term

Is identical with Ref.: 7090121016V, Type: 4420819
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121122V, Type: 4420920
Nasal Sampling Line with O2 Supply, Pediatric, Long Term

Is identical with Ref.: 7090121017V, Type: 4420820
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121123V, Type: 4420921
Airway Adapter Sampling Line, Adult/Pediatric, Short Term

Is identical with Ref.: 7090121022V, Type: 4420821

Ref.: 7090121124V, Type: 4420922
Airway Adapter Sampling Line, Infant, Short Term

Is identical with Ref.: 7090121023V, Type: 4420822

Ref.: 7090121125V, Type: 4420923
Airway Adapter Sampling Line, Adult/Pediatric, Long Term

Is identical with Ref.: 7090121022V, Type: 4420823
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121126V, Type: 4420924
Airway Adapter Sampling Line, Infant, Long Term

Is identical with Ref.: 7090121023V, Type: 4420824
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h


Ref.: 7090121140V, Type: 4420925
Airway Adapter Sampling Line, Adult/Pediatric, Long Term

Is identical with Ref.: 7090121022V, Type: 4420825
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h High Humidity


Ref.: 7090121141V, Type: 4420926
Airway Adapter Sampling Line, Infant, Long Term

Is identical with Ref.: 7090121023V, Type: 4420826
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h High Humidity

Tag Layout:



Tag for short term sampling lines



Tag for long term sampling lines